

Sep 30, 2023

Plasmid Construction Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldm5nol5b/v1>

Xuefeng Ren^{1,2}, Annan Cook^{1,2}

¹University of California, Berkeley; ²Aligning Science Across Parkinson's



Annan SI Cook

University of California, Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldm5nol5b/v1>

Protocol Citation: Xuefeng Ren, Annan Cook 2023. Plasmid Construction Protocol. protocols.io
<https://dx.doi.org/10.17504/protocols.io.n92ldm5nol5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 25, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 88511

Keywords: ASAPCRN, pcag mammalian expression vector, plasmid construction protocol, plasmid construction protocol this protocol detail, insertion of gene, plasmid, gene, insertion

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

Abstract

This protocol details plasmid construction in general terms for the insertion of genes into the pCAG mammalian expression vector.

Attachments



[852-2204.pdf](#)

31KB

Materials

Materials

- DNA sequences encoding PI3KC3-C1 subunits and VPS15 mutants (codon-optimized and synthesized)
- pCAG vector
- Restriction enzymes (e.g., New England Biolabs enzymes)
- Gibson assembly kit
- Sanger sequencing service



Design Plasmid Constructs

- 1 This protocol uses the pCAG vector and NEB restriction enzymes and Gibson Assembly kits.
- 2 Obtain the DNA sequences for the genes of interest from your preferred DNA synthesis vendor.

Digestion or Gibson Assembly

- 3 Decide whether to use restriction digestion or Gibson assembly for subcloning.
- 4 **Restriction Digestion:**
 - 4.1 Digest both the pCAG vector and the DNA fragments (genes/tags) using appropriate restriction enzymes.
 - 4.2 Purify the digested fragments using a DNA purification kit.
- 5 **Gibson Assembly:**

Follow the manufacturer's protocol for the Gibson assembly kit to assemble the DNA fragments into the pCAG vector.



Ligation

- 6 If using restriction digestion, perform a ligation reaction to insert the digested DNA fragments into the linearized pCAG vector.
- 7 Use a DNA ligation kit and follow the manufacturer's instructions.

Transformation

- 8 Transform the ligated plasmids into competent *E. coli* cells.



- 9 Plate the transformed cells on selective agar plates containing the appropriate antibiotic for plasmid selection.

Bacterial Culture

- 10 Incubate the plates  Overnight at an appropriate temperature for *E. coli* growth (e.g.,  37 °C).



Plasmid Isolation

- 11 Pick colonies that have grown on the plates and inoculate them into liquid culture with the same antibiotic.
- 12 Grow the cultures to obtain a sufficient amount of plasmid DNA.
- 13 Isolate the plasmid DNA using a plasmid purification kit.
- 14 Sanger sequence the results.