

Nov 27, 2022

## Plasmid Expansion

DOI

[dx.doi.org/10.17504/protocols.io.8epv5jrjl1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5jrjl1b/v1)

Ali Albalakhi<sup>1,2</sup>, Ning Xia<sup>1,2</sup>

<sup>1</sup>Massachusetts General Hospital; <sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Jackson Schumacher

Massachusetts General Hospital

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5jrjl1b/v1>

**Protocol Citation:** Ali Albalakhi, Ning Xia 2022. Plasmid Expansion. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5jrjl1b/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working.

**Created:** November 27, 2022

**Last Modified:** December 12, 2024

**Protocol Integer ID:** 73277

**Keywords:** ASAPCRN, protocol for plasmid expansion, plasmid expansion, protocol

## Abstract

This is the protocol for plasmid expansion.

## Attachments



Plasmid\_Expansion\_AL..

49KB

## Troubleshooting



- 1 Thaw the competent cells on ice
- 2 Chill approximately 5 µg (2 µl) of Plasmid DNA in 1.5 mL microcentrifuge tube
- 3 Add 50 µl of competent cells to the DNA. Mix gently by pipetting up and down or flicking the tube 4–5 times to mix the cells and DNA. Do not vortex.
- 4 Place the mixture on ice for 30 minutes. Do not mix.
- 5 Heat shock at 42°C for 30 seconds\*. Do not mix.
- 6 Add 950 µl of room temperature media\* to the tube.
- 7 Place tube at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.
- 8 Warm selection plates to 37°C.
- 9 Spread 50–100 µl of the cells and ligation mixture onto the plates.
- 10 Incubate overnight at 37°C.
- 10.1 Please note: For the duration and temperature of the heat shock step as well as for the media to be used during the recovery period, please follow the recommendations provided by the competent cells' manufacturer.

## Plasmids and antibiotics

---

11



A	B
PLVX mCherry (HLA-A2)	Amp 50µg/mL
Ploc MLANA (MART1)	Blasticidin (50µg/mL)
psPAX2	Amp 50µg/mL
VSV.G	Amp 50µg/mL

## Preparation of Liquid Media - LB-Ampicillin Broth and LB-Blasticidin Broth:

- 12 Add 20g of LB broth Lennox into 1 liter of distilled water in an Erlenmeyer flask
- 13 Autoclave the LB broth in the Flask on a wet cycle for 30min
- 14 Storage of LB in 4°C is recommended and is good for up to 1 year
- 15 Ensure that the LB broth has come to room temperature before you add ampicillin
- 16 The final concentration of ampicillin in the broth should be 50µg/mL
- 17 LB broth containing ampicillin should be stored at 4°C for up to 1 month

## Preparation of Solid Media

- 18 LB-Amp Agar and LB-Blasticidin Agar: Each Petri Dish takes about 10mL of LB-Agar so scale the volume accordingly.  
  
1 LB agar tablet makes 50mL
- 19 Take 1 tablet and dissolve it in 50mL of distilled water in an autoclaved flask
- 20 Autoclave the LB agar for 30min on a wet cycle



- 21 Once the cycle is complete check to make sure that all the agar is dissolved
- 22 Allow the LB agar to cool until it is comfortable to the touch (50°C), a water bath set to 50°C is useful for this step
- 23 While the agar is cooling label the plates using sharpie for each antibiotic
  - 23.1 a. Ampicillin = Red
  - b. Blastcidin = Blue
- 24 The final concentration of ampicillin and Blastcidin should be 50µg/mL
- 25 Make sure that the LB agar has cooled to 50°C as excessive heat will degrade the antibiotics