

Sep 28, 2023

Plasmid Extraction (Plasmid Isolation)

DOI

dx.doi.org/10.17504/protocols.io.yxmvm3yq6l3p/v1

NUS iGEM¹

¹National University of Singapore



NUS iGEM

National University of Singapore

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm3yq6l3p/v1>

Protocol Citation: NUS iGEM 2023. Plasmid Extraction (Plasmid Isolation). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvm3yq6l3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 27, 2023



Last Modified: September 28, 2023

Protocol Integer ID: 88484

Keywords: Plasmid Extraction, Plasmid, Buffer P1, Buffer P2, Buffer N3, Buffer PB, Buffer PE, Plasmid Isolation, Isolation, plasmid isolation, plasmid extraction, plasmid from the cell, plasmid, singapore igem team, isolation, extraction

Abstract

2023 NUS-Singapore iGEM team followed this protocol to isolate the plasmid from the cells.

Protocol materials

⊗ Buffer N3 **Qiagen Catalog #19064**

⊗ PB buffer **Qiagen Catalog #19066**

⊗ Buffer P2 **Qiagen Catalog #19052**

⊗ Buffer PE **Qiagen Catalog #19065**

⊗ Buffer P1 **Qiagen Catalog #19051**

⊗ RNase A **Qiagen Catalog #19101**




⊗ Buffer PE **Qiagen Catalog #19065**

Troubleshooting

Safety warnings





- ! ■ Proper lab PPE must be worn at all times.

Before start

Glycerol stock ( 900 µL of cell stock in  300 µL of 100% glycerol solution) may be prepared and keep it in the  -80 °C fridge before the plasmid extraction procedure.









Cell Culture












- 1 Incubate the cells containing the plasmid of interest in a Falcon tube with  5 mL of LB media and  5 μ L of the appropriate antibiotics  Overnight at  37 °C before starting the plasmid extraction.

Plasmid Extraction

5m

- 2 Take out the Falcon tube with the cultured cells from the incubator.
- 3 Centrifuge the Falcon tube at  5000 rpm, 4°C, 00:05:00 .
- 4 Discard the supernatant in the Falcon tube and keep the cell pellet.
- 5 Add  250 μ L of  Buffer P1 **Qiagen Catalog #19051** (with  RNase A **Qiagen Catalog #19101** added, kept in the  4 °C fridge) into the Falcon tube and resuspend the cell pellet.
- 6 Transfer the whole solution into a Eppendorf tube.
- 7 Add  250 μ L of  Buffer P2 **Qiagen Catalog #19052** into the Eppendorf tube to lyse the cells.
- 8 Shake the Eppendorf tube for mixing.
- 9 Add 350 μ L of  Buffer N3 **Qiagen Catalog #19064** into the Eppendorf tube.
- 10 Shake the Eppendorf tube for mixing, a cloudy solution shall be observed.



- 11 Centrifuge the Eppendorf tube at  13 rpm, 00:10:00 . 10m
- 12 Transfer the supernatant into a Mini Prep tube and discard the cell debris.
- 13 Centrifuge it at  13 rpm, 00:01:00 . 1m
- 14 Discard the flow-through and place the Mini Prep tube back into the same tube.
- 15 Add  500 μL of  PB buffer **Qiagen Catalog #19066** into the Mini Prep tube.
- 16 Centrifuge it at  13 rpm, 00:01:00 .
- 17 Discard the flow-through and place the Mini Prep tube back into the same tube.
- 18 Add  700 μL of  Buffer PE **Qiagen Catalog #19065** into the Mini Prep tube.
- 19 Centrifuge it at  13 rpm, 00:01:00 .
- 20 Discard the flow-through and place the Mini Prep tube back into the same tube.
- 21 Centrifuge it at  13 rpm, 00:01:00 to remove the  Buffer PE **Qiagen Catalog #19065** residual.
- 22 Transfer the Mini Prep column into the newly labelled Eppendorf tube.
- 23 Add  50 μL of DI water into the Eppendorf tube.



- 24
- Centrifuge it at

13 rpm, 00:01:00

, ensuring that the direction of the Eppendorf tube's cap is the same as the direction of spinning to avoid breaking.
- 25
- Discard the Mini Prep column, the solution left in the Eppendorf tube is the plasmid.
- 26
- Use NanoDrop

Equipment	
NanoDrop™ One/OneC Microvolume UV-Vis Spectrophotometer	NAME
UV-Vis Spectrophotometer	TYPE
Thermo Scientific	BRAND
ND-ONE-W	SKU

-
- to measure the concentration (under "dsDNA" mode) and the purity of the plasmid.
- 27
- Store the isolated plasmid at room temperature.