

Oct 11, 2022

Plasmid Miniprep for 2022



Forked from [Plasmid Miniprep](#)

DOI

dx.doi.org/10.17504/protocols.io.dm6gpjwdpgzp/v1

Team Fudan iGEM¹

¹Fudan University

2022



Team Fudan iGEM

Fudan University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpjwdpgzp/v1>

Protocol Citation: Team Fudan iGEM 2022. Plasmid Miniprep for 2022. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gpjwdpgzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 10, 2022

Last Modified: October 11, 2022

Protocol Integer ID: 71085





Keywords: plasmid miniprep, plasmid extraction, purification

Abstract

Plasmid extraction and purification, based on Vazyme Mini-prep Kit manual.

Troubleshooting



- 1 Pellet bacteria culture in a 1.6mL centrifugation tube at 13,000rpm for 1 minute. Decant supernatant. The pellet in the tube can be stored in -20 freezer if time not enough for the following steps.
- 2 Add 250uL of resuspension solution (P1 with RNase) and vortex. After adding RNase into the P1, P1 is stored in 4 degree refrigerator.
- 3 Add 250uL of lysis solution (P2) and invert gently 4-6 times. P2 contains NaOH, and must keep tightly closed.
- 4 Add 350uL of neutralization solution (P3) and invert gently 4-6 times. Do not vortex, which will break bacteria genomic DNA and genomic DNA fragments will contaminate purified plasmids.
- 5 Spin at 13,000rpm for 10 minutes, at 4 degree.
- 6 Take the clear supernatant, without any white protein cloud. Bind the DNA to the column by decanting.
- 7 Centrifuge at 13,300 rpm for 1 minute.
 00:01:00
- 8 Wash the column with 600uL of wash solution (PW2) and centrifuge for 1 minute. Remind adding pure ethanol into PW2 as indicated on the bottle.
 00:01:00
- 9 Discard the flow-through, wash again with 600uL of PW2, and centrifuge for 1 minute. Discard the flow-through and centrifuge the column again for 2 minutes.
 00:03:00
- 10 Keep the cap open and wait until no ethanol smell. Move the column into a new centrifugation tube. Elute the purified DNA with 50uL of elution buffer (TE; 10 mM Tris-Cl, 1 mM EDTA). Then, centrifuge for 1 minutes. Collect the flow-through, which is the plasmid DNA. Because pET28 based is low-copy plasmid, pre-warm TE to 50 degree helps the elution.
 00:01:00



3m