

Nov 14, 2023

PlatereaderPCAoxidationAssay

DOI

dx.doi.org/10.17504/protocols.io.bp2l6xm6dlqe/v1

Lev Tsypin¹, Dianne K Newman², Allen W Chen², Scott Saunders³

¹Stanford University; ²California Institute of Technology; ³UT Southwestern



Lev Tsypin

Stanford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6xm6dlqe/v1>

Protocol Citation: Lev Tsypin, Dianne K Newman, Allen W Chen, Scott Saunders 2023. PlatereaderPCAoxidationAssay. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l6xm6dlqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 14, 2023

Last Modified: November 14, 2023

Protocol Integer ID: 90896

Keywords: plate reader

Funders Acknowledgements:

NIH

Grant ID: 1R01AI127850-01A1

Doren Family Foundation

NSF GRFP

Abstract

Protocol for measuring phenazine-1-carboxylic acid oxidation in a plate reader

Troubleshooting

Pre-grow cultures for assay

- 1 Two days before, streak frozen culture stock on LB agar plate, grow over night at 30 °C; inoculate liquid cultures for next step with a mixed patch of cells.
- 2 Grow 5 mL cultures of cells overnight at 30 °C in LB for 17 hours +/- 10 minutes.
- 2.1 Depending on the experiment, tubes are either slanted shaking at 250 rpm or standing non-shaking with caps sealed with Parafilm

Prepare cells for assay

- 3 Wash entire volume into basal medium with no terminal electron acceptor and no PCA (essentially just buffer and salts): 20 mM potassium phosphate buffer (pH 7-7.1); 1 mM sodium sulfate; 10 mM ammonium chloride; and 1× freshwater salt solution (17.1 mM sodium chloride, 1.97 mM magnesium chloride, 0.68 mM calcium chloride, and 6.71 mM potassium chloride)
- 3.1 In 1 mL aliquots, wash 3x into basal medium by spinning for 2 min at 6000x(g), aspirating with vacuum trap, and resuspending by pipetting.
- 3.2 Measure OD600 and normalize all cultures to OD600 = 0.2-1, target OD600 = 0.5.
- 3.3 For slow growing cultures, like *menAubiC-tlKO*, this amounts to spinning down 1 mL of the overnight culture and resuspending it in 250 µL, then adjusting the volume. May need to scale up if inoculating a lot of wells.
- 4 Bring cultures into the anaerobic chamber.
- 5 Transfer washed cultures to anoxic microcentrifuge tubes (tubes that have been in chamber for at least three days).
- 6 Let stand for at least 1 hour. To test that this is fine, can track parallel culture with resazurin to see when it turns pink or clear.

Set up plate (BRAND PureGrade™ S, Cat. No. 781671) inside anoxic chamber



- 7 Do this as much as possible while cells are incubating in the anoxic chamber
- 8 Prepare reduced PCA calibration wells
- 8.1 Prepare 500 μM solution by diluting PCAred stock in basal medium.
- 8.2 In one row of the plate, prepare calibration (in μM): 250, 200, 175, 150, 125, 100, 75, 50, 25, 10, 5, 0.

This corresponds to the following volumes of the 500 μM PCAred per well (in μL): 100, 80, 70, 60, 50, 40, 30, 20, 10, 4, 2, 0.
- 8.3 Bring total volume in each well to 200 μL with basal medium
- 9 Desired concentration of PCA in PCA wells is 200 μM (40 μL of 1 mM stock).
- 10 Desired concentration of cells is 40 μL of OD600 = 0.5 culture for a target of OD600 = 0.1 in the wells.
- 11 Other components depend on the experiment: e.g., 2 μL of 1 M sodium nitrate for a 10 mM experimental concentration.
- 12 Note: all solutions prepared in the same basal medium.
- 13 Final volume in each well is 200 μL .
- 14 Order of preparing wells (except calibration): basal medium, then PCA, then TEA, then cells (always last).

Plate reader protocol (BioTek Synergy 4 or HTX)

- 15 Incubate at 30 $^{\circ}\text{C}$ throughout experiment. If using HTX model, set the temperature limits +/- 1 $^{\circ}\text{C}$ to prevent condensation



- 16 Medium shaking throughout experiment
- 17 Measure absorbance at 440 and 600 nm (PCA and cells, respectively), as well as fluorescence (360/40 ex and 528/20 em). Test that sensitivity is such that the calibration curve spans the dynamic range (different on each instrument). Include pathlength correction for absorbance measurements.
- 18 Measure every 5 minutes