

Oct 01, 2024

pMOA 189F-682R

DOI

dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1

Roey Angel¹, Eva Petrova¹

¹Institute of Soil Biology and Biogeochemistry, Biology Centre CAS

Anaerobic and Molecular Microbiology Lab, Biology Centre CAS

Tech. support email: eva.petrova@bc.cas.cz



Eva Petrova

Soil and Water Research Infrastructure

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1>

Protocol Citation: Roey Angel, Eva Petrova 2024. pMOA 189F-682R. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1>

Manuscript citation:

Holmes AJ, Costello A, Lidstrom ME, Murrell JC. Evidence that particulate methane monooxygenase and ammonia monooxygenase may be evolutionarily related. FEMS Microbiol Lett. 1995;132:203–8.

<http://www.sciencedirect.com/science/article/pii/037810979500311R>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 24, 2018

Last Modified: October 01, 2024

Protocol Integer ID: 12426

Keywords: PCR, pmoA, PMMO, Methanotrophs, Methane oxidising bacteria, Particulate methane monooxygenase, particulate methane monooxygenase gene, alpha subunit of the particulate methane monooxygenase gene, amoA gene of ammonia, general assay for methane, pcr assay, oxidising bacteria, amoA gene, methane, ammonia, pmoa, pcr, general assay, bacteria, assay, using primer, alpha subunit, 682r, primer

Abstract

A PCR assay targeting the alpha subunit of the particulate methane monooxygenase gene (*pmoA*). This is a general assay for methane oxidising bacteria using primers 189F–682R.

Note: these primers might occasionally also amplify the *amoA* gene of ammonia-oxidising bacteria.

189F GGN GAC TGG GAC TTC TGG (target seq. 189 - 207) Holmes et al. (1995),
682R GAA SGC NGA GAA GAA SGC (target seq. 664 - 682) Holmes et al. (1995).

Fragment size: 508bp.

Troubleshooting



PCR mixture

1

Reagent	Final. conc.	1 tube (50µl)	1 tube(25 µl)	96 tubes (25µl x100)	96 tubes (12.5µl x100)
PCR H ₂ O		37.9	18.95	1255	627.5
10X Buffer (I for RT and II for DNA)	1x	5	2.5	500 (Promega)	250 (Promega)
dNTP mixture (2.0 mM each	0.2mM	5	2.5	250	125
MgCl ₂ (50 mM)*	3mM	1.5*	0.75*	250	125
BSA (20µg/µl)	0.8 µg/µl	2	1	100	50
189f (25µM)	0.4 µM	0.8	0.4	40	20
682r (25 µM)	0.4 µM	0.8	0.4	40	20
AccuPrime™ Taq	?	1	0.5	15	7.5
Template		1	0.5	0.5 × 100	0.25 × 100
Final volume		50	25	2500	1250

* Buffer contains 1.5 mM MgCl₂ at final concentration

PCR program



2

1. $94^{\circ}\text{C} - 4'$

2. x **10** {

a. $94^{\circ}\text{C} - 1'$

b. $62^{\circ}\text{C} - 1^{\circ}\text{C} - 1'$

c. $68^{\circ}\text{C} - 1'$

}

3. x **25** {

a. $94^{\circ}\text{C} - 1'$

b. $52^{\circ}\text{C} - 1'$

c. $68^{\circ}\text{C} - 1'$

}

4. $68^{\circ}\text{C} - 10'$