


Aug 26, 2024

Pollen germination on wheat stigmas

DOI

dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1



Marina Millán Blázquez¹

¹John Innes Centre

Uauy Lab



Marina Millán Blázquez

John Innes Centre

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1>

Protocol Citation: Marina Millán Blázquez 2024. Pollen germination on wheat stigmas . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2024

Last Modified: August 26, 2024

Protocol Integer ID: 106332

Keywords: pollen germination, stigma, wheat, aniline blue, pollen germination on wheat stigma, aniline blue staining of pollen tube growth, pollen tube growth, wheat stigma, pollination, min after pollination, blue staining

Abstract

Aniline blue staining of pollen tube growth on the wheat stigma and through the style 4 h and 30 min after pollination.

Troubleshooting

Hand pollination

- 1 Please refer to [wheat-training.com](https://www.wheat-training.com/wp-content/uploads/Wheat_growth/pdfs/How-to-cross-wheat-pdf.pdf) for detailed explanations on how to perform hand pollinations. Link to pdf: https://www.wheat-training.com/wp-content/uploads/Wheat_growth/pdfs/How-to-cross-wheat-pdf.pdf

Carpel dissection and fixation

- 2 Using a pair of tweezers, dissect carpels 4.5 h after pollination to allow sufficient time for pollen tube emergence.
- 3 Store samples in a fixative solution of 95% ethanol and absolute acetic acid (75% v/v) and kept at 4 °C until sample preparation for fluorescence microscopy.

Aniline blue staining of pollinated stigmas

- 4 On the day of sampling, prepare a solution of 0.1% aniline blue in 0.1 M K_3PO_4 . You can prepare a stock solution of 1% aniline blue dissolved in 1x PBS. The stock solution should be kept in the fridge at 4 °C. Use tin foil to avoid exposure to light.
- 5 Wash fixed samples three times for 5 minutes in sterile water and transferred to 0.1% aniline blue solution and kept overnight at 4 °C. Use tin foil to avoid exposure to light.
- 6 Without washing the samples, dissect out the ovary using a sharp razor blade. Try not to damage the stigma. Leave the remaining stigmatic tissue to dry at 45 °C in a hot plate for a few minutes until most of the aniline blue solution is evaporated. Cover the hot plate with an opaque lid to avoid exposure to light.
- 7 Use Vectashield (catalogue No. H-1000-1, 2BSCIENTIFIC LTD) as an antifade mounting media to preserve fluorescence.
- 8 Happy microscopy.