

Aug 16, 2018 Version 2

## Pollen metabarcoding V.2

DOI

[dx.doi.org/10.17504/protocols.io.sn8edhw](https://dx.doi.org/10.17504/protocols.io.sn8edhw)

Tomasz Suchan<sup>1</sup>

<sup>1</sup>Institute of Botany, Polish Academy of Sciences

Molecular Biogeography...



**Tomasz Suchan**

W. Szafer Institute of Botany, Polish Academy of Sciences

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.sn8edhw](https://dx.doi.org/10.17504/protocols.io.sn8edhw)

**Protocol Citation:** Tomasz Suchan 2018. Pollen metabarcoding. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.sn8edhw>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
correcting some typos

**Created:** August 16, 2018

**Last Modified:** August 16, 2018

**Protocol Integer ID:** 14784

## Guidelines

Perform reactions in small batches until you are confident that there is no cross-contamination among the samples. Including isolation blanks and PCR blanks is crucial for the quality control.



## Materials

### MATERIALS

✂ Q5 Hot Start High-Fidelity DNA Polymerase - 100 units **New England Biolabs Catalog #M0493S**

✂ Water, nuclease free

✂ dNTP mix (25 mM of each)

✂ Phire Plant Direct PCR Kit **Thermo Fisher Scientific Catalog #F130WH**

### Before start

Prepare 5  $\mu$ M primer solutions:

ITS2 primers used in the 1st PCR:

ITS2-4R GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTNNNNNNNTCCTCCGCTTATTGATATGC

ITS2-S2F ACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNNNATGCGATACTTGGTGTGAAT

Indexing primers used in the 2nd PCR (xxxxxxxx - index):

AATGATACGGCGACCACCGAGATCTACACxxxxxxxxACACTCTTTCCCTACACGACGC

CAAGCAGAAGACGGCATACGAGATxxxxxxxxGTGACTGGAGTTCAGACGTGTGC

## Pollen extraction

- 1 Vortex the butterfly in 50 µl of water with 0.1% SDS.

### Safety information

Add blank sample at this step (= "isolation blank")  
<http://null>

- 2 Evaporate water in speedvac.
- 3 Add 5 µl of the Phire Plant Direct sample buffer.
- 4 Spin max speed for 2 min.

## 1st PCR

- 5 Prepare the mix:
  - 🧴 14 µL molecular grade water
  - 🧴 25 µL Phire Plant Direct PCR mix
  - 🧴 5 µL ITS2-S2F primer
  - 🧴 5 µL ITS2-4R primer
- 6 Add 1 µl of the sample to 49 µl of the mix. Use water instead of the sample for the blanks.

### Safety information

Add another blank sample at this step (= "PCR blank")

- 7 Run the PCR program: initial denaturation at 98°C for 5 min; 20 cycles of denaturation at 98°C for 40 s, annealing at 49°C for 40 s and elongation at 72°C for 40 s; followed by a



final extension step at 72°C for 5 min.

## Purification

- 8 Perform AMPure purification with ratio 1x. Elute in 10 µl.

### Protocol



NAME

**AMPure purification protocol**

CREATED BY

Tomasz Suchan

**PREVIEW**

- 8.1 Remove beads from the fridge, equilibrate to the room temperature, mix thoroughly but do not vortex.
- 8.2 Add desired ratio of AMPure beads to the DNA to be purified and mix well by pipetting.
- 8.3 Incubate 5 minutes.
- 8.4 Place on the magnetic rack.
- 8.5 Let it stand for 5 minutes on the rack, aspirate and discard supernatant.
- 8.6 Add 200 µl of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.
- 8.7 Repeat the wash: add 200 µl of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.



- 8.8 Wait until the pellet dries out completely. DO NOT remove the plate/tubes from the magnetic rack, as the dried AMPure might pop out from the tubes!
- 8.9 Add desired volume of Tris 10 mM or water, wetting the dried AMPure pellets (add 1  $\mu$ L to the final volume to avoid pipetting out the beads)
- 8.10 Remove from the magnetic rack.
- 8.11 Resuspend by pipetting or vortexing.
- 8.12 Incubate 10 minutes, incubating in 37°C can improve DNA yield.
- 8.13 Place on the magnetic rack.
- 8.14 Let it stand for 5 minutes, pipette out and save supernatant. The eluted DNA is in the supernatant, do not discard it!

## 2nd PCR

- 9 Prepare the mix:

🧪 4.82  $\mu$ L molecular grade water

🧪 2  $\mu$ L Q5 reaction buffer

🧪 0.08  $\mu$ L dNTPs (25 mM each)

🧪 0.1  $\mu$ L Q5 Hot Start polymerase

- 10 Add 1  $\mu$ L of the template to 7  $\mu$ L of the mix.

- 11 Add 2  $\mu$ L of each 5  $\mu$ M primer (forward and reverse).


### Safety information

Use different pair of indexed primers for each sample.



- 12 Run the PCR program: 30 s denaturation at 98°C; 12 cycles of denaturation at 98 °C for 10 s, combined annealing and extension at 72°C for 30 s (shuttle PCR); the final extension at 72°C for 5 min.
- 13 Check profiles on a gel.

#### Expected result



## Pooling

- 14 Pool all the samples together.
- 15 Perform AMPure purification with ratio 1x. Elute in 100 µl.

#### Protocol



NAME

**AMPure purification protocol**

CREATED BY

Tomasz Suchan

**PREVIEW**

- 15.1 Remove beads from the fridge, equilibrate to the room temperature, mix thoroughly but do not vortex.
- 15.2 Add desired ratio of AMPure beads to the DNA to be purified and mix well by pipetting.
- 15.3 Incubate 5 minutes.
- 15.4 Place on the magnetic rack.
- 15.5 Let it stand for 5 minutes on the rack, aspirate and discard supernatant.
- 15.6 Add 200  $\mu$ l of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.
- 15.7 Repeat the wash: add 200  $\mu$ l of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.
- 15.8 Wait until the pellet dries out completely. DO NOT remove the plate/tubes from the magnetic rack, as the dried AMPure might pop out from the tubes!
- 15.9 Add desired volume of Tris 10 mM or water, wetting the dried AMPure pellets (add 1  $\mu$ l to the final volume to avoid pipetting out the beads)
- 15.10 Remove from the magnetic rack.
- 15.11 Resuspend by pipetting or vortexing.
- 15.12 Incubate 10 minutes, incubating in 37°C can improve DNA yield.
- 15.13 Place on the magnetic rack.



- 15.14 Let it stand for 5 minutes, pippete out and save supernatant. The eluted DNA is in the supernatant, do not discard it!
- 16 Check the concentration using Qubit and the profile using Tapestation/Fragment Analyzer. Calculate molarity from that and proceed to the sequencing. Add 15% PhiX to the sequencing run.