



Oct 11, 2022

## Post-Surgical Dissection of Fallopian Tubes

DOI

[dx.doi.org/10.17504/protocols.io.14egn75eqv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn75eqv5d/v1)

Stephen Fisher<sup>1</sup>, Marielena Grijalva<sup>1</sup>, Rong Guo<sup>1</sup>, Sarah A Johnston<sup>1</sup>, Hieu Nguyen<sup>1</sup>, John Renz<sup>2</sup>, Jean G Rosario<sup>1</sup>, Steven Rudich<sup>2</sup>, Brian Gregory<sup>1</sup>, Junhyong Kim<sup>1</sup>, Kate O'Neill<sup>1</sup>

<sup>1</sup>University of Pennsylvania; <sup>2</sup>Gift of Life Donor Program



Stephen Fisher

University of Pennsylvania

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn75eqv5d/v1>

**Protocol Citation:** Stephen Fisher, Marielena Grijalva, Rong Guo, Sarah A Johnston, Hieu Nguyen, John Renz, Jean G Rosario, Steven Rudich, Brian Gregory, Junhyong Kim, Kate O'Neill 2022. Post-Surgical Dissection of Fallopian Tubes. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn75eqv5d/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 05, 2022

**Last Modified:** October 11, 2022

**Protocol Integer ID:** 68275

**Keywords:** dissection of the fallopian tube, fallopian tube, fallopian tube the muscle layer, dissection, dissection pin, tube, preparation for 10x visium, muscle layer, 10x visium, 10x multiomic, pathology review, pathology, muscular layer, inner mucosa

**Funders Acknowledgements:**







NIH

Grant ID: U54HD104392

## Abstract

This protocol describes dissection of the Fallopian tubes in preparation for 10X Visium, 10X Multiomics, pathology review, and biobanking. The Fallopian tubes have an inner mucosa surrounded by a muscular layer that is coated by serosa. When sectioning the Fallopian tube the muscle layer may retract or protrude and distort orientation. This can be addressed by using a stabilizing hand or dissection pins.

## Materials

- Ice
- Clean cutting surface
- Disposable scalpels
- Kimwipes
- Disposble rulers (with centimeter measurements)
- Ice bucket
- Biohazard bags and container
-  Nuclease-free Water
-  DMEM/F-12 **Thermo Fisher Catalog #11320033**
-  RPMI 1640 Medium **Thermo Fisher Catalog #11875085**
-  MACS® Tissue Storage Solution **Miltenyi Biotec Catalog #130-100-008**
-  DPBS (10X), no calcium, no magnesium **Thermo Fisher Scientific Catalog #14200075** - dilute to 1X with nuclease-free water
-  Marking dye **Cancer Diagnostics Catalog #03000P**

### Equipment

VWR® Pour-Boat Weighing Dishes

NAME

flat-bottomed dish

TYPE

VWR International

BRAND

10803-166


SKU

## Troubleshooting

### Before start

Dilute phosphate buffered saline (Gibco; 14200-075) to 1X with nuclease-free water.

- 1 After removing uterus and Fallopian tubes from ice, dry each Fallopian tube and use marking dye to label the superior edge (12 o'clock), inferior edge (6 o'clock), anterior surface (3 o'clock) and posterior surface (9 o'clock).
- 2 Measure the Fallopian tube from cornua to fimbriated end.



**Measurement of Fallopian Tube.** Left Fallopian tube is displayed from isthmus (left) to fimbriated end (right). Black line made with marking ink at 3 o'clock and red line made at 6 o'clock where the Fallopian tube attaches to the mesosalpinx.

- 3 Using a disposable scalpel, divide one Fallopian tube at the cornua and excise mesosalpinx.
- 4 Weigh Fallopian tube in a disposable weigh boat on ice.
- 5 Place on ice the uterus with the remaining contralateral Fallopian tube.
- 6 Remove the first Fallopian tube from ice and place on a cutting surface.
- 7 Divide Fallopian tube into four sections: isthmus, ampulla, infundibulum, and fimbriae.
- 8 Working one section at a time, slice tissue into four sections.
- 9 Using disposable weigh boats, weigh each tissue piece.
- 10 Tissue can be processed with protocols **OCT-Embedded Tissue Preparation**, **Tissue Fixation Preparation**, or **Snap-Frozen Tissue Preparation** depending on the desired



downstream processing.

- 11 Remove the uterus with the remaining contralateral Fallopian tube from ice and repeat above steps.