

Jul 10, 2024

Postnatal astrocyte labeling by electroporation (PALE)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8j44dg2w/v1>

Shiyi Wang¹

¹Duke University



Shiyi Wang

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8j44dg2w/v1>

Document Citation: Shiyi Wang 2024. Postnatal astrocyte labeling by electroporation (PALE). protocols.io
<https://dx.doi.org/10.17504/protocols.io.5jyl8j44dg2w/v1>

License: This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: May 23, 2023

Last Modified: August 13, 2024



Document Integer ID: 82328

Keywords: ASAPCRN, postnatal astrocyte labeling by electroporation, postnatal astrocyte labeling, astrocyte labeling, electroporation, postnatal, pale

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020607

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

Postnatal astrocyte labeling by electroporation (PALE)

1. Late P0/early P1 mice were sedated by hypothermia until anesthetized, and 1 μ l of plasmid DNA mixed with Fast Green Dye was injected into the lateral ventricle of one hemisphere using a pulled glass pipette (Drummond).
2. For shRNA knockdown experiments in wild-type CD1 mice, the 1 μ l of DNA contained 1 μ g of pGLAST-PBase and 1 μ g of pPB-shRNA-mCherryCAAX was injected.
3. To label astrocytes in WT and LRRK2 G2019S^{ki/ki} mice, the 1 μ l of DNA contained 1 μ g of pGLAST-PBase and 1 μ g of pPB-mCherry-CAAX was injected per mouse.
4. For PALE-mediated overexpression of phospho-mimetic EZRIN in shRNA knockdown experiments, 0.5 μ g pGLAST-PBase, 0.5 μ g pPB-shRNA-mCherryCAAX, and 1 μ g pZac2.1-GfaABC1D-Ezrin T567D-BioID2-HA were injected in a total volume of 1 μ l.
5. For phospho-dead EZRIN overexpression in WT and LRRK2 G2019S^{ki/ki} mice, 0.5 μ g pGLAST-PBase, 0.5 μ g of pZac2.1-gfaABC1D-mCherry-CAAX and 1 μ g of pZac2.1-GfaABC1D-Ezrin T567A-BioID2-HA were injected in a total volume of 1 μ l.
6. Following DNA injection, electrodes were oriented with the positive terminal above the frontal cortex and the negative terminal below the chin, and 5 discrete 50 ms pulses of 100 V spaced 950 ms apart were applied.
7. Pups were recovered on a heating pad, returned to their home cage, and monitored until collection at P21.
8. The number of mice used for each experiment is indicated in the figure legends.
9. All animals appeared healthy at the time of collection. Brain sections were examined for the presence of electroporated cells before staining.