



May 14, 2018

## Pour plating protocol for *Emiliana huxleyi* (single colonies)

DOI

[dx.doi.org/10.17504/protocols.io.p4ydqwx](https://dx.doi.org/10.17504/protocols.io.p4ydqwx)

Inbal Nussbaum<sup>1</sup>, Daniella Schatz<sup>1</sup>

<sup>1</sup>Weizmann Institute of Science

Protist Research to Opti...



**Daniella Schatz**

Weizmann Institute of Science

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.p4ydqwx](https://dx.doi.org/10.17504/protocols.io.p4ydqwx)

**Protocol Citation:** Inbal Nussbaum, Daniella Schatz 2018. Pour plating protocol for *Emiliana huxleyi* (single colonies). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.p4ydqwx>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 14, 2018

**Last Modified:** May 14, 2018

**Protocol Integer ID:** 12152

### Abstract

A protocol for generating single colonies of *Emiliana huxleyi* embedded in soft-agarose.

- 1 Autoclave 75 mL DDW with 0.75 g agarose
- 2 For 300 mL (= 7 plates): heat 225 mL FSW in a 50 °C water bath.
- 3 Cool autoclaved agar to 50 °C, mix with heated FSW, add F/2 stock and antibiotics if necessary. Swirl/mix well.
- 4 For each sample, pour 40 mL of FSW-agar mixture to a 50 mL conical tube.
- 5 Move to 32 °C bath.
- 6 Prepare the cells at desired cell concentration. For each sample prepare 0.5ml cells.
- 7 Lay out plates.
- 8 When you are ready, remove FSW-agar mix from 32 °C bath and allow to cool 1-2 degreesCritical: the mixture will start to solidify at 27-28 °C.
- 9 Add cell suspension to each tube, mix gently, then pour into plates.
- 10 Let solidify to a loose consistency
- 11 Transfer to growth room. Single colonies should appear within 10 days (without selection).