



Oct 10, 2025

Pre-MALDI Autofluorescence Microscopy of Lung Tissues

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g75z23lwz/v1>

Brittney Gorman¹, Christopher Anderton¹

¹Pacific Northwest National Laboratory

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Brittney Gorman

Pacific Northwest National Lab

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g75z23lwz/v1>

Protocol Citation: Brittney Gorman, Christopher Anderton 2025. Pre-MALDI Autofluorescence Microscopy of Lung Tissues. [protocols.io https://dx.doi.org/10.17504/protocols.io.q26g75z23lwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g75z23lwz/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 25, 2025

Last Modified: October 10, 2025

Protocol Integer ID: 124931

Keywords: Autofluorescence, Zeiss, MALDI, optical imaging, HuBMAP, Lung, PNNL, autofluorescence microscopy images of lung tissue section, autofluorescence microscopy image, brightfield image of each tissue section, lung tissue, lung tissue section, visualization of tissue morphology, imaging, tissue morphology, additional imaging modality, brightfield image, tissue section, tissue, image

Funders Acknowledgements:

Hubmap

Grant ID: U54HL165443

Abstract

Obtain autofluorescence microscopy images of lung tissue sections. This will produce an RGB plus brightfield image of each tissue section that enables visualization of tissue morphology, and these images can be co-registered with additional imaging modalities.

Guidelines

Carefully handle all slides with gloves

Materials

Zeiss Axio Zoom Microscope

Autofluorescence Microscopy Imaging

- 1 If tissue-mounted slides are frozen, allow them to equilibrate to room temperature (~20 °C) before opening the vacuum-sealed bag (~30 min); otherwise, proceed directly to step 2.
- 2 Place the microscope slide(s) onto the stage of the Zeiss Axio Zoom Microscope
- 3 Open Zeiss Zen software
- 4 Set up the method for acquiring autofluorescence (RGB-channels) and transmitted light images

Note

Magnification: 100 x
Objective: PalnNeoFluar Z 2.3x

Blue: Chroma DAPI (ex. 353 nm, em. 465 nm)

Light Source Intensity: 100%

Exposure time: 350 ms

- Captures NAD(P)H excitation present in cells and elastin in the extracellular matrix

Green: Chroma GFP (ex. 488 nm, em. 509 nm)

Light Source Intensity: 100%

Exposure time: 500 ms

- Captures Flavins in metabolically active zones (increased in the vasculature)

Red: Chroma DSRED (ex. 545 nm, em. 572 nm)

Light Source Intensity: 100%

Exposure time: 1300 ms

- Captures lipopigments often found in muscle

Transmitted Light:

Light Source Intensity: 100%

Exposure time: 1450 ms

- General morphological assessment of the tissue

- 5 Define imaging regions to include each individual tissue section within a single region
- 6 Perform manual focus adjustment at ~6 'support points' per tissue section



- 6.1 Switch to the green channel, press 'Live Scan', and perform focus adjustments at each support point
- 7 Press 'Start Scan' to begin image acquisition. The microscope will automatically scan each tissue section and interpolate the focus between support points
- 8 The method will acquire a tiled brightfield and (auto)fluorescence image.
- 9 Export each image region as a Tiff and PNG using Zen
(The .czi files can be visualized using Zen or QuPath software)
- 9.1 Under the processing tab use the 'Split Scenes' method to create individual .czi files for each region
- 9.2 Under the processing tab, use the 'Image export' method to export as a Tiff and PNG
Resize: 50%