

May 20, 2019

Preparation of a cell-free expression system from Escherichia coli

DOI

dx.doi.org/10.17504/protocols.io.2xkgfkw

Renuka Kudva¹

¹Stockholm University



Renuka Kudva

Stockholm University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.2xkgfkw

Protocol Citation: Renuka Kudva 2019. Preparation of a cell-free expression system from Escherichia coli. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.2xkgfkw>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol is based on a published paper, it has been tried and tested under different conditions. All these have been mentioned. Please feel free to get in touch with any comments.

Created: May 16, 2019

Last Modified: May 20, 2019

Protocol Integer ID: 23244

Keywords: Cell-free protein production, E.coli, S30, S135, in vitro transcription and translation, membrane-free E.coli cell-extract,



Abstract

This protocol describes the procedure make an S30 derived E.coli cell lysate that lacks membrane components. Modified from a protocol published in <https://doi.org/10.1091/mbc.e11-07-0590>
The procedure takes three days in total, including preparation time.

Guidelines

Start up ultra-centrifuges 30 min prior to use to stabilise vacuum. Work on ice. Plot growth curves for the strain used in the medium used for preparation (see materials.)



Materials

MATERIALS

- ✕ Magnesium acetate
- ✕ Sodium Hydroxide **Fisher Scientific Catalog #BP359500**
- ✕ L-aminoacids **Merck MilliporeSigma (Sigma-Aldrich) Catalog #LAA21-1KT**
- ✕ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**
- ✕ NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**
- ✕ Tris acetate **Bio Basic Inc. Catalog #TD0101.SIZE.500g**
- ✕ BD Bacto™ Yeast Extract **Becton Dickinson (BD) Catalog #212750**
- ✕ Sucrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S7903**
- ✕ Potassium acetate **Merck Millipore (EMD Millipore) Catalog #1.04820.1000**
- ✕ Tryptone (pancreatic digest of casein) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T9410**
- ✕ Phospho(enol)pyruvic acid mono potassium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #860077**
- ✕ Adenosine 5'-triphosphate (ATP) disodium salt hydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A2383**
- ✕ Pyruvate kinase from rabbit muscle **Catalog #83330**

Cells can be cultured in LB medium (buffered with phosphates), TB or S30 medium. Prepare 4 litres i.e. 1 litre each in 4 5 litre erlenmeyer flasks. Autoclave.

S30 medium

9 g/l bacto-tryptone

0.8 g/l yeast extract

5.6 g/l NaCl

1 mM NaOH

Prepare a mix of the 20 amino acids 1 mM each. Aliquot and store at -20°C.

Prepare stock solutions of the acetates: 1M Tris Acetate at pH 7.5, 1M Magnesium acetate at pH 7.5, 4 M Potassium acetate at pH 7.5

Prepare stock solutions: 0.25 M ATP pH neutralised with KOH, 0.2 M PEP. Make aliquots and store at -20°C.

Strains tested: E.coli KC6, MRE600, MC4100.

If using a rich medium, do a growth curve for the strain being used under those growth conditions, and harvest at early exponential phase.

S30 buffer (pH 7.5)



10 mM Tris-acetate
14 mM magnesium acetate
60 mM potassium acetate
1 mM DTT

Sucrose cushion (pH 7.5)
10 mM Tris acetate
14 mM Magnesium acetate
60 mM potassium acetate
1 mM DTT
1.44 M sucrose

Safety warnings

- ⚠ Use ultracentrifuges carefully, ensure to fill tubes up to requisite volumes. Handle chemicals according to local safety requirements. Wear eye-goggles and a mask while weighing yeast extract and peptone. Use fume hoods where indicated.
Handle sharps with care. Dispose off responsibly.
Handle liquid nitrogen with cryo-protective gear including eye goggles.
Handle autoclaves carefully.

Before start

Wash all flasks and bottles with milliQ water. Prepare all buffers in milliQ water.

