

Aug 21, 2018

## Preparation of Cas9 Protein

DOI

[dx.doi.org/10.17504/protocols.io.rczd2x6](https://dx.doi.org/10.17504/protocols.io.rczd2x6)



Ashwin Banfal<sup>1</sup>, Brett T. Staahl<sup>2</sup>, Steve Lin<sup>2</sup>

<sup>1</sup>University of California, San Francisco; <sup>2</sup>University of California, Berkeley

The Center for Genome ...



Meredith Triplet

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.rczd2x6](https://dx.doi.org/10.17504/protocols.io.rczd2x6)

**Protocol Citation:** Ashwin Banfal, Brett T. Staahl, Steve Lin 2018. Preparation of Cas9 Protein. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rczd2x6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 28, 2018

**Last Modified:** August 21, 2018

**Protocol Integer ID:** 13433

## Abstract

This is a protocol from the **Doudna Lab** for preparation of Cas9 protein.

## Attachments



Cas9\_purification\_Br...

135KB

## Guidelines

The protocol workflow is as follows:

### 1. Expression of Cas9 in Rosetta2 E. coli

- Transformation ([Steps 1-5](#))
- Starter culture ([Steps 6-8](#))
- Expression ([Steps 9-21](#))

### 2. Purification of Cas9


- Sonication ([Steps 22-29](#))
- HisTrap ([Steps 30-40](#))
- TEV ([Steps 41-45](#))
- Heparin ([Steps 46-53](#))
- BPTrap/Superdex 200 ([Steps 54-65](#))

## Materials

### MATERIALS

 Please see before starting in Guidelines section for materials.

## Safety warnings

 Please refer to the SDS (Safety Data Sheet) for safety warnings and hazard information.

## Before start

### Materials

<b>TY media</b>	<b>1 L</b>
Tryptone	10 g
Yeast extract	6 g
NaCl	10 g
Dissolve in 400 mL Rinse flasks thoroughly with MQ water Add 600 mL MQ water, media Autoclave (45', no dry time)	

### Buffers

<b>Lysis</b>	<b>1L, pH7.5</b>
20mM HEPES	20 ml of 1M
1M KCl	74.55 g
10% glycerol	100 ml
5mM TCEP	1.43 g
10mM imidazole	680 mg

<b>Elution</b>	<b>500mL, pH7.5</b>
20mM HEPES	10 ml of 1M
0.1M KCl	50 ml of 1M
10% glycerol	50 ml
5mM TCEP	717 mg
300mM imidazole	10.2 g




<b>Ion Exchange buffer A</b>	<b>1L, pH7.5 (low salt)</b>
20mM HEPES	20 ml of 1M
300mM KCl	300 ml of 1M
10% glycerol	100 ml
1mM TCEP	287 mg

<b>Ion Exchange buffer B</b>	<b>1L, pH7.5 (high salt)</b>
20mM HEPES	20 ml of 1M
1M KCl	74.55 g
10% glycerol	100 ml
1mM TCEP	287 ml

<b>Gel Filtration</b>	<b>1L, pH7.5</b>
20mM HEPES	20 ml of 1M
150mM KCl	150 ml of 1M
10% glycerol	100 ml
1mM TCEP	287 mg


## Expression of Cas9 in Rosetta2 E. coli - Transformation


- 1 Mix 1  $\mu$ L of 10 ng/ $\mu$ L DNA + 10  $\mu$ L Rosetta 2 cells (pRAR-rare codons).


 1  $\mu$ L 10 ng/ $\mu$ L DNA

 10  $\mu$ L Rosetta 2 cells


- 2 Heat shock: 10' ice, 45" 42°C, 5' ice.

 00:10:00 ice


 00:00:45 42°C

 00:05:00 ice

- 3 Recover: Add 300  $\mu$ L LB, 50' 37°C on shaker.

 300  $\mu$ L LB

 37 °C Shaker

 00:50:00 Shaker


- 4 Spin down (6000 xg, 1'), and remove 150  $\mu$ L LB.

 00:01:00 Spin down

- 5 Plate 100  $\mu$ L, 30  $\mu$ L on 2 AMP plates.

## Expression of Cas9 in Rosetta2 E. coli - Starter culture

- 6 Pick 6 colonies to add to 20 mL LB + Amp (Cf=100 ug/mL) starter culture.

 20 mL LB + Amp (Cf=100 ug/mL) starter culture

- 7 Grow overnight.

- 8 Next day: Store XL1 blue 20% glycerol stocks in -80°C.

 -80 °C Storage

## Expression of Cas9 in Rosetta2 E. coli - Expression

- 9 Add AMP to the flasks (Cf=100 ug/mL).




10 Add 20 ml starter culture in the morning,


 20 mL starter culture

11 Grow to OD600=0.6 (150 rpm, ~4 hrs, 37°C).

 37 °C Growing

 04:00:00 Growing

12 Set shaker temp to 16°C.


 16 °C Shaker temperature

13 \*Save 1 ml for SDS-PAGE: **uninduced**.


14 Add 0.5 mM IPTG.


15 Grow overnight (150 rpm, at least 16 hrs, 16°C).

 16 °C Growing

 16:00:00 Growing

16 Harvest cells before lunch: Spin down (4000 rpm, 20', 4°C).


 4 °C Spin down

 00:20:00 Spin down


#### Note

**\*Keep the cells cold from here.**

17 Completely resuspend pellets in 20 ml Lysis Buffer + 1mM PMSF per bottle.

 20 mL Lysis Buffer + 1mM PMSF

18 Use additional 20 ml Lysis Buffer to wash all bottles.

 20 mL Lysis Buffer

19 Transfer to 50 ml-tubes.

20 \*Save 20 µl for SDS-PAGE: **cell pellet**



21 Freeze the cells at -80°C.

🧊 -80 °C Freezing cells

## Purification of Cas9 - Sonication

22 Thaw cells in water.

23 Push beaker with a stir bar to the bottom of ice bucket.

24 Add ~50 ml Lysis Buffer and 1 protease inhibitor tablet/50 ml cell. Dissolve the tablets.

🧴 50 mL Lysis Buffer

25 Transfer the cells to the beaker.

26 Sonicate while stirring: 8' total, 10" on, 10" off, level 4. Lysate should turn fluid and thin. Repeat 8 min cycle, if necessary.

⌚ 00:08:00 total

⌚ 00:00:10 on/off

27 Transfer to Oak Ridge tubes (30 ml/tube). \*Balance within +/- 0.01g.

28 Pellet cell debris (18,000 rpm, 30', 4°C). Prepare HisTrap purification.

### HisTrap

CV=5mL, FR=2mL/min

Buffer A: Lysis buffer, Buffer B: Elution buffer

🧊 4 °C Pellet cell debris

⌚ 00:30:00 Pellet cell debris

29 \*Save 20 µl supernatant for SDS-PAGE: **cell lysate**



## Purification of Cas9 - HisTrap

30 Filter sample.

### Note

#### **HisTrap**

CV=5ml, FR=2ml/min

Buffer A: Lysis buffer, Buffer B: Elution buffer

31 Set-up Akta:

1. Turn on UV lamp.
2. Wash Akta with water: pump wash, and flow 50% A, 50% B.
3. Pump wash.
4. Flow 50% A, 50% B.
5. Flow 100% A, until steady baseline. Set baseline to zero.

32 Apply sample with pump, in the cold room:

1. Fill tubing with water, connect column.
2. Wash: 5CV water.
3. Equilibrate: 5CV lysis buffer.
4. Apply sample to the column. Collect flow-through: Ni-sample loading.  
Re-apply flow-through to the column, if necessary.

33 Purification:

In the cold room:

First wash: 15CV lysis buffer. Collect flow-through: Ni-wash1.

Collect flow-through: **Ni-wash1**. \*Save 20µl for SDS-PAGE.

34 Wash2: 3CV lysis buffer + 0.1% Triton X.

Collect flow-through: **Ni-wash2**.








\*Save 20 µl for SDS-PAGE.


- 35 Purification using Akta:  
First, attach column to Akta (FR=0.9 ml/min).
- 36 Wash: 5CV lysis buffer, until steady UV baseline.
- 37 Elute: 25CV elution buffer, 100% gradient.
- 38 Collect fractions (Fraction size=2.5 ml).
- 39 Clean-up:
1. Wash: 5CV water
  2. Store: 5CV 20% ethanol
  3. Detach column.
  4. Wash Akta with water.
- 40 Pool appropriate fractions. Rinse collection tubes with elution buffer.  
\*Save 10 µl for SDS-PAGE: **Ni-elution**

## Purification of Cas9 - TEV

- 41 Concentrate with 30K Amicon to ~3 ml volume. Collect flow-through: Ni-Amicon FT.
- 42 Add ~1.25 mg TEV/3 ml sample, and mix with pipet.  
 1.25 mg TEV
- 43 Incubate on shaker, at 4°C overnight.  
 4 °C Incubation on shaker  
 16:00:00 overnight incubation
- 44 Next day: \*Save 10 µl for SDS-PAGE: **postTEV**



45 Add 25 ml of buffer A to sample.

 25 mL Buffer A

## Purification of Cas9 - Heparin

46 **Heparin**

CV=5 mL, FR=1 mL/min

Buffer A: IEX low salt buffer, Buffer B: IEX high salt buffer

47 Set-up Akta (FR=2 ml/min).

48 Apply sample with pump, in the cold room:

1. Equilibrate: 10CV buffer A.

2. Apply sample to the column. Collect flow-through: Hep-sample loading.

49 Purification:

1. Attach column to Akta (FR=0.9 ml/min).

2. Wash: 15CV buffer A, until steady baseline.

3. Collect flow-through: **Hep-wash**. \*Save 20 µl for SDS-PAGE.

4. Elute: 20CV buffer B, 100% gradient.

5. Collect fractions (Fraction size=2 ml).

50 Clean-up (FR=2ml/min):


1. Wash: buffer B, until steady conductivity baseline for over 30 ml.

2. Store: 5CV 20% ethanol.

3. Detach column.

4. Wash Akta with water.



- 51 Pool appropriate fractions. Rinse collection tubes with buffer B.  
\*Save 10  $\mu$ l for SDS-PAGE: **Hep-elution**
- 52 Concentrate sample to ~1 ml. Collect flow-through: Hep-Amicon FT.
- 53 Leave concentrator on ice, at 4°C overnight.  
 4 °C Concentrator on ice




## Purification of Cas9 - MBPTrap/Superdex 200

### 54 **MBPTrap/Superdex 200**

CV=120 mL, FR=0.5-1 mL/min  
Buffer A: Gel filtration buffer

#### **Previous day:**

1. Set-up Akta (FR=2 mL/min):
  2. Prepare columns:
    - a. MBPTrap wash: 5CV water, 5CV gel filtration buffer (FR=2 mL/min), in the cold room.
    - b. Attach MBPTrap and Superdex200 to Akta.
    - c. Equilibrate overnight: 2CV gel filtration buffer (FR=0.3 mL/min).
- \*Set end timer volume=240 mL, max pressure alarm=0.25 MPa.

- 55 Concentrate sample to ~500  $\mu$ l. Rinse concentrator with 500  $\mu$ l gel filtration buffer.  
 500  $\mu$ L Gel filtration buffer
- 56 Spin down sample (11000 rpm, 4', 4°C). Transfer supernatant to a new 1.5 ml-tube.  
 4 °C Spin down  
 00:04:00 Spin down
- 57 To inject sample, first, clean injection port, sample loop with water, gel filtration buffer.
- 58 Now, draw up sample with 1 ml-syringe, and replace the syringe on the injection port.



59 [LOAD] Inject sample to the loop with syringe.

60 Purification

1. [INJECT] Apply sample to the column (~3 ml).

2. [LOAD] Flow ~1CV gel filtration buffer. Cas9 elutes at 65-85 ml.

3. Collect fractions (Fraction size=2 ml).

61 Clean-up (FR=2 ml/min):

1. Detach MBPTrap.

2. Wash MBPTrap: 10CV 10mM maltose + gel filtration buffer, in the cold room.

3. Store MBPTrap: 5CV 20% ethanol.

4. Wash Superdex200: 2CV water (Long-term storage in 20% ethanol).

5. Detach Superdex200.

6. Wash Akta with water.

62 Pool appropriate fractions. Rinse collection tubes. \*Save 5 µl for SDS-PAGE: **Cas9**

63 Concentrate to desired concentration (~40 uM).

*Cas9	MW=158,441.4 g/mol Molar absorption coefficient=120,575 L/mol·cm
-------	---

64 Aliquot 20 µl each into 1.5 ml-tubes.

65 Snap-freeze in LN2, and store at -80°C.

 -80 °C Storage