

Jul 16, 2018

## Preparation of chemical competent E. coli cells (w/ calcium chloride)

DOI

[dx.doi.org/10.17504/protocols.io.qe4dtgw](https://dx.doi.org/10.17504/protocols.io.qe4dtgw)

Dennis DD Dienst<sup>1</sup>, Hanna Eriksson

<sup>1</sup>Department of Chemistry - Microbial Chemistry, Ångström Laboratory, Uppsala, Sweden

CyanoWorld



Hanna Eriksson

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.qe4dtgw](https://dx.doi.org/10.17504/protocols.io.qe4dtgw)

**Protocol Citation:** Dennis DD Dienst, Hanna Eriksson 2018. Preparation of chemical competent E. coli cells (w/ calcium chloride). [protocols.io](https://dx.doi.org/10.17504/protocols.io.qe4dtgw) <https://dx.doi.org/10.17504/protocols.io.qe4dtgw>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 24, 2018

**Last Modified:** July 16, 2018

**Protocol Integer ID:** 12476



## Materials

### MATERIALS

⊗ Glycerol **Catalog #G5516**

⊗ Calcium Chloride Dihydrate **Fisher Scientific Catalog #C79**

⊗ Manganese(II) chloride tetrahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3634**

⊗ 1.5 mL Eppendorf tubes

⊗ Magnesium chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich)**




## Preparation of buffer CCMB 80

- 1
  - 10 mL of 1M K-Acetate stock
  - 11.8 g  $\text{CaCl}_2 \cdot \text{H}_2\text{O}$
  - 2 g  $\text{MgCl}_2 \cdot 6 \text{H}_2\text{O}$
  - 4 g  $\text{MnCl}_2 \cdot 4 \text{H}_2\text{O}$
  - 100 mL glycerol
  - adjust pH to 6.4 w/ 0.1M HCl
  - add ddH<sub>2</sub>O to 1 L
- filter buffer through syringe filter (e.g. Acrodisc Syringe Filter, 0.2  $\mu\text{M}$  Supor Membrane, Pall Life Science)  
→ work in the Clean Bench

### Safety information

HCl is corrosive! Wear goggles, gloves and don't inhale! Work in the hood!  
[https://en.wikipedia.org/wiki/Hydrochloric\\_acid](https://en.wikipedia.org/wiki/Hydrochloric_acid)

 1000 mL Buffer CCMB 80

## Autoclave material

- 2
  - 255 mL LB medium
  - >100 × 1.5 mL Eppendorf tubes (e.g. in a big flask or beaker)
  - 1 mL Erlenmeyer flask (closed w/ aluminium foil)
  - 1 x > 10 mL test tube (if not using disposables)
- store sterile Eppendorf tubes at -20°C

## Overnight (o/n) culture

- 3
  - Inoculate 5 mL of LB medium w/ seed cells (e.g. DH5 $\alpha$ ) in a test tube
  - Incubate at 37 °C and >200 rpm

## Overday (o/d) culture

- 4
  - Inoculate 250 mL fresh LB medium w/ 5 mL o/n culture (1:50) in 1L Erlenmeyer flask
  - Incubate at 37 °C and >200 rpm until OD<sub>600</sub> ~ 0.3



#### Note

Work in the Clean Bench!

## CaCl<sub>2</sub> treatment

- 5
  - Cool down centrifuge including inlets.
  - transfer o/d cultures to five sterile 50 mL tubes (Falcon)
  - Centrifuge at 3000 \*rcf and 4 °C for 10 min
  - gently resuspend each pellet in 16 mL **ice cold** CCMB
  - incubate on ice for 20 min
  - Centrifuge o/d culture at 3000 \*rcf and 4 °C for 10 min
  - gently resuspend each pellet in 2 mL **ice cold** CCMB
  - incubate on ice for 20 min

🧊 4 °C keep cells on ice

## Prepare aliquots

- 6
  - transfer each 100 µL of cell suspension in 1.5 mL Eppendorf tube
  - store at -80°C

🧊 4 °C keep tubes on ice

#### Note

**Option:** Quick-freeze aliquots instantly in liquid nitrogen

#### Note

This yields slightly more than 100 1,5 ml Eppendorf tubes. Tip: mark the tubes beforehand. This makes the tubes easier to find in the ice when transforming the cells afterwards. To reduce workload, mark the whole freezer box with the type of competent cells and only make a line or dot on the Eppendorf tubes, preferably with a thick marker.