

Sep 28, 2023

Preparation of Competent Cells (10 β E. coli Strain)

DOI

dx.doi.org/10.17504/protocols.io.q26g7p5w9gwz/v1

NUS iGEM¹

¹National University of Singapore



NUS iGEM

National University of Singapore

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7p5w9gwz/v1>

Protocol Citation: NUS iGEM 2023. Preparation of Competent Cells (10 β E. coli Strain). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7p5w9gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 27, 2023

Last Modified: September 28, 2023




Protocol Integer ID: 88455

Keywords: Competent Cells, Transformation, preparation of competent cell, competent cell, cell, singapore igem team, preparation


Abstract

2023 NUS-Singapore iGEM team followed this protocol to make competent cells that would be used for transformation.

Guidelines


This protocol demonstrates the process of making 5 tubes of competent cells from a  5 mL cell culture. Generally, within our protocol, every  5 mL of cultured cells can be transformed into 5 tubes of competent cells, with each tube containing  50 μ L of cells.

Materials

1.  NEB 10-beta Competent E.coli (High Efficiency) - 6×0.2 ml **New England Biolabs Catalog #C3019I**
2. LB Media
3. MgCl₂ Solution
4. CaCl₂ Solution
5. 100% Glycerol Solution







Troubleshooting

Safety warnings





-  Proper lab PPE must be worn at all times.
- Thermal gloves shall be worn when handling cell stock from the -80°C fridge.



Cell Culture from Cell Stock


- 1 Prepare a Falcon tube with  5 mL of LB media.
- 2 Prepare an ice box.
- 3 Take out a tube of  Sample cell stock from the  -80 °C fridge and put it into the ice box.
- 4 In the biosafety cabinet (BSC), use an inoculation loop to inoculate some competent cells into the Falcon tube with  5 mL of LB media.
- 5 Incubate the cells in an incubator at  37 °C for  Overnight

Refresh Cell Culture

- 6 Prepare a new Falcon tube and add  10 mL of LB media into the tube.
- 7 Add  100 µL of the pre-cultured cells into this new Falcon tube to refresh the cells.
- 8 Incubate the cells at  37 °C for  02:00:00 .

2h

Pre-Cell Washing

- 9 Pre-cool the centrifuge machine to  4 °C .
- 10 Take out the Falcon tube from the incubator, ensuring that the optical density (OD) of the cultured cells is 0.6OD to 0.8OD.






- 11 Place the Falcon tube in ice for 00:30:00 . 30m
- 12 Prepare 10 mL of 0.1M MgCl₂ and 10 mL of 0.1M CaCl₂, put both tubes in ice for 00:30:00 . 30m

Cell Washing

- 13 Centrifuge the Falcon tube with cultured cells in the pre-cooled centrifuge machine at 5000 rpm for 00:05:00 . The temperature in the centrifuge machine must be kept at 4 °C the whole time in the "Cell Washing" section. 5m
- 14 Discard the supernatant and keep the cell pellet.
- 15 Add a small amount of MgCl₂ solution prepared in the earlier step into the Falcon tube to resuspend the cell pellet. Then, pour the rest of the MgCl₂ solution into the Falcon tube.
- 16 Centrifuge the Falcon tube again at 5000 rpm for 00:05:00 . 5m
- 17 Discard the supernatant and keep the cell pellet.
- 18 Add a small amount of CaCl₂ solution prepared in the earlier step into the Falcon tube to resuspend the cell pellet. Then, pour the rest of the CaCl₂ solution into the Falcon tube.
- 19 Place the Falcon tube in ice for 01:00:00 . 1h
- 20 Centrifuge the Falcon tube with the cells at 5000 rpm for 00:05:00 . 5m
- 21 Discard the supernatant and keep the cell pellet.

Storage



- 22 Prepare a  1 mL mixed solution composed of 20% glycerol and 80% 0.1M CaCl₂ solution, and put it into the ice box to cool its temperature down.
- 23 Add  250 µL of glycerol-CaCl₂ solution into the Falcon tube with the cell pellet and resuspend the cell pellet.
- 24 Split the cells in the Falcon tube into 5 new Eppendorf tubes with each tube containing  50 µL of cells.
- 25 The competent cells can be stored in a -80°C fridge or used immediately.