

Feb 26, 2019

Preparation of feeder-free iPSCs culture

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.x83fryn>

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x83fryn>

Protocol Citation: Celeste Karch, Rita Martinez, Jacob Marsh 2019. Preparation of feeder-free iPSCs culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x83fryn>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2019

Last Modified: February 26, 2019

Protocol Integer ID: 20475

Keywords: free ipscs culture, preparation

Attachments



IPSC CORTICAL

DIFFER...

179KB

Guidelines

This protocol is part of the **IPSC CORTICAL DIFFERENTIATION** collection.

This method should be performed using sterile technique.

Materials

Please refer to the attached full manuscript for required materials.

Safety warnings

⚠ Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



- 1 Coat appropriate sized tissue culture dish with Matrigel. Swirl plate to distribute evenly across surface of dish.
- 2 Allow Matrigel to set for  01:00:00 at room temperature or in a  37 °C humidified chamber.

Note

This will provide enough cells for a full 96 well plate of neural aggregates.

Note

We typically thaw 1 vial of iPSC into 1 well of a 6 well (1mL Matrigel per well). Volumes below are based on 6 well plates.

- 3 Prepare mTeSR1 feeder-free medium by aseptically adding  100 mL of thawed 5x supplement to  400 mL Basal Medium.
- 4 Add  5 mL penicillin/streptomycin if needed. Mix thoroughly and warm to room temperature.
- 5 Thaw iPSCs ($1-2 \times 10^6$ cells/mL) in a  37 °C water bath.
- 6 Add  1 mL of thawed iPSC to  9 mL of DMEM/F12 in a 15mL conical tube (Do not mix). Centrifuge iPSC at 750 rpm for  00:03:00 .
- 7 Carefully aspirate medium from iPSC pellet. Add  2 mL of mTeSR1 supplemented with Rock inhibitor (10 μ M final) - **avoid breaking up clumps of suspended cells.**



- 8 Aspirate Matrigel from tissue culture dish. Add  2 mL of suspended iPSCs. Place plate in  37 °C , 5% CO₂ and 95% humidified chamber, making a "T" motion to ensure even distribution of cells.
- 9 Change medium with  2 mL of fresh mTesR1 medium daily.

Note

As iPSC become more confluent, increase mTesR to 3mL per well.