

Apr 05, 2024

Preparation of Free Floating Coronal Mouse Brain Sections

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8pzk9g2w/v1>

madalynn.erb Erb¹

¹Van Andel Research Institute

Madalynn's Workspace



madalynn.erb Erb

Van Andel Research Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8pzk9g2w/v1>

Protocol Citation: madalynn.erb Erb 2024. Preparation of Free Floating Coronal Mouse Brain Sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8pzk9g2w/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 26, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94228

Keywords: ASAPCRN, floating coronal mouse brain sections protocol, floating mouse brain section, coronal mouse brain sections protocol, mouse brain section, mouse brain, immunostaining, mouse, situ hybridization

Disclaimer

The [protocols.io](https://www.protocols.io) team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

Protocol to prepare free floating mouse brain sections for immunostaining or *in situ* hybridization.



Collect Mouse Brain Tissue

- 1 Deeply anesthetize mice via intraperitoneal injection of 2X Avertin solution
- 2 Perform transcardial perfusion using chilled saline solution
 - 2.1 0.9% NaCl kept on ice
 - 2.2 Use approximately 60mL saline solution per mouse
- 3 Switch from saline solution to chilled PFA
 - 3.1 4% paraformaldehyde in 0.1M phosphate buffer (PB) pH 7.4
 - 3.2 Use approximately 60mL PFA per mouse
- 4 Remove brain immediately after PFA perfusion
- 5 Incubate brain in PFA for 24 hours at 4°C
- 6 Transfer brains to 30% Sucrose / 0.1M Phosphate Buffer (PB) – keep at 4°C for ≥ 24 hours
- 7 Tissue should be completely saturated with sucrose before sectioning
 - 7.1 Brains will sink to the bottom of the vial when saturated

Section Tissue

- 8 Use a Leica SM2010R Microtome
 - 8.1 Blade: Leica 16cm, knife angle set at 0 degree
 - 8.2 Cut thickness: 35µm
- 9 Adjust the microtome platform so that it is level with the blade and lock it into place
- 10 Chill the microtome platform by covering it with crushed dry ice
- 11 Apply 5-10 drops of 30% Sucrose / 0.1M Phosphate Buffer (PB) solution to the chilled microtome platform and wait for it to solidify
- 12 Use the microtome to gently shave the solidified sucrose to make a flat surface
- 13 Use a razor blade to remove any spinal cord from the brain making a flat surface perpendicular to the rostral / caudal axis
- 14 Apply 2-3 drops of sucrose to the existing sucrose platform and quickly position the brain with the olfactory bulbs pointing upwards
 - 14.1 Apply 2-3 more drops of sucrose to the top of the brain to securely freeze it to the microtome platform
- 15 Gently cover the brain in crushed dry ice to freeze the tissue
- 16 Adjust the microtome platform so that the rostral / caudal axis is perpendicular to the blade and the dorsal / ventral axis is level with the blade
- 17 Collect sections in a 24-well plate prefilled with cryoprotectant solution (0.1M PB +30% Sucrose + 30% Ethylene Glycol)



- 18 Seal plate with parafilm and store tissue stored at -20°C

Recipes

- 19 0.2M Phosphate Buffer (PB)

	A	B
	1L	H ₂ O
	22.72 g	Na ₂ HPO ₄
	5.52 g	NaH ₂ PO ₄

- 20 4% PFA

	A	B
	500mL	8% paraformaldehyde
	500mL	0.2M PB

- 21 Saline

	A	B
	1L	H ₂ O
	9 g	NaCl

- 22 30% Sucrose

	300g	Sucrose
	Fill to 1L	0.1M PB

1. Stir and heat to 60°C until dissolved



23 Cryoprotectant Solution

	A	B
	300g	Sucrose
	300mL	Ethylene glycol
	Fill to 1L	0.1M PB

Completely dissolve sucrose in 0.1M PB before adding ethylene glycol. After adding ethylene glycol, fill to 1L with 0.1M PB.