

May 08, 2020

Preparation of Gelatin Slide Coating Solution

DOI

dx.doi.org/10.17504/protocols.io.bfzvjp66

Allen Institute for Brain Science¹

¹Allen Institute

BICCN / BICAN

Allen Institute for Brain S...



Allen Institute

Allen Institute

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bfzvjp66

Protocol Citation: Allen Institute for Brain Science 2020. Preparation of Gelatin Slide Coating Solution. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bfzvjp66>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 05, 2020

Last Modified: March 17, 2021

Protocol Integer ID: 36629

Keywords: RP0194, glass microscope slides



Abstract

This protocol describes the preparation of gelatin solution to be used for coating glass microscope slides.

Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Attachments



RP0194_Gelatin_slide...

250KB

