

Sep 25, 2019

🌐 Preparation of LB Agar

DOI

dx.doi.org/10.17504/protocols.io.7pjhmkn

NUS iGEM¹

¹National University of Singapore



NUS iGEM

National University of Singapore

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7pjhmkn

Protocol Citation: NUS iGEM 2019. Preparation of LB Agar. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7pjhmkn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 25, 2019

Last Modified: September 11, 2023

Protocol Integer ID: 28107

Materials

MATERIALS

 Water refers to sterilized deionized water

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

STEP MATERIALS

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

 Water refers to sterilized deionized water

Protocol materials

 Water refers to sterilized deionized water

 Water refers to sterilized deionized water

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

 Water refers to sterilized deionized water

1 Weigh  32 g of LB agar powder.

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

2 Add the powder into  1 L of deionized water.

 Water refers to sterilized deionized water

3 Autoclave entire bottle of LB media.