

Sep 26, 2019

Version 2

## Preparation of LB Media V.2

DOI

[dx.doi.org/10.17504/protocols.io.7pdhmi6](https://dx.doi.org/10.17504/protocols.io.7pdhmi6)

NUS iGEM<sup>1</sup>

<sup>1</sup>National University of Singapore



**NUS iGEM**

National University of Singapore

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.7pdhmi6>

**Protocol Citation:** NUS iGEM 2019. Preparation of LB Media. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7pdhmi6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 26, 2019

**Last Modified:** September 12, 2023

**Protocol Integer ID:** 28101

**Keywords:** preparation

## Materials

### MATERIALS

⊗ Water refers to sterilized deionized water

⊗ Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo Fisher Catalog #12795027**

### STEP MATERIALS

⊗ Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo Fisher Catalog #12795027**

⊗ Water refers to sterilized deionized water

## Protocol materials

⊗ Water refers to sterilized deionized water

⊗ Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo Fisher Catalog #12795027**

⊗ Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo Fisher Catalog #12795027**


⊗ Water refers to sterilized deionized water


⊗ Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo Fisher Catalog #12795027**


⊗ Water refers to sterilized deionized water


## Troubleshooting



1 Weigh  25 g of Luria Broth Base powder.

 Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo**  
**Fisher Catalog #12795027**

2 Add the powder into  1 L of water.

 Water refers to sterilized deionized water

3 Autoclave entire bottle of LB media.