


May 23, 2022

🌐 Preparation of LRRK1 RCKW cryo-EM grids

DOI

dx.doi.org/10.17504/protocols.io.b3rqm5w

Preparation of LRRK1 RCKW cryo-EM grids



Mariusz Matyszewski¹, David Snead¹

¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093



Mariusz Matyszewski

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b3rqm5w>

Protocol Citation: Mariusz Matyszewski, David Snead 2022. Preparation of LRRK1 RCKW cryo-EM grids. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b3rqm5w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 12, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 56848

Keywords: cryo-EM, LRRK2, structural biology, ASAPCRN, Irrk1 rckw grids for cryo, preparation of Irrk1 rckw cryo, Irrk1 rckw grid, Irrk1 rckw cryo, em grids protocol, cryo

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000519

MJFF

Grant ID: 18321

Abstract

Protocol used to create LRRK1 RCKW grids for cryo-EM used in Snead, Matyszewski, Dickey et al.

Materials

LRRK1 Buffer:

- [M] 20 millimolar (mM) HEPES pH 7.4
- [M] 80 millimolar (mM) NaCl
- [M] 0.5 millimolar (mM) TCEP
- [M] 2.5 millimolar (mM) MgCl₂
- [M] 20 micromolar (μM) GDP

Note: please change salt as needed to maintain final salt of 80 mM NaCl

Troubleshooting



Safety warnings

! For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).
Take proper precautions while freezing grids.

Before start

Decide which protein concentration to use, and create the proper LRRK1 buffers in order to obtain the right salt concentration (80 mM NaCl).



Preparing Sample

10m


- 1 Spin down purified LRRK1 RCKW.  10000 rcf, 4°C, 00:10:00 , (can be faster)
Leave protein on ice afterwards.

Note



For best results, reduce the amount of time between spinning and freezing samples.

Freezing Grids

20s

- 2 Plasma clean grids.
We used UltrAuFoil Holey Gold 1.2/1.3 300 mesh grids and plasma cleaned them in the Solarus II (Gatan) using the QuantiFoil Au preset.
- 3 Dilute samples to desired concentration in the **LRRK1 buffer**. Make sure final salt is at 80 mM NaCl.
For best results, make  10 μL samples, good for freezing 2 grids. This is to minimize time spent outside of storage buffer, reducing aggregation.


Note

In our publication, we used concentrations ranging from  2 micromolar (μM) to  6 micromolar (μM) . Lower concentrations were used for samples collected on using a tilted stage.

- 4 Apply protein to grids and plunge freeze.
We used a Vitrobot (FEI) to blot away excess sample and plunge freeze



Note

Each Vitrobot is slightly different, so use your preferred settings. We used  3.5 μL of sample and a 4 second blot at force 20, but it is unlikely to work on other Vitrobots.

- 5 Store grids in liquid nitrogen until ready for imaging.