

Oct 15, 2024

Preparation of Sylgard® 184 PDMS

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpzk7plzp/v1>



Guillermo Pérez-Mateluna¹

¹Pontificia Universidad Católica de Chile

Wenzel Lab



Guillermo Pérez-Mateluna

Pontificia Universidad Católica de Chile

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpzk7plzp/v1>

Protocol Citation: Guillermo Pérez-Mateluna 2024. Preparation of Sylgard® 184 PDMS. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpzk7plzp/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 25, 2024

Last Modified: October 15, 2024

Protocol Integer ID: 108261

Keywords: PDMS, polydimethylsiloxane, pdm, silicon wafer, pdms this protocol, safety methods for the preparation, preparation, service providers of the laboratory, laboratory

Abstract

This protocol outlines the main steps, materials, equipment, and safety methods for the preparation of Polydimethylsiloxane (PDMS) on silicon wafers (or similar). It is intended for all users, visitors, and service providers of the laboratories of the Institute of Biological and Medical Engineering (IIBM).

Guidelines

Usar guantes de nitrilo

Antiparras protectoras para los ojos

Delantal de laboratorio



Materials

Materiales

 SLYGARD 184 Biosearch Technologies

Wafer with chip

Cup

Petri dish

Scalpel

Fork for mixing the mixture

#10 scalpel blade

Equipos

Drying oven

Vacuum generator located in Laboratory No. 2 ID EQ_184

Balance ID: EQ_179

Safety warnings

 Todas los insumos que toquen PDMS, tales como, recipientes, agitadores y papeles deben ser eliminadas en el basurero rojo de sustancias peligrosas

Before start

El PDMS debe realizarse en el laboratorio 2 en la zona delimitada para trabajar con sustancias peligrosas



- 1 Place a weighing dish or plastic cup on the balance, press TARE, and prepare the required amount of Sylgard® 184 PDMS by pouring the necessary quantities of base and curing agent in a 10:1 ratio (for example, 5 g of base + 0.5 g of curing agent).
- 2 Mix both components thoroughly with a plastic fork until the compounds are fully homogenized (transparent).




2.1



- 3 After mixing the compounds, bubbles may form, so the PDMS should be degassed by placing it in a vacuum chamber located in Laboratory No. 2.


3.1



- 4 Place the container with PDMS inside a vacuum chamber, close the valve, turn on the vacuum pump, and allow it to degas until all the bubbles disappear.
- 5 Slowly open the vacuum chamber valve to let air in and release the negative pressure.
- 6 Pour the mixture onto the silicon wafer (or another container/mold) that will be used. Keep the silicon wafer inside its container and avoid spilling PDMS near the edge of the wafer. Ensure no PDMS remains stuck inside the container.

The wafer with the chip should never be handled directly, not even with gloves, as it can easily get dirty or scratched. Always handle it inside its case.

6.1



- 7 Place the wafer, still inside its case with the lid on, into the drying oven at 80°C for

45m


 00:45:00 .


Ensure that the container is properly leveled to guarantee a uniform thickness of the PDMS.

- 8 Allow the wafer and the PDMS to cool in a clean environment with the case lid on.

- 9 When the PDMS is cold (the next day), it should be carefully cut along the defined surface to avoid damaging or scratching the wafer. Using a scalpel with a #10 blade, gently lift the PDMS with tweezers to detach it from the wafer and the plate.

9.1





10 Now the PDMS can be mounted on a microscope slide, well plate, or Petri dish.