

May 07, 2018

Preparation of the sgRNA-Barcode Amplicon Library

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.pufdntn

Shiqi Xie¹, Anne Cooley¹, Gary Hon¹

¹Green Center for Reproductive Biology, Department of Obstetrics and Gynecology, UT Southwestern Medical Center



Shiqi Xie

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.pufdntn

External link: <https://doi.org/10.1371/journal.pone.0198635>

Protocol Citation: Shiqi Xie, Anne Cooley, Gary Hon 2018. Preparation of the sgRNA-Barcode Amplicon Library. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.pufdntn>

Manuscript citation:

Xie S, Cooley A, Armendariz D, Zhou P, Hon GC (2018) Frequent sgRNA-barcode recombination in single-cell perturbation assays. PLoS ONE 13(6): e0198635. doi: [10.1371/journal.pone.0198635](https://doi.org/10.1371/journal.pone.0198635)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 29, 2018

Last Modified: May 07, 2018

Protocol Integer ID: 11879



Guidelines

This protocol is used to assess the sgRNA-to-barcode linkage in Mosaic-seq. The 3kb fragment can be amplified from either plasmids or genomic DNA extract from cells after lentiviral infection. The fragments are then circularized and used to amplify a smaller amplicon (about 270bp) which is compatible with Illumina sequencing.

PCR to get the 3kb fragment

- 1 Prepare the following mixture:

Reagent	Vol.(μ l)
NEBNext mixture	25
PCR1 primers +/- (10 μ M)	1/1
H ₂ O	Up to 50
Plasmid / gDNA	indicated
Total	50

Note

Use 1000ng genomic DNA or 100ng plasmid DNA for the PCR

- 2 PCR cycles

98°C	30sec	
98°C	10sec	
53°C	30sec	22 or 12 cycles
72°C	3min	
72°C	5min	

Note

22 cycles for gDNA, 12 cycles for plasmids

- 3 Gel extract the PCR product:
 - Run the PCR products on 1% agarose gel. Include vector control in another lane.
 - Excise 3kb fragment and perform MinElute PCR cleanup. Elute in 2 × 10 μ L of EB.

PNK treatment to phosphorylate 5' ends

- 4 Set up the following reaction:

Reagent	Vol.(μ l)
PCR1 product in water	39
10X T4 PNK buffer	5
10 mM ATP	5
T4 PNK	1



5 Incubate at 37°C for 30 min.

🔥 37 °C 30 min

6 Perform MinElute PCR cleanup. Elute in 2×10ul. Spec the DNA concentration by Qubit.

Self-ligation with Quick Ligase

7 Setup the following reaction:

Reagent	Vol.(μl)	
50 ng PCR1 product in water	78	
2X Quick ligase buffer	80	
Quick ligase	2	

8 Incubate at 25°C for 2 hours.

🔥 25 °C 2 hours

9 Perform MinElute PCR cleanup. Elute in 2×10 μl EB.

PlasmidSafe digestion

10 Setup the following reaction:

Reagent	Vol.(μl)
DNA in water	42
25 mM ATP	2
10X Reaction buffer	5
Plasmid-Safe DNase	1

11 Incubate at 37C for 30 min.

🔥 37 °C 30 min

12 Perform MinElute PCR cleanup. Elute in 2×10ul. Spec the DNA concentration by using Qubit.

2nd PCR

13 Setup the following reaction:



Reagent	Vol.(μ l)
KAPA HiFi HS MasterMix	25
BarLib P5 Stagger mix (10 mM)	1
BarLib P7 N7XX (10 mM)	1
Circulated DNA	18
Water	to 50

14 Perform PCR by using the following program:

95°C	3min	
98°C	20sec	
60°C	15sec	16 cycles
72°C	1min	
72°C	3min	

15 Run samples on 1% agarose gel. Excise band at 400bp

16 Spec DNA concentration by Qubit.

17 Run the samples on Tapestation (D1000 tape) to check the final size of the library.