

Jul 24, 2019 Version 2

## Prepare 2L of B-Broth for culturing bacteria V.2

DOI

[dx.doi.org/10.17504/protocols.io.5spg6dn](https://dx.doi.org/10.17504/protocols.io.5spg6dn)

Cancer Research UK / Wellcome Gurdon Institute media kitchen<sup>1</sup>

<sup>1</sup>Cancer Research UK / Wellcome Trust Gurdon Institute



Cristian Riccio

University of Cambridge

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.5spg6dn](https://dx.doi.org/10.17504/protocols.io.5spg6dn)

**Protocol Citation:** Cancer Research UK / Wellcome Gurdon Institute media kitchen 2019. Prepare 2L of B-Broth for culturing bacteria. [protocols.io https://dx.doi.org/10.17504/protocols.io.5spg6dn](https://dx.doi.org/10.17504/protocols.io.5spg6dn)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 24, 2019

**Last Modified:** July 24, 2019

**Protocol Integer ID:** 26159


### Abstract

Prepare 2L of B-Broth for culturing bacteria

## Materials


### MATERIALS

 double distilled water (ddH<sub>2</sub>O)

 Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical **Fisher Scientific Catalog #S/3160/65**

 Bacto™ Tryptone Life Technologies **VWR International (Avantor) Catalog #211705**

## Safety warnings

 Make sure your know how to use the autoclave.



1



2

Ingredients	Quantity	
Bacto tryptone	20g	
NaCl		10g
Double distilled H2O	up to 2L	

3

1	Dissolve Bacto tryptone and NaCl in 1L double distilled H2O in a 2L measuring cylinder
2	Make up to final volume using double distilled H2O
3	Aliquot in 100 ml bottles, label and autoclave