

Oct 11, 2022

🌐 Prepare agarose pad for microscopy (2022 iGEM)

🔗 Forked from [Fluorescence microscopy of Chlamydomonas reinhardtii - mCherry | Chlorophyll](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5jze5l1b/v1>

igem¹


¹Fudan University

2022



Team Fudan iGEM

Fudan University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5jze5l1b/v1>

External link: <https://doi.org/10.1371/journal.pone.0192433>

Protocol Citation: igem 2022. Prepare agarose pad for microscopy (2022 iGEM). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8epv5jze5l1b/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 10, 2022

Last Modified: October 11, 2022

Protocol Integer ID: 71105

Keywords: agarose pad for microscopy, fluorescence microscopy, microscopy, agarose pad, agarose

Abstract

This protocol describe the steps to perform fluorescence microscopy for *E. coli*

Guidelines

Cells should be cultured to a late log phase to increase cell number for microscopy.

Before start

- Culture cells in liquid media


Cell growth

- 1 Grow cells to late-log phase density (e.g. $OD_{600} > 0.6$).

Agarose pads

- 2 Prepare a 1% agarose in distilled water, by adding sufficient low-melting agarose to water and dissolve agarose by heating in a microwave for 30-60 seconds.

16h



- 3 Prepare a glass slide with a double adjacent slides barrier (see below), forming a wall around the pad slide. Add slides on top of the poured agarose solution (1-1.2 mL) to allow formation of a smooth surface. It takes time for 1% agarose solution to cool down to ~40 degree, and the solution is viscous but not runny.



- 4 Wait for 45-60 minutes. Carefully remove adjacent and top slides, leaving the bottom slide.
- 5 Cut the pad into 15×15 mm square pierces, and carefully transfer the pad pierces to other cleared slide (left one agarose pad on the existing glass slide, if it was clean).
- 6 Add 2-5 μ L bacteria culture onto the agarose pad, and cover it was a 18×18 mm #1.5 square coverslip. Seal with wax if needed.



Microscopy GFP

5m

- 7
 - Place the slide onto the microscope stage, with the coverslip facing the objective.
 - Use a laser at 488 nm to excite GFP.
 - Record images.
- 8 Analyze the images. For our microscope with 150x objective, the pixel size is 0.0806 μ m.

5m