

Jan 30, 2025

Version 1

Prepare gDNA for CRISPR Screen V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk9rb1v5r/v1>

Christopher Sansam¹

¹Oklahoma Medical Research Foundation



Christopher Sansam

Oklahoma Medical Research Foundation

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk9rb1v5r/v1>

Protocol Citation: Christopher Sansam 2025. Prepare gDNA for CRISPR Screen. protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlk9rb1v5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 29, 2025

Last Modified: January 30, 2025

Protocol Integer ID: 119249

Keywords: gdna for crispr screen, critical step for crispr screening application, crispr screening application, using quickextract dna extraction solution, dna extraction, quickextract dna extraction solution, extracting gdna, preparation of genomic dna, crispr screen, final dna pellet, phase separation efficiency during organic extraction, purity dna, phase lock gel technology, gdna, quality dna suitable for downstream analysis, quality dna, genomic dna, dna, organic extraction, rnase treatment, extraction, such as pcr amplification, pcr amplification, pcr, chloroform extraction, phase separation efficiency

Funders Acknowledgements:

Oklahoma Center for Adult Stem Cell Research (OCASCR)

Grant ID: January 2022

National Institutes of Health

Grant ID: R01GM121703

Presbyterian Health Foundation

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

The preparation of genomic DNA (gDNA) is a critical step for CRISPR screening applications, ensuring high-quality DNA suitable for downstream analysis such as PCR amplification and sequencing. This protocol details a streamlined method for extracting gDNA from approximately 30 million cells using QuickExtract DNA Extraction Solution. The process involves cell lysis, RNase treatment, phenol-chloroform extraction, and isopropanol precipitation to ensure the isolation of high-purity DNA. Phase Lock Gel technology is employed to enhance phase separation efficiency during organic extraction. The final DNA pellet is resuspended in EB buffer and quantified before proceeding to PCR.

Materials

Reagents and Chemicals

Quick Extract DNA Extraction Solution (Lucigen, Catalog #QE0905T)

A specialized buffer designed for rapid and efficient DNA extraction from cells. It facilitates cell lysis and initial DNA release for downstream processing.

PureLink™ RNase A (20 mg/mL) (Thermo Fisher, Catalog #12091021)

An enzyme that degrades RNA contaminants in the extracted genomic DNA, ensuring RNA-free samples for accurate downstream applications.

UltraPure™ Phenol:Chloroform:Isoamyl Alcohol (25:24:1, v/v) (Thermo Fisher Scientific, Catalog #15593031)

A reagent used for organic phase separation, effectively removing proteins and other cellular contaminants from the DNA extract.

Chloroform

Used in the extraction process to enhance the separation of DNA from proteins and lipids, ensuring high-purity DNA recovery.

Phase Lock Gel Heavy 2 mL (Quantabio, Catalog #2302830)

A gel-based barrier that simplifies phenol-chloroform extraction by preventing interphase contamination and improving DNA recovery.

Ammonium Acetate (7.5 M solution)

Used in DNA precipitation to remove proteins and enhance nucleic acid purity by reducing the solubility of DNA.

Isopropanol (100%)

Facilitates DNA precipitation, allowing efficient recovery of genomic DNA from the aqueous phase. The use of 100% isopropanol ensures optimal DNA yield.

Ethanol (70%)

Used for washing the DNA pellet, removing residual salts and impurities, and replacing isopropanol for easier redissolution.

EB Buffer (10 mM Tris-Cl, pH 8.5)

A low-salt buffer used for DNA resuspension, providing stable conditions for DNA storage and downstream applications.

Equipment and Consumables

Eppendorf Thermomixer R (1.5 mL block)

- Maintains precise temperature and agitation settings for enzymatic digestion and incubation steps.

Microcentrifuge (capable of 16,000 x g spin speeds at 4°C and room temperature)

Essential for phase separation, DNA precipitation, and pellet washing.

Sherlock Tube Cap (USA Scientific)

Used to securely cap tubes during heating and vortexing steps.

Pipettes and Pipette Tips

Used for precise liquid handling during all stages of DNA preparation.

1.5 mL Microcentrifuge Tubes

Used for sample processing, extraction, and DNA storage.



Protocol materials

⊗ Quick Extract DNA Extraction Solution **Lucigen Catalog #QE0905T**

⊗ PureLink® RNase A (20 mg/mL) **Thermo Fisher Catalog #12091021**

⊗ UltraPure™ Phenol:Chloroform:Isoamyl Alcohol (25:24:1, v/v) **Thermo Fisher Scientific Catalog #15593031**

⊗ Phase Lock Gel Heavy 2 ml **Quantabio Catalog #2302830**



Lyse Cells in QuickExtract Buffer

1h 30m

- 1 Add  500 μL
 Quick Extract DNA Extraction Solution **Lucigen Catalog #QE0905T** to approximately 30 million cell pellet.
- 2 Mix by vortexing vigorously for 15 seconds.
- 3 Cap with Sherlock tube cap (USA Scientific). Heat in Eppendorf Thermomixer R (1.5 ml block) at  65 °C  1000 rpm for  01:00:00
- 4 Vortex vigorously for 15 seconds and centrifuge at full speed in microcentrifuge for 5 seconds.
- 5 Heat at  98 °C for 5 minutes.
- 6 Vortex vigorously for 15 seconds and centrifuge at full speed in microcentrifuge for 5 seconds.
- 7 Incubate on ice for 1 minute.
- 8 Add  5 μL
 PureLink® RNase A (20 mg/mL) **Thermo Fisher Catalog #12091021** .
Incubate in Thermomixer at  37 °C  1000 rpm for  00:30:00
- 9 Keep on ice.

1h

30m

Phenol-Chloroform Extraction

15m



- 10 Add  500 μL
 UltraPure™ Phenol:Chloroform:Isoamyl Alcohol (25:24:1, v/v) **Thermo Fisher Scientific Catalog #15593031**
- 11 Vortex vigorously for 1 minute.
- 12 Centrifuge at full speed in microcentrifuge ( 16873 rcf, 25°C, 00:05:00) for 5 minutes. 5m
- 13 Prepare a phase lock tube ( Phase Lock Gel Heavy 2 ml **Quantabio Catalog #2302830**) by centrifuging at 12,000 to 16,000 x g for 15 seconds.
- 14 In hood, carefully transfer the QuickExtract lysate into the the phase lock tube by decanting or pipetting.
- 15 Centrifuge  16000 x g, 25°C, 00:05:00 to separate the phases. 5m
- 16 Add  500 μL to upper phase. Mix thoroughly without vortexing or disrupting the Phase Lock Gel barrier.
- 17 Centrifuge  16000 x g, 25°C, 00:05:00 to separate the phases. 5m
- 18 Carefully transfer by decanting or pipetting the upper aqueous layer to a new 1.5 ml centrifuge tube.

Isopropanol precipitation

10m

- 19 Add ammonium acetate ($[M]$ 7.5 Molarity (M)) to a final concentration of $[M]$ 2.5 Molarity (M) (i.e. 1/2 volume ~  250 μL).
- 20 Add 0.7X volume of Isopropanol (~  525 μL) and vortex briefly to mix.



- 21 Incubate at room temperature for  00:10:00
- 22 Centrifuge the sample immediately at 10,000–15,000 x g for 15–30 min at 4°C.
- 23 Carefully decant the supernatant without disturbing the pellet.
- 24 To wash the DNA pellet, add 1 ml of room-temperature 70% ethano. This step removes co-precipitated salts and replaces isopropanol with the more volatile ethanol, facilitating easier redissolution of the DNA.
- 25 Centrifuge at 10,000–15,000 x g for 5–15 min at 4°C.
- 26 Carefully remove the supernatant without disturbing the pellet.
- 27 Quick spin in microcentrifuge to bring residual Ethanol to bottom of tube.
- 28 Using a p200 pipette, remove residual Ethanol. Leave pellet.
- 29 After briefly air drying (< 2 minutes), resuspend in  100 µL EB buffer (10 mM Tris-Cl, pH 8.5).
- 30 Let DNA solution sit at room temperature overnight to rehydrate. Vortex well.
- 31 Quantify DNA and proceed to PCR.

10m