

Mar 07, 2023

Version 1

## Prepare Samples for Miseq V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbyknyvpk/v1>



FishFloorUCL<sup>1</sup>

<sup>1</sup>University College London

Talia Pittman

FishFloorUCL



Talia Pittman

University College London (UCL)

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.81wgbyknyvpk/v1>

**Protocol Citation:** FishFloorUCL 2023. Prepare Samples for Miseq. protocols.io

<https://dx.doi.org/10.17504/protocols.io.81wgbyknyvpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 02, 2023

**Last Modified:** March 22, 2023

**Protocol Integer ID:** 77981

**Keywords:** miseq, sequencing, dna sequence multiple sample, fish with mosaic genotype, samples for miseq, miseq nano v2 chip 300 cycle, miseq multiplexes sequence, genotyping fish, miseq nano v2 chip, dna sequence, mosaic genotype, multiple sample, dna, genotype, miseq, sample

## Abstract

(MiSeq nano V2 chip 300 cycle), MiSeq multiplexes sequences, allowing you to DNA sequence multiple samples at once. This is particularly useful for genotyping fish with mosaic genotypes, such as F0 knockouts.

@FishFloorUCL

## Protocol materials

 Qubit™ 1X dsDNA High Sensitivity (HS) Assay Kit **Invitrogen - Thermo Fisher Catalog #Q33231**

 ExoSAP-IT® PCR Product Cleanup Reagent **Thermo Fisher Catalog #78201.1.ML**

## Design Miseq Primers

### 1 Design primers surrounding your site of interest


- Optimal amplicon size **200bp**
- As far as possible, the site of interest should be in the middle of the designed amplicon. (This is so that it has the best chance of being read by both the forward and reverse reads)
- If possible, try and avoid designing primers in non-coding regions of DNA since these regions are more likely to have variation and SNPs, which might affect the efficacy of your primer.

*I usually use the **NCBI Primer-Blast** tool. Enter ~300bp of DNA sequence around your region of interest under 'PCR Template' and set the 'Range' for the forward and reverse primer so that your site of interest is more or less in the middle of the PCR product.*

If you previously used a set of primers to check your TALEN/CRISPR efficiency using Melt Curve Analysis, or if you have your primer sequence from **CHOPCHOP**, you can use the same sequence to design the primers for MiSeq. Anecdotally this seems to work well for most people on their first try.

### 2 Add miseq tags/overhangs to your primer sequence.

- Universal forward (F) tag: 5' **TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG** 3'
- Universal reverse (R) tag: 5' **GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG** 3'



*These overhang sequences correspond to the Illumina nextera transposase overhangs and will be compatible with Illumina kits using Nextera XT*



- 3 (Optional) Test primers (including tag) for secondary structure formation (e.g., dimers and hairpins) using this online [IDT OligoAnalyzer tool](#)

NB. It is unlikely that you will have high scores due to the length of the primers and amplicon restrictions. Just pick the best

- 4 Order primers - standard desalted primers are fine for this application.

- 5 If the primers arrive as dry pellets, reconstitute the primers in autoclaved distilled or nuclease-free water to create a stock solution and working solution.

- [M] 100 micromolar ( $\mu\text{M}$ ) **stock solution** by adding nuclease-free water. If  $x =$  nMoles of your primer (check the tube), resuspend the dry pellet in  $10x \mu\text{L}$  of water
- [M] 10 micromolar ( $\mu\text{M}$ ) **working solution**

## (Optional) Test Primers

11m 30s

6

11m 30s

### Note

Testing your primers is optional. Although, it is recommended, especially if you are having trouble getting your PCR to work or if you have a lot of samples and want to ensure your primers work before running the entire experiment.

For optimal Miseq results, it is best to use a high-fidelity polymerase. However, at this stage, you can use a cheaper polymerase. Feel free to use a DNA polymerase and PCR protocol that works for you and follow the instructions with the kit. I usually use standard *Taq* DNA Polymerase.

 Taq DNA Polymerase recombinant **Invitrogen - Thermo**  
**Fisher Catalog #10342**

### Run a PCR

I have designed a handy [master mix calculator](#) for *Taq*, Platinum *Taq* and Pre-mix Schier PCR mix.

Here is an example of a PCR mix for a single well with a  10  $\mu\text{L}$  total volume

| A                                           | B                  |
|---------------------------------------------|--------------------|
| Autoclaved, distilled (nuclease free) water | 6.46 $\mu\text{L}$ |
| 10X PCR Buffer, -Mg                         | 1 $\mu\text{L}$    |



| A                                 | B       |
|-----------------------------------|---------|
| 50 mM MgCl <sub>2</sub>           | 0.3 µL  |
| 10 mM (each) dNTP Mix             | 0.2 µL  |
| 10 µM forward primer              | 0.5 µL  |
| 10 µM reverse primer              | 0.5 µL  |
| <i>Taq</i> DNA Polymerase         | 0.04 µL |
| Template DNA (from 100µL HotShot) | 1 µL    |

Note, *Taq* is unstable and should be taken out of the freezer for the shortest amount of time possible. It is sensible to prepare PCR mix on ice.

Suggested PCR programme:

1. 🔥 95 °C , ⌚ 00:05:00
2. 🔥 95 °C , ⌚ 00:00:30 (*Denaturation step*)
3. Temperature gradient\* around 🔥 60 °C , ⌚ 00:00:30 (*Annealing step*)
4. 🔥 72 °C , ⌚ 00:00:30 (*Extension step*)
5. Repeat steps 2-4, 25-35 times
6. 🔥 72 °C , ⌚ 00:05:00 (*Final extension step*)
7. 🔥 10 °C hold

\*For your test run, you should use a gradient PCR to test the ideal temperature for your primers (Note. By default, CHOPCHOP designs primers for an optimal melting temperature of 60°C). We usually run a gradient PCR for five different temperatures, with 60°C being optimal, then varying 1°C-1.5°C above and below.

## 7 Run on a 1% agarose gel

### Expected result

Assuming a 200bp amplicon, you should expect a final PCR product of 267bp on your gel. If you don't see a band in the expected region, you may need to re-design your primers.

## Prepare samples for miseq

11m 30s

## 8 Amplify up PCR products to analyse by MiSeq.

11m 30s

## Note

Note **you can only submit up to 96 samples per amplicon** on one MiSeq run.

**We recommend using a high-fidelity Taq polymerase.** We recommend using  Taq DNA Polymerase Platinum **Invitrogen - Thermo Fisher** or Phusion® High-Fidelity DNA Polymerase from New England Biolabs (M0530), but any high-fidelity polymerase will do.

I have designed a handy [master mix calculator](#).

Here is an example of a PCR mix for a single well with  10 µL total volume

| A                                              | B      |
|------------------------------------------------|--------|
| Autoclaved, distilled (or nuclease free) water | 7 µL   |
| 10X PCR Buffer, -Mg                            | 1 µL   |
| 50 mM MgCl <sub>2</sub>                        | 0.3 µL |
| 10 mM (each) dNTP Mix                          | 0.2 µL |
| 10 µM forward primer                           | 0.2 µL |
| 10 µM reverse primer                           | 0.2 µL |
| Platinum Taq DNA Polymerase                    | 0.1 µL |
| Template DNA                                   | 1 µL   |

Note, *Taq* is unstable and should be taken out of the freezer for the shortest amount of time possible. It is sensible to prepare PCR mix on ice.

Suggested PCR programme:

1.  95 °C ,  00:05:00 (*Initial Denaturation Step*)
2.  95 °C ,  00:00:30 (*Denaturation step*)
3.  60 °C \*  00:00:30 (*Annealing step*)
4.  72 °C ,  00:00:30 (*Extension step*)
5. Repeat steps 2-4, 25-35 times
6.  72 °C ,  00:05:00 (*Final extension step*)
7.  10 °C hold

\*Or  $T_m$  of your given primer pair

- 9 **Run on a 1% gel** to confirm a single band at the anticipated weight. This is just to check the PCR worked. To save time, just run a couple of your samples.

10 **Clean up your samples**

30m

Traditional column-based PCR clean-up kits would work fine, but we recommend using ExoSAP,



ExoSAP-IT<sup>®</sup>; PCR Product Cleanup Reagent **Thermo**  
**Fisher Catalog #78201.1.ML**

. Although it is more expensive (~£30 per plate), it is much more practical, especially for many samples.

1. Transfer  5  $\mu\text{L}$  of your PCR product into a new plate
2. Add  2  $\mu\text{L}$  of ExoSAP reagent to each well
3. (Vortex optional) Spin down & run the ExoSAP programme thermocycler.
  -  00:15:00 at  37 °C – Degrade the remaining primers and nucleotides
  -  00:15:00 at  80 °C – Inactivate ExoSAP

- 11 **Check the concentration of your clean PCR product.** This can be done using Nanodrop or Qubit. I would recommend Qubit due to its superior accuracy and reliability. We recommend only testing a few of your samples to save time and money.

We usually use the following kit from Thermo Fisher, which comes with a working solution.




Qubit<sup>™</sup> 1X dsDNA High Sensitivity (HS) Assay Kit **Invitrogen - Thermo**  
**Fisher Catalog #Q33231**

Setup and label the lids of the required number of assay tubes for your **samples** plus **two standards** (Use thin-wall, clear, 0.5-mL Qubit Assay PCR tubes (Cat. No. Q32856)).

1. Add  190  $\mu\text{L}$  Qubit<sup>™</sup> 1X dsDNA working solution to each **standard** tube
2. Add  10  $\mu\text{L}$  of each Qubit standard to the appropriate **standard** tube
3. Add  199  $\mu\text{L}$  Qubit<sup>™</sup> 1X dsDNA working solution to each **sample** tube
4. Add  1  $\mu\text{L}$  of each sample to the appropriate tube

All tubes should have a final volume of  200  $\mu\text{L}$



Read concentration using the Qubit machine and note them down.

## 12 **Prepare your samples** for miseq:

You need to submit a minimum of  6 µL **at a final concentration of 15-25 ng/µL**

Label your plate/tubes with your name, the date and '*For Miseq*' and leave the samples with the **Payne lab in the Cancer Insititute or in the Fish Facility B03 fridge in the anatomy building**, whichever is most convenient for you. Once you have dropped off your samples, email Alexandra Lubin, letting her know.

**Note**

You could save time and plastic by combining **clean-up** with dilution to **15-25 ng/μL**.

**Step 1:** Measure the concentration of the PCR product before clean-up

**Step 2:** In a new plate/tubes, add your PCR product, ExoSAP and nuclease-free water (if necessary). Volumes will differ according to the concentration of your samples. Calculate how much you need to dilute your PCR product according to the concentration you measured. Add  2 μL ExoSAP for every  5 μL you add of your PCR product and top up with nuclease-free water where needed. The following table shows some examples:

|  | A             | B                 | C            | D                         |
|--|---------------|-------------------|--------------|---------------------------|
|  | Concentration | PCR product (vol) | ExoSAP (vol) | nuclease-free water (vol) |
|  | < 15-25 ng/μL | 5 μL              | 2 μL         | 0 μL                      |
|  | 30-50 ng/μL   | 5 μL              | 2 μL         | 3 μL                      |
|  | 45-75 ng/μL   | 2.5 μL            | 1 μL         | ~6 μL                     |

**Step 3:** (Vortex optional) Spin down & run the ExoSAP programme thermocycler.

-  00:15:00 at  37 °C – Degrade the remaining primers and nucleotides
-  00:15:00 at  80 °C – Inactivate ExoSAP

The resulting plate/tubes will be ready to submit for miseq.