

Nov 02, 2020

Version 1

Preparing biological samples for metabarcoding V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bkyrkxv6>

Tim Regan¹

¹The Roslin Institute, University of Edinburgh



Tim Regan

The Roslin Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bkyrkxv6>

Protocol Citation: Tim Regan 2020. Preparing biological samples for metabarcoding. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bkyrkxv6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 07, 2020



Last Modified: November 02, 2020

Protocol Integer ID: 41713

Keywords: Metabarcoding, metagenomics, DNA extraction, preparing biological sample, pcr inhibition during amplicon library preparation, amplicon library preparation, preparation of biological sample, dna from lysate, biological sample, using metapolyzyme, dna elution, efficient dna release, concentration estimates of dna elution, amplicon, qiagen dneasy kit, pcr inhibition, enzymatic lysi, qiagen buffer atl, metabarcoding, sample, preparation, dna,

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the preparation of biological samples (specifically from a marine environment e.g. hatchery or RAS unit) for amplicon sequencing. Starting with a biological sample stored in Qiagen buffer ATL, or similar, it begins with a bead beating process to homogenise the sample. Enzymatic lysis using Metapolyzyme and Proteinase K are employed to ensure efficient DNA release. The Qiagen DNeasy kit is used to column extract DNA from lysates. Following concentration estimates of DNA elutions, samples are diluted >1:10 to avoid PCR inhibition during amplicon library preparation.

Guidelines

In every step following enzymatic digestion of samples (and in general), ensure samples are kept at 4C to maximise sample stability.

Freeze DNA samples if not being used for >1 week following extraction.

Otherwise, storing DNA at 4C in fridge is preferable.



Materials

MATERIALS

✕ Buffer AL **Catalog #19075**

✕ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

✕ Buffer ATL **Qiagen Catalog #19076**

✕ Proteinase K, 100mg **Promega Catalog #V3021**

✕ PBS

✕ Ethanol 70%

✕ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

✕ UltraPure®; DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977015**

✕ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

Centrifuge.

Bead beater.

Incubator (for 37C and 56C).

Pipettes and tips.

Protocol materials

✕ Proteinase K, 100mg **Promega Catalog #V3021**

✕ UltraPure®; DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977015**

✕ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

✕ Buffer AL **Catalog #19075**

✕ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

✕ Ethanol 70%

✕ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

✕ Buffer ATL **Qiagen Catalog #19076**

✕ PBS

✕ Buffer ATL **Qiagen Catalog #19076**

✕ Buffer AL, Lysis buffer **Qiagen Catalog #19076**

Safety warnings

❗ Refer to manufacturer's MSDS information for each reagent used to ensure appropriate and safe use.

Before start

Ensure leaving time for samples to thaw if frozen. Avoid leaving samples thaw for too long as this may lead to degradation.



Bead beating

45m

- 1 Starting with biological sample (filter, swab, water, biofilm, tissue etc.) stored in Qiagen Buffer ATL (or similar), transfer up to 1 mL to Matrix A bead tube.

30m

 Buffer ATL **Qiagen Catalog #19076** Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

- 2 Perform bead beating in a disruptor at at 5.0 M/s (speed) for 00:00:40 x2 (ensure tube looks homogenous).

15m

Enzymatic digestion

2h 15m

- 3 Add 5 μL of Metapolyzyme to each tube and vortex briefly.

2h 15m

Incubate samples at 37 °C for 02:00:00 .

 MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

- 4 Add 20 μL of Proteinase K to each tube, vortex for 00:00:10 , then incubate at

16h

56 °C Overnight .

 Proteinase K, 100mg **Promega Catalog #V3021**

DNA extraction

- 5 Vortex samples for 00:00:15 and centrifuged 13000 x g for 00:01:00 .

- 6 Transfer the supernatant from each tube (up to 900 μL) into a new tube and centrifuged at 13000 x g, 00:01:00 .



- 7 Transfer up to $600\ \mu\text{L}$ of bead-free supernatant to a new $2\ \text{mL}$ tube.
- 8 Premix 70% ethanol and Qiagen lysis buffer AL 1:1 to add to sample at a ratio of 1:1:1 e.g. for 10 samples of $500\ \mu\text{L}$ each, premix $550\ \mu\text{L}$ of buffer AL and $550\ \mu\text{L}$ of 70% ethanol and add $1\ \text{mL}$ of ethanol/buffer AL mixture to each sample.
-  Buffer AL, Lysis buffer **Qiagen Catalog #19076**
- 9 Hereafter, the manufacturer's protocol for the Qiagen DNeasy Blood and Tissue kit is followed with some modifications:
- Load $< 600\ \mu\text{L}$ of lysate mixture (ATL, AL and EtOH) at a time into the column
 - Spin at $6000 \times g$, 00:01:00 and discard flow-through.
 - Repeat as necessary until all lysate is loaded on column e.g. mixture of $1500\ \mu\text{L}$ may take x3 initial spins and flow through discarding to complete column binding.
-  QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**
- 10
- Place the DNeasy Mini spin column in a new $2\ \text{mL}$ collection tube (provided), add $500\ \mu\text{L}$ Buffer AW1, and centrifuge at $6000 \times g$, 00:01:00 (8000 rpm).
 - Discard flow-through and collection tube.
- 11
- Place the DNeasy Mini spin column in a new $2\ \text{mL}$ collection tube (provided), add $500\ \mu\text{L}$ Buffer AW2, and centrifuge for at $20000 \times g$, 00:03:00 to dry the DNeasy membrane.
 - Discard flow-through and collection tube.
- 12 Perform final elution in $100\ \mu\text{L}$ of AE buffer.

Preparing concentration for library preparation

- 13 Check approximate concentration of extracted DNA using a Nanodrop.
- 14 Prepare 1:10 dilution of each extraction for PCR (to avoid PCR inhibition).



Perform further dilution of sample to a maximum final concentration of ~ [M] 1 ng/μl -

[M] 10 ng/μl

- 15 Use ~  50 ng of DNA in a  20 μL per sequencing library PCR reaction (see amplicon library PCR protocol).