

Sep 28, 2020 Version 2

## Preparing biological samples for metabarcoding V.2

DOI

[dx.doi.org/10.17504/protocols.io.bms8k6hw](https://dx.doi.org/10.17504/protocols.io.bms8k6hw)

Tim Regan<sup>1</sup>

<sup>1</sup>The Roslin Institute, University of Edinburgh



Tim Regan

The Roslin Institute

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bms8k6hw](https://dx.doi.org/10.17504/protocols.io.bms8k6hw)

**Protocol Citation:** Tim Regan 2020. Preparing biological samples for metabarcoding. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bms8k6hw> Version created by **Tim Regan**

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it works for us to generate metabarcoding libraries from biological samples (stored in e.g. Qiagen buffer ATL) for amplicon sequencing using MinION.

**Created:** September 28, 2020

**Last Modified:** September 28, 2020

**Protocol Integer ID:** 42560

**Keywords:** Metabarcoding, metagenomics, DNA extraction,



## Disclaimer

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol describes the preparation of biological samples (specifically from a marine environment e.g. hatchery or RAS unit) for amplicon sequencing. Starting with a biological sample stored in Qiagen buffer ATL, or similar, it begins with a bead beating process to homogenise the sample. Enzymatic lysis using Metapolyzyme and Proteinase K are employed to ensure efficient DNA release. The Qiagen DNeasy kit is used to column extract DNA from lysates. Following concentration estimates of DNA elutions, samples are diluted >1:10 to avoid PCR inhibition during amplicon library preparation.

## Guidelines

In every step following enzymatic digestion of samples (and in general), ensure samples are kept at 4C to maximise sample stability.

Freeze DNA samples if not being used for >1 week following extraction.

Otherwise, storing DNA at 4C in fridge is preferable.



## Materials

### MATERIALS

⊗ Buffer AL **Catalog #19075**

⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

⊗ Buffer ATL **Qiagen Catalog #19076**

⊗ Proteinase K, 100mg **Promega Catalog #V3021**

⊗ PBS

⊗ Ethanol 70%

⊗ MetaPolyzyme **Sigma Aldrich Catalog #MAC4L-5MG**

⊗ UltraPure®; DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977015**

⊗ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

### STEP MATERIALS

⊗ Buffer ATL **Qiagen Catalog #19076**

⊗ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

⊗ MetaPolyzyme **Sigma Aldrich Catalog #MAC4L-5MG**

⊗ Proteinase K, 100mg **Promega Catalog #V3021**

⊗ Buffer AL, Lysis buffer **Qiagen Catalog #19076**

⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

Centrifuge.

Bead beater.

Incubator (for 37C and 56C).

Pipettes and tips.

## Protocol materials

⊗ Proteinase K, 100mg **Promega Catalog #V3021**

⊗ Buffer AL **Catalog #19075**

⊗ Buffer ATL **Qiagen Catalog #19076**

⊗ Proteinase K, 100mg **Promega Catalog #V3021**

⊗ PBS

⊗ UltraPure<sup>®</sup>; DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977015**

⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

⊗ Ethanol 70%

⊗ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

⊗ Buffer ATL **Qiagen Catalog #19076**

⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

⊗ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

⊗ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

⊗ Buffer AL, Lysis buffer **Qiagen Catalog #19076**

⊗ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

⊗ Buffer ATL **Qiagen Catalog #19076**

⊗ Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

⊗ MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

⊗ Proteinase K, 100mg **Promega Catalog #V3021**

⊗ Buffer AL, Lysis buffer **Qiagen Catalog #19076**

⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

## Safety warnings

❗ Refer to manufacturer's MSDS information for each reagent used to ensure appropriate and safe use.


## Before start

Ensure leaving time for samples to thaw if frozen. Avoid leaving samples thaw for too long as this may lead to degradation.



## Bead beating


45m

- 1 Starting with biological sample (filter, swab, water, biofilm, tissue etc.) stored in Qiagen Buffer ATL (or similar), transfer up to  1 mL to Matrix A bead tube.

30m

 Buffer ATL **Qiagen Catalog #19076**

 Lysing Matrix A 2 mL tube **MP Biomedicals Catalog #SKU 116910050-CF**

- 2 Perform bead beating in a disruptor at at 5.0 M/s (speed) for  00:00:40 x2 (ensure tube looks somewhat homogenous).

15m


## Enzymatic digestion



2h 15m

- 3 Add  5  $\mu\text{L}$  of Metapolyzyme to each tube and vortex briefly.

2h 15m

Incubate samples at  37 °C for  02:00:00 .

 MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**




- 4 Add  20  $\mu\text{L}$  of Proteinase K to each tube, vortex for  00:00:10 , then incubate at



16h







 56 °C  Overnight .

 Proteinase K, 100mg **Promega Catalog #V3021**




## DNA extraction


- 5 Vortex samples for  00:00:15 and centrifuge a  13000 x g for  00:01:00 .








- 6 Transfer the supernatant from each tube (up to  900  $\mu\text{L}$  ) into a new tube and centrifuged at  13000 x g, 00:01:00 .

- 7 Transfer up to  600  $\mu\text{L}$  of bead-free supernatant to a new  2 mL tube.
- 8 Premix 70% ethanol and Qiagen lysis buffer AL 1:1 to add to sample at a ratio of 1:1:1 e.g. for 10 samples of  500  $\mu\text{L}$  each, premix  550  $\mu\text{L}$  of buffer AL and  550  $\mu\text{L}$  of 70% ethanol and add  1 mL of ethanol/buffer AL mixture to each sample.

 Buffer AL, Lysis buffer **Qiagen Catalog #19076**

- 9 Hereafter, the manufacturer's protocol for the Qiagen DNeasy Blood and Tissue kit is followed with some modifications:
- Load <  600  $\mu\text{L}$  of lysate mixture (ATL, AL and EtOH) at a time into the column
  - Spin at  6000 x g, 00:01:00 and discard flow-through.
  - Repeat as necessary until all lysate is loaded on column e.g. mixture of  1500  $\mu\text{L}$  may take x3 initial spins and flow through discarding to complete column binding.

 QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

- 10
- Place the DNeasy Mini spin column in a new  2 mL collection tube (provided), add  500  $\mu\text{L}$  Buffer AW1, and centrifuge at  6000 x g, 00:01:00 (8000 rpm).
  - Discard flow-through and collection tube.
- 11
- Place the DNeasy Mini spin column in a new  2 mL collection tube (provided), add  500  $\mu\text{L}$  Buffer AW2, and centrifuge for at  20000 x g, 00:03:00 to dry the DNeasy membrane.
  - Discard flow-through and collection tube.
- 12 Perform final elution in  100  $\mu\text{L}$  of AE buffer.



## Preparing concentration for library preparation

- 13 Check approximate concentration of extracted DNA using a Nanodrop.
- 14 Prepare 1:10 dilution of each extraction for PCR (to avoid PCR inhibition).



Perform further dilution of sample to a maximum final concentration of ~ [M] 1 ng/μl -

[M] 10 ng/μl

- 15 Use ~  50 ng of DNA in a  20 μL per sequencing library PCR reaction (see amplicon library PCR protocol).