



Jan 02, 2025

Preparing *C. elegans* for image acquisition

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzjr8lx1/v1>

Justin Donnelly¹

¹Arcadia Science

Arcadia Science



Arcadia Science

Arcadia Science

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzjr8lx1/v1>

Protocol Citation: Justin Donnelly 2025. Preparing *C. elegans* for image acquisition. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzjr8lx1/v1>

Manuscript citation:

Avasthi P, Borges AL, Cheveralls K, Donnelly J, Mets DG, Reiter T. (2025). An experimental and computational workflow to characterize nematode motility behavior. <https://doi.org/10.57844/arcadia-b89a-7e76>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 07, 2024

Last Modified: March 13, 2025

Protocol Integer ID: 109298

Keywords: c. elegans, worms, nematodes, worm, nematode, imaging, microscopy, phenotype, phenotyping, image acquisition, imaging, images on an upright widefield microscope, image, synchronized l1, upright widefield microscope, motility analysis, based motility analysis

Abstract

This protocol describes a procedure to prepare synchronized *C. elegans* for image-based motility analysis. It describes methods for maturing synchronized L1s to adulthood, transferring them to NGM plates without OP50 for imaging, and acquiring the images on an upright widefield microscope.

Materials

M9 buffer

6 cm Petri dishes with NGM(OP50)

6 cm Petri dishes with NGM (no OP50)

15 mL centrifuge tubes

Dissection scope

Nikon upright widefield microscope with Kinetix sCMOS camera



Maturation of synchronized L1s

50m

- 1 Once synchronized worms have hatched into L1 larvae, plate  200–300 μL synchronized worm solution per plate onto NGM (OP50) plates under a flame. If you resuspended worms in  1 mL M9 buffer, this should correspond to ~four plates per strain.
- 2 Cover plates and allow buffer to evaporate. Tap the plate gently to remove any bubbles. This process takes approximately  00:20:00 –  00:30:00 .
- 3 Incubate plates at room temperature for  96:00:00 .

50m

4d

Replating worms for imaging

50m

- 4 After incubation, worms on synchronized plates should have matured to young adults. Confirm this using the dissection scope.
- 5 For each strain, detach adults from plates using  0.5 mL sterile M9 buffer per plate under a flame. Swirl and tilt plate to ensure that buffer passes over the entirety of the plate, then use a pipette to transfer the buffer to a sterile 15 mL centrifuge tube. This process should remove the vast majority of worms from the plate.
- 6 Set tube upright in a rack to allow worms to settle under 1 G. This process takes approximately  00:20:00 –  00:30:00 .
- 7 After worms have settled, remove ~half of remaining supernatant.
- 8 Resuspend worms in remaining supernatant, then transfer to a fresh NGM plate **without** OP50. If there are too many worms for one plate, transfer half of resuspended worms to one plate, then plate an additional NGM plate after  00:15:00 .

50m

15m

**Note**

During buffer evaporation, worms can become "knotted" around each other. This effect is especially pronounced when bubbles are present. To disperse worm knots, it is best to sharply strike the plate against the bench. If necessary, you can also disperse worms by gently prodding them with a flame-sterilized inoculation loop.

9 Incubate plates on bench for  01:00:00 to allow worms to habituate.

1h

10 Repeat this process with additional strain(s), staggering the replating times to allow for consistent habituation time between strains.

Acquire imaging data

11 For each plate, image 25 unique FOVs (if sufficient worms are available) for 30 s each on an upright microscope. Use the widest field of view possible.

Image acquisition parameters:

- 30 s acquisition at 25 fps (40 ms exposure time)
- Plan Apo D 4×/0.20
- Nikon upright widefield microscope with Kinetix sCMOS camera