

Feb 13, 2023

Preparing Combined Indexed Primer Plates (IDT Ultramers) for the Illumina MiSeq - IDT UDIs

DOI

dx.doi.org/10.17504/protocols.io.6qpvr63wbvmk/v1

André M Comeau¹, Alessi Kwawukume¹

¹Integrated Microbiome Resource (IMR), Dalhousie University

Integrated Microbiome R...



André M Comeau

Integrated Microbiome Resource (IMR), Dalhousie University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr63wbvmk/v1>

Protocol Citation: André M Comeau, Alessi Kwawukume 2023. Preparing Combined Indexed Primer Plates (IDT Ultramers) for the Illumina MiSeq - IDT UDIs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr63wbvmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 07, 2022

Last Modified: February 13, 2023

Protocol Integer ID: 59158

Keywords: Illumina, primers, unique dual index, idt working primer stocks of illumina udi primer, illumina udi primer, idt working primer stock, idt ultramer, idt udis the preparation, preparing combined indexed primer plate, combined indexed primer plate, indexed primer plate, use in imr pcr prep, imr pcr prep, idt udi, primer, illumina miseq, pcr, illumina, preparation

Abstract

The preparation of diluted combined (F+R) IDT working primer stocks of Illumina UDI primers for use in IMR PCR preps.



Materials

The following materials list contains those consumables used specifically at the IMR to complete the present protocol.

- IDT stock primer DWP
- Eppendorf (or similar) DWPs
- Microplates sealing film Bio-Rad
- Tips ClipTip 20
- Tips ClipTip 300
- UltraPure water
- Reservoirs

Troubleshooting

Order Primers

- 1 Use our Excel template ( [Illumina-IDT-UDI-8bp-customfusio...](#)) to copy existing 16S/18S/ITS primers or to design your own custom gene primers with the proper Illumina indices and Nextera adapter orientations. We order **IDT "Ultramers"** for such long primers (~80-90 nt) as their coupling efficiency is one of the highest available (critical for obtaining high proportions of full-length oligos in the mix you obtain). Order the fusion primers at  4 nanomolar (nM) scale in deep-well plates (DWP); there will be 4 x 96-well plates, with **one unique Forward primer + one unique Reverse primer combined in each well**, arranged as follows:






F1R1 Plate (P1) = UDIs 1-96 (ie: i5_1+i7_1 in well A1, i5_2+i7_2 in well A2, etc.)

F1R2 Plate (P2) = UDIs 97-192

F2R1 Plate (P3) = UDIs 193-288





F2R2 Plate (P4) = UDIs 289-384

Prepare Archival Stocks

- 2 Once arrived, do a short spin of the plate in case lyophilized material was dislodged, then add  400 μL of PCR-grade water to each well of each plate containing the primers in order to reconstitute them at a concentration of  10 micromolar (μM) (1/10th the typical 100 μM working stock concentration for primers). Mix well by pipetting up and down at least 3 times and seal the plate with Bio-Rad film. Alternatively, the plate is sealed with Bio-Rad film and mixed well by vortexing it on a benchtop vortex for  00:00:30 and then doing a short spin at approx.  500 rpm, 00:00:30 . We have found that these primers usually need a significant incubation time for the lyophilized pellets to re-suspend well – we typically leave them overnight at  4 $^{\circ}\text{C}$ before continuing.

1m

Prepare Combined Working Stocks

- 3 Prepare the combined  1 micromolar (μM) working stock **F1R1 Primer Plate** by pipetting  228 μL of PCR-grade water into each well of an empty 96-well DWP from a sterile reservoir. Working by column and changing tips each time, transfer  12 μL of reconstituted **F1R1 stock** primer into each well of each column, mixing well by pipetting. Once complete, the resulting plate will have enough primer for 30 PCR plates (8 μL combined F+R per rxn \times 30 = 240 μL). Seal the plate with PCR film and store at  -20 $^{\circ}\text{C}$.



- 4 Prepare the combined [1M] 1 micromolar (μM) working stock **F1R2 Primer Plate** by repeating Step 3, but using the above reconstituted **F1R2 stock** deep-well primer plate.
- 5 Prepare the combined [1M] 1 micromolar (μM) working stock **F2R1 Primer Plate** by repeating Step 3, but using the above reconstituted **F2R1 stock** deep-well primer plate.
- 6 Prepare the combined [1M] 1 micromolar (μM) working stock **F2R2 Primer Plate** by repeating Step 3, but using the above reconstituted **F2R2 stock** deep-well primer plate.
- 7 Once all aliquoting is complete, seal the DWPs with PCR film and archive at $-20\text{ }^{\circ}\text{C}$ until new aliquots are required (minimized freeze-thaw cycles).

(Optional) Prepare Blocking Primer Stocks

- 8 **Optional:** For the generation of 18S V4 amplicons from microbiome samples containing substantial non-target host DNA (ex: human, mouse, etc.), order (ex: from PNA Bio) a custom PNA mammalian blocking primer (elongation arrest in the V4 region) with the sequence: 5'-**TCTTAATCATGGCCTCAGTT**-3' (courtesy of Laura Parfrey and Matt Lemay, UBC). Once arrived, prepare an archival stock of [1M] 100 micromolar (μM) and a working stock of [1M] 10 micromolar (μM) using PCR-grade water.