


Jan 10, 2022 Version 2

🌐 Preparing Gene of interest for GateWay cloning (2 step PCR process) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.b3muqk6w>



Johannes Wolfram JWD. Debler¹

¹Curtin University



Johannes Wolfram JWD Debler

Curtin University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b3muqk6w>

External link: https://www.embl.de/pepcore/pepcore_services/cloning/cloning_methods/gateway/2step/

Protocol Citation: Johannes Wolfram JWD. Debler 2022. Preparing Gene of interest for GateWay cloning (2 step PCR process). protocols.io <https://dx.doi.org/10.17504/protocols.io.b3muqk6w> Version created by [Johannes Wolfram JWD Debler](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 10, 2022

Last Modified: January 10, 2022

Protocol Integer ID: 56724

Keywords: Gateway, attB1, attB2, gateway recombination cloning, gene specific pcr with primer, utilizing universal gateway primer, universal gateway primer, gateway cloning, pcr product via primer, step gateway pcr experiment, preparing gene, step pcr process, specific pcr, second pcr, short extension of the first pcr product, gateway attb1, first pcr product, binding primer sequence, primer extension, pcr product, pcr, gene, primer sequence, attb2, attb2 site, gene of interest, primer, gateway attachment site, full attb1

Abstract

GateWay recombination cloning is achieved by flanking your gene of interest with GateWay attachment sites. In our case attB1 and attB2. Those sites are added to the PCR product via primers with 5' extensions. Since those primes create 31 bp and 30 bp 5' primer extensions respectively, plus about 20 bp of actual binding primer sequence it becomes expensive fast if you need 2 x ~50 bp primers for every GOI. We therefore use a 2 step PCR process to attach GateWay attB1 and attB2 sites. We first run a gene specific PCR with primers carrying short 5' extesions, and then a second PCR utilizing universal GateWay primers which bind to the short extension of the first PCR product to create the full attB1 and attB2 sites.

This protocol has been adapted from: **2-STEP GATEWAY PCR EXPERIMENTS**

Materials

MATERIALS

☒ Q5 High-Fidelity PCR Kit - 200 rxns **New England Biolabs Catalog #E0555L**

☒ SuperFi Polymerase **Thermo Fisher Scientific**

Use your favourite High Fidelity Polymerase or Mastermix. We have used Q5 and SuperFi, both work fine.

Primer design

- 1 There are a few things we need to keep in mind when designing the primers in the wider context of **Gateway cloning**. The idea is that you create an "entry clone" containing your gene of interest. Ususally **INCLUDING** the start codon, but **EXCLUDING** the stop codon. This is because one of the points of Gateway recombination is that you use the same entry clone to shuttle your gene of interest into different destination vectors with different properties (different N- or C- terminal tags). If you included the stop codon you would prevent yourself from being able to attach C-terminal tags.

Gene specific forward primer (PCR 1): 5'-AA AAA GCA GGC TNN -(15 to 20 bp template specific sequence)-3'

The 'NN' here can be any base. **Don't use AA, AG or GA though, as that would introduce an in-frame stop codon!** They are inserted to keep the reading frame if a destination vector with an N-terminal tag is used. I use 'CC' in all of my primer design.


Gene specific reverse primer (PCR 1): 5'-A GAA AGC TGG GTN -(15 to 20 bp template specific sequence)-3'

The 'N' here can be any base, just be careful to not accidentally introduce a stop codon. I use 'A' in all of my primer design.

attB1 adapter primer (PCR 2): 5'-GGG GAC AAG TTT GTA CAA AAA AGC AGG CT-3'

attB2 adapter primer (PCR 2): 5'-GGG GAC CAC TTT GTA CAA GAA AGC TGG GT-3'

This is what we are trying to achieve:



PCR 1 - gene specific PCR

2

Use your favourite High Fidelity Polymerase. I have used NEB Q5 and Thermo Fisher SuperFi Mastermix as listed below with great success.

10 ul total					
5x Q5 Buffer	2 ul		98° C	1 min	
2 mM dNTPs	1 ul		98° C	10 s	
10 uM F primer	0.5 ul		50° C	20 s	} 40 x
10 uM R primer	0.5 ul		72° C	30 s (per kb)	



	Q5 polymerase	0.1 ul		72°C	2 min	
	H2O	5.4 ul		10°C	hold	
	template	0.5 ul				

	10 ul total					
	2x SuperFi Mastermix	5		98°C	30 s	
	10 uM F primer	0.5 ul		98°C	10 s	
	10 uM R primer	0.5 ul		50°C	20 s	} 40 x
	H2O	3.5 ul		72°C	30 s (per kb)	
	template	0.5 ul		72°C	2 min	
				10°C	hold	

PCR 2 - GateWay attB1 and attB2 adapter primers

- 3 Use PCR 1 product as input for PCR 2.

	10 ul total			98°C	1 min	
--	--------------------	--	--	------	-------	--



	5x Q5 Buff er	2 ul		98° C	10 s	
	2 mM dNT Ps	1 ul		45° C	30 s	} 40 x
	Q5 poly mer ase	0.1 ul		72°C	30 s (per kb)	
	10 uM attB1 prim er	0.5 ul		72°C	2 min	
	10 uM attB 2 prim er	0.5 ul		10°C	hold	
	PCR 1 prod uct	1 ul				
	Cres ol red	1.7 ul				
	H2O	3.2 ul				

	10 ul total			98° C	1 min	
	2x Sup erFi Mast ermi x	5 ul		98° C	10 s	
	10 uM attB1	0.5 ul		45° C	20 s	} 40 x



primer					
10 uM attB 2 primer	0.5 ul		72°C	30 s (per kb)	
PCR 1 product	1 ul		72°C	2 min	
H2O	3 ul		10°C	hold	

Run a gel

- 4 Run PCR 2 product on a gel to make sure you have a crisp band (and to know if the PCR has worked).

Gel extraction

- 5 Extract the PCR 2 band with your gel extraction kit of choice.

Proceed to BP recombination step

- 6 Gateway BP recombination of attB tailed PCR products into pDONR/zeo