



Feb 13, 2023

🌐 Preparing Indexed Primer Plates (IDT Ultramers) for the Illumina MiSeq - Nextera Dual Indices

DOI

dx.doi.org/10.17504/protocols.io.3byl4bq1jvo5/v1

André M Comeau¹, Alessi Kwawukume¹

¹Integrated Microbiome Resource (IMR), Dalhousie University

Integrated Microbiome R...



André M Comeau

Integrated Microbiome Resource (IMR), Dalhousie University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4bq1jvo5/v1>

Protocol Citation: André M Comeau, Alessi Kwawukume 2023. Preparing Indexed Primer Plates (IDT Ultramers) for the Illumina MiSeq - Nextera Dual Indices. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4bq1jvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 06, 2022

Last Modified: February 13, 2023

Protocol Integer ID: 59117

Keywords: Illumina, primers, dual index, illumina dual index primer, idt working primer stock, diluted idt working primer stock, idt ultramer, preparing indexed primer plate, indexed primer plate, nextera dual indices the preparation, use in imr pcr prep, imr pcr prep, nextera dual indice, primer, illumina miseq, pcr, illumina, imr

Abstract

The preparation of diluted IDT working primer stocks of Illumina Dual Index primers for use in IMR PCR preps.



Materials

The following materials list contains those consumables used specifically at the IMR to complete the present protocol.

IDT stock primer DWP
PCR microplates 96-well Bio-Rad
PCR microplates sealing film Bio-Rad
Tips ClipTip 20
Tips ClipTip 200
Tips ClipTip 300
UltraPure water
Reservoirs






Troubleshooting

Order Primers

- 1 Use our Excel template ( [Illumina-CDI-8bp-customfusionpri...](#)) to copy existing 16S/18S/ITS primers or to design your own custom gene primers with the proper Illumina indices and Nextera adapter orientations. We order **IDT "Ultramers"** for such long primers (~80-90 nt) as their coupling efficiency is one of the highest available (critical for obtaining high proportions of full-length oligos in the mix you obtain). Order the fusion primers at  4 nanomolar (nM) scale in deep-well plates (DWP); one set per 96-well plate, arranged as follows, leaving blank rows in between sets:

S502	S503	S505	S506	S507	S508	S510	S511	←	Forward Set 1		
S513	S515	S516	S517	S518	S520	S521	S522	←	Forward Set 2		
N701	N702	N703	N704	N705	N706	N707	N710	N711	N712	N714	N715
								↖	Reverse Set 1		
N716	N718	N719	N720	N721	N722	N723	N724	N726	N727	N728	N729
								↖	Reverse Set 2		

Prepare Archival Stocks

- 2 Once arrived, do a short spin of the plate in case lyophilized material was dislodged, then add  400 µL of PCR-grade water to each well containing the primers in order to reconstitute them at a concentration of  10 micromolar (µM) (1/10th the typical 100 µM working stock concentration for primers). Mix well by pipetting up and down at least 3 times and seal the plate with Bio-Rad film. Alternatively, the plate is sealed with Bio-Rad film and mixed well by vortexing it on a benchtop vortex for  00:00:30 and then doing a short spin at approx.  500 rpm, 00:00:30 . We have found that these primers usually need a significant incubation time for the lyophilized pellets to re-suspend well – we typically leave them overnight at  4 °C before continuing.

1m

Prepare Working Stocks

- 3 Prepare the **1** micromolar (μM) working stock **Forward Set 1 Primer Plate (F1)** by pipetting **90** μL of PCR-grade water into each well of the 96-well PCR plate from a sterile reservoir. Rotate the deep-well primer plate 90° clockwise and align it so that the 8 occupied wells (= 8 different indices) of row 1 line up with the 8 rows of the new plate. Working by column and keeping the same set of tips, transfer **10** μL of reconstituted primer into each well of each column, mixing well by pipetting. Once complete, each column of the resulting plate will have enough primer for two complete 96-well plate PCRs (the Pure+1/10th plates; 12 columns \times 4 μL \times 2 plates = 96 μL required). Seal the plate with PCR film and store at **-20** °C .
- 4 Prepare the **1** micromolar (μM) working stock **Forward Set 2 Primer Plate (F2)** by repeating Step 3, but using row 3 of the reconstituted deep-well primer plate.
- 5 Prepare the **1** micromolar (μM) working stock **Reverse Set 1 Primer Plate (R1)** by pipetting **63** μL of PCR-grade water into each well of the 96-well PCR plate from a sterile reservoir. Align the deep-well primer plate horizontally (normal orientation) so that the 12 occupied wells (= 12 different indices) of row 5 line up with the 12 columns of the new plate. Working by row and keeping the same set of tips, transfer **7** μL of reconstituted primer into each well of each row, mixing well by pipetting. Once complete, each row of the resulting plate will have enough primer for two complete 96-well plate PCRs (the Pure+1/10th plates; 8 rows \times 4 μL \times 2 plates = 64 μL required). Seal the plate with PCR film and store at **-20** °C .
- 6 Prepare the **1** micromolar (μM) working stock **Reverse Set 2 Primer Plate (R2)** by repeating Step 5, but using row 7 of the reconstituted deep-well primer plate.
- 7 Once all aliquoting is complete, seal the deep-well plate with PCR film and archive at **-20** °C until new aliquots are required (minimized freeze-thaw cycles).

(Optional) Prepare Blocking Primer Stocks

- 8 **Optional:** For the generation of 18S V4 amplicons from microbiome samples containing substantial non-target host DNA (ex: human, mouse, etc.), order (ex: from PNA Bio) a custom PNA mammalian blocking primer (elongation arrest in the V4 region) with the sequence: 5'-**TCTTAATCATGGCCTCAGTT**-3' (courtesy of Laura Parfrey and Matt Lemay, UBC). Once arrived, prepare an archival stock of **100** micromolar (μM) and a working stock of **10** micromolar (μM) using PCR-grade water.

