

Aug 09, 2023

Preparing mitochondrial samples for immunoblot analysis

DOI

dx.doi.org/10.17504/protocols.io.n2bvj38enk5/v1

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis

WEHI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj38enk5/v1>

Protocol Citation: Louise Uoselis 2023. Preparing mitochondrial samples for immunoblot analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj38enk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 86215








Keywords: ASAPCRN, mitochondrial samples for immunoblot analysis protocol, mitochondrial samples for immunoblot analysis, preparation of mitochondrial sample, preparing mitochondrial sample, mitochondrial sample, immunoblot analysis, immunoblot analysis protocol

Abstract

Protocol for preparation of mitochondrial samples for immunoblot analysis.

Troubleshooting



- 1 Thaw mitochondrial stocks on ice, and aliquot out the desired amount of mitochondria.
- 2 Centrifuge each aliquot for  00:10:00  10000 x g, 4°C 10m
- 3 Carefully aspirate the supernatant from each sample.
- 4 Add a volume of 1x SDS sample buffer (5% w/v SDS, 10% v/v glycerol, 100 mM DTT, 50 mM Tris-Cl pH 6.8) equal to the amount of mitochondria (in ug) to each sample. Eg. If each sample contains  20 µg of mitochondria, add  20 µL of 1x SDS sample buffer.
- 5 Vortex samples for ~5 seconds to mix, and then boil at  99 °C with shaking at max speed for  00:10:00 . 10m
- 6 Allow samples to cool to room temperature, quickly centrifuge to collect liquid to the bottom of the tube, and vortex for ~3 seconds to ensure the samples are homogenous.
- 7 Samples can now be directly loaded onto an SDS-PAGE gel, or stored at  -20 °C for later use.