

Feb 09, 2023

Preparing multiplexed WGS/MetaG libraries with the Illumina DNA Prep kit for the Illumina NextSeq or MiSeq

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvob31bl4o/v1>

André Comeau¹, Gina Filloramo¹

¹Integrated Microbiome Resource (IMR), Dalhousie University

Integrated Microbiome R...



André M Comeau

Integrated Microbiome Resource (IMR), Dalhousie University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvob31bl4o/v1>

Protocol Citation: André Comeau, Gina Filloramo 2023. Preparing multiplexed WGS/MetaG libraries with the Illumina DNA Prep kit for the Illumina NextSeq or MiSeq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvob31bl4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 07, 2022

Last Modified: February 09, 2023

Protocol Integer ID: 59164

Keywords: illumina dna prep, libraries with the illumina dna prep, kit for the illumina nextseq, illumina nextseq, dna prep, illumina nextera flex, illumina protocol, illumina protocol 10000000254 16 v10, preparing multiplexed wg, metag, multiplexed wg, udi indices at the imr, dna, illumina, kit with cd, udi indice, kit

Abstract

The preparation of (meta)genomic libraries using the Illumina Nextera Flex (now "DNA Prep") kit with CD/UDI indices at the IMR.

Based upon Illumina protocol 10000000254 16 v10 (Aug.2021).



Materials

The following materials list contains those consumables+quantities used specifically at the IMR to complete the present protocol for 192 metagenomes (2 plates) on the NextSeq2000.

Pre-Library Sample Quantification

Quant-iT 1X HS dsDNA Kit (1000 samples)

Optical plate

PCR microplates 96-well Bio-Rad

PCR microplates sealing film Bio-Rad

Tips ClipTip 20

Tips ClipTip 200

Tips ClipTip 300

Reservoir

UltraPure water

Illumina Library Preparation

Nextera Flex DNA Sample Prep Kit (96 samples)

IDT for Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 Indexes, 96 Samples)

IDT for Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 Indexes, 96 Samples)

PCR microplates 96-well Bio-Rad

PCR microplates sealing film Bio-Rad

Manual p1000 tips

Tips ClipTip 20

Tips ClipTip 200

PCR tube strips (Axygen, 8 × 0.2 mL)

15 mL Falcon tubes

50 mL Falcon tubes

Ethanol

Reservoir

UltraPure water

Post-Library Quantification and Pooling

Quant-iT 1X HS dsDNA Kit (1000 samples)

Optical plate

PCR microplates sealing film Bio-Rad

Tips ClipTip 20

Tips ClipTip 300

PCR tube strips (Axygen, 8 × 0.2 mL)

1.5 mL Eppendorf tubes

Reservoir

Loading and Sequencing



PhiX Control Kit v3 (at least for 5 stocks)
NextSeq 2000 P3 Kit 300 cycle
Quant-iT 1X HS dsDNA Kit (1000 samples)
Qubit Tubes
1.5 mL Eppendorf tubes
Various manual pipette tips



Prepare and Quantify gDNA

- 1 Using two 96-well PCR plates, prepare  30 μL of 1:10 dilution using nuclease-free water (NFW) for up to 192 gDNA samples.
- 2 Quantify  5 μL of the **1:10 diluted gDNA** using the Invitrogen Quant-iT dsDNA HS assay kit on the fluorescent 96-well plate reader.

Note

Sample concentrations are required for calculating equal molar DNA ratios for library preparation input. Our standard practice is to add the appropriate volume of each sample so that  5 ng of gDNA are added to each sample well.

Perform Illumina DNA Prep Kit Library Preparation

7m

- 3 Follow the library preparation instructions as described for the **Illumina DNA Prep** kit (https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/illumina_prep/illumina-dna-prep-reference-guide-1000000025416-10.pdf) with the below modifications.

Note

We were originally doing up to 96 samples (1 plate) on the NextSeq550 and so using the **Nextera DNA CD Indexes (96 Indexes)** for indexing. However, we have now increased to up to 192 samples (2 plates) on the NextSeq2000 and so are now using the **IDT for Illumina DNA/RNA UD Indexes Set A and B, Tagmentation (96+96 Indexes)** for indexing.

- 3.1 **Amplify Tagmented DNA:** use 12 PCR cycles (based on 5 ng of input DNA).
- 3.2 **Final library clean-up:** wash beads with 80% ethanol as described, remove excess ethanol and allow beads to dry on the magnet for  00:05:00 . Remove the plates from the magnet and add  60 μL of Resuspension Buffer (RSB) to the beads. Pipette mix until thoroughly resuspended and incubate at room temperature for  00:02:00 .

7m



Place the plates back on the magnet until clear and transfer  55 μL of the supernatant into new 96-well plates.

Quantify and Pool Final Library

- 4 Quantify  5 μL of each **tagmented, amplified and purified** sample using the Invitrogen Quant-iT dsDNA HS assay kit on the fluorescent 96-well plate reader.
- 5 Pool **equal molar amounts** of up to 192 samples in a 1.5 mL microcentrifuge tube.
- 6 **Dilute** pooled library with RSB to the recommended loading concentration specific to the intended sequencer per Illumina's guidelines.