


Nov 07, 2018

Preparing primary T cells for fluorescence microscopy

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.vede3a6



Bulent Arman Aksoy¹, Pinar Aksoy¹, Jeff Hammerbacher¹

¹Medical University of South Carolina

Hammer Lab

Tech. support phone: +18437924527 email: arman@hammerlab.org



Bulent Arman Aksoy

Medical University of South Carolina

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.vede3a6

Protocol Citation: Bulent Arman Aksoy, Pinar Aksoy, Jeff Hammerbacher 2018. Preparing primary T cells for fluorescence microscopy. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.vede3a6>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 07, 2018

Last Modified: November 07, 2018

Protocol Integer ID: 17573



Materials

MATERIALS

⊗ PBS Invitrogen - Thermo Fisher

⊗ BD Cytotfix/Cytoperm™ Fixation/Permeabilization Solution Kit **Becton Dickinson (BD) Catalog #554714**

⊗ Pierce™ 16% Formaldehyde (w/v) Methanol-free **Thermo Fisher Scientific Catalog #28906**

⊗ Alexa Fluor™ 488 Phalloidin **Catalog #A12379**

⊗ ProLong™ Glass Antifade Mountant with NucBlue™ Stain **Thermo Fisher Scientific Catalog #P36983**



⊗ Cytology Funnels for Shandon CytoSpin™ Centrifuge Biomedical Polymers **VWR International (Avantor) Catalog #80094-254**





Cell Culture

- 1 Culture primary T cells using standard tissue culture techniques



Fixation, permeabilization, and staining

- 2 Collect 200K to 1 million cells and spin them at 300g for  00:05:00
- 3 Discard the supernatant, resuspend the cells in 1 mL PBS, spin them at 300g for  00:05:00 , and discard the supernatant.

- 4 Fix the cells:

If the cells are too clumpy and sticking to each other, first resuspend them in  750 μ L PBS and mix them thoroughly until there are no clumps and add  250 μ L of 16% formaldehyde on top.

If there are no cell clumps, simply resuspend cells in  1 mL of Cytofix solution.

- 5 Incubate the tubes at  4 °C for  00:30:00 on a nutating mixer, e.g.

Equipment

Fisherbrand™ Nutrating Mixers - Variable Speed

NAME

Mixer/Shaker

TYPE

Fisherbrand

BRAND





88-861-043

SKU

<https://www.fishersci.com/shop/products/nutating-mixers-variable-speed/p-6707018#?keyword=nutating+mixer>

LINK



- 6 Spin the cells at 300g for  00:05:00 , discard the supernatant, resuspend them in  1 mL of 1X CytoPerm Wash Buffer.
- 7 Optionally, add your stains, dyes (e.g. phalloidin stain) into the permeabilization buffer.
- 8 Cover the tubes with aliminum foil and incubate at  25 °C (Room temperature) for  00:30:00 or for the specific duration required for the stain on a nutating mixer, e.g.

Equipment

Fisherbrand™ Nutrating Mixers - Variable Speed

NAME

Mixer/Shaker

TYPE

Fisherbrand

BRAND





88-861-043

SKU




<https://www.fishersci.com/shop/products/nutating-mixers-variable-speed/p-6707018#?keyword=nutating+mixer>

LINK



- 9 Spin the cells at 300g for  00:05:00 , discard the supernatant, resuspend them in  1 mL of PBS
- 10 Spin the cells at 300g for  00:05:00 , discard the supernatant, resuspend them in  200 μ L of PBS.

Slide preparation

- 11 Cytospin  200 μ L of the stained cell suspension onto a glass slide or sample cover slip at 300-500g for  00:05:00
- 12 Mount the sample in between a glass slide and cover slip using one drop (at least 10 μ L) of ProLong Anti-fade reagent and keep the slide within a dark container for at least  12:00:00