



Mar 27, 2025

Prevention of Rhodanese Aggregation Assay

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkdbrxg5r/v1>

Patricia Yuste-Checa¹, Andreas Bracher¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkdbrxg5r/v1>

Protocol Citation: Patricia Yuste-Checa, Andreas Bracher, F Ulrich Hartl 2025. Prevention of Rhodanese Aggregation Assay. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlkdbrxg5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 20, 2025

Last Modified: March 27, 2025

Protocol Integer ID: 124797

Keywords: ASAPCRN, prevention of rhodanese aggregation assay, rhodanese aggregation assay, prevention of rhodanese aggregation, extracellular chaperone clusterin, rhodanese aggregation, clusterin, assay, prevention

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol is based on Weber F, Hayer-Hartl M. *Methods Mol Biol.* (2000) and details how to efficiently monitor the prevention of rhodanese aggregation by the extracellular chaperone Clusterin using a spectrophotometer.

Materials

Buffers and reagents:

- Rhodanese from bovine liver (Merck, R1756, MW 33,296 Da). Prepare a ~ 1M 200 micromolar (μ M) stock dilution in 1M 20 millimolar (mM) MOPS-NaOH pH 7.4 , 1M 50 millimolar (mM) NaCl, divide into small aliquots, flash-freeze in liquid nitrogen and store at -80°C .
- Denaturation buffer: 1M 6 Mass Percent guanidinium-HCl, 1M 5 millimolar (mM) DTT (in water or PBS, depending on the reaction buffer).
- PBS 1x
- Clusterin (see protocol "Clusterin purification from HEK293E cells", [dx.doi.org/10.17504/protocols.io.bvvkn64w](https://doi.org/10.17504/protocols.io.bvvkn64w)) or any other chaperone.



Prevention of Rhodanese Aggregation Assay

- 1 Thaw and lyophilize a rhodanese aliquot using a SpeedVac concentrator.
- 2 Resuspend lyophilized rhodanese with denaturation buffer to a final concentration of **1M 60 micromolar (μM)** (D-Rhod).

Note

For example, a **20 μL** aliquot at **8 μL** pellet in **80.8 μL** buffer. **1 μL** of this stock will be used for each reaction.

- 3 Incubate at **25 °C** for **01:00:00** . 1h
- 4 Make Clusterin stocks at **1M 30 micromolar (μM)** and **1M 90 micromolar (μM)** , for 1:1 and 1:3 D-Rhod:Clu molar ratios, respectively.
- 5
 - Prepare assay in a spectrophotometer microcuvette. Final volume **120 μL** .
 - Use a spectrophotometer with temperature control of the cuvette compartment.
 - Equilibrate temperature at **25 °C** .
- 5.1 Add **2 μL** Clusterin stock at **1M 30 micromolar (μM)** or **1M 90 micromolar (μM)** (final concentration 0.5 μM or 1.5 μM) to the bottom of the cuvette.

- 5.2 Add **117 μL** PBS buffer and mix gently with the pipette. Equilibrate at **25 °C** in the cuvette compartment.


Note

NOTE: Other reaction buffers can be used to test the prevention of aggregation. For example, to test the prevention of aggregation of rhodanese by Clusterin at **pH 5.2** use **1M 20 millimolar (mM)** Na-acetate **pH 5.2** , **1M 150 millimolar (mM)** NaCl, **1M 2 millimolar (mM)** CaCl₂ buffer.

- 5.3 Add  1 μL of denatured rhodanese stock at  60 micromolar (μM) (final concentration 0.5 μM) to the reaction mix.
- 5.4 Immediately mix by vortexing shortly, place the cuvette in the spectrophotometer, press autozero and start measurement at 320 nm wavelength every  00:00:02 .

Note

NOTE: Rhodanese rapidly aggregates once diluted in buffer. Clusterin at molar ratio D Rhod:Clu 1:3 almost completely prevents rhodanese aggregation.



Protocol references

Weber F, Hayer-Hartl M. Prevention of rhodanese aggregation by the chaperonin GroEL. *Methods Mol Biol.* 2000;140:111-5. doi: 10.1385/1-59259-061-6:111. PMID: 11484477.