


Aug 20, 2025

prime-seq 2

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn97kpl5d/v1>



Felix Pförtner¹, Eva Briem¹, Wolfgang Enard¹, Daniel Richter¹

¹Ludwig-Maximilians-Universität München

Felix Pförtner: The first two authors contributed equally to this work

Eva Briem: The first two authors contributed equally to this work



Daniel Richter

LMU

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn97kpl5d/v1>

Protocol Citation: Felix Pförtner, Eva Briem, Wolfgang Enard, Daniel Richter 2025. prime-seq 2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn97kpl5d/v1>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 02, 2024

Last Modified: August 20, 2025

Protocol Integer ID: 113399

Keywords: bulk RNA sequencing, RNA-seq, transcriptomics, SCRB-seq, prime-seq, efficient bulk rna, best bulk rna, seq protocol, increased final unique molecular identifier, sequencing cost, seq option, seq, final unique molecular identifier, optimized protocol, addressing dnase digestion, transcription, dnase digestion

Funders Acknowledgements:

DFG

Grant ID: 238187445

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This is the optimized version of the [prime-seq.protocol](#).

We systematically optimized the powerful and cost-efficient bulk RNA-seq protocol prime-seq. Notable efficiency enhancements were achieved by addressing DNase digestion, reverse transcription, adapter ligation, and amplification.

The optimized protocol increased final unique molecular identifiers (UMIs) by 60 % at equal sequencing costs or, conversely, reduced sequencing costs by 38 % at equal counts. This further improves one of the most cost-efficient bulk RNA-seq protocols available, making it currently one of the best bulk RNA-seq options.

Guidelines

- All reagents and plastic-ware can be found in the 'Materials' section.
- Use only RNase free supplies and clean all surfaces and tools with RNase Away prior to working
- Make sure all steps involving cell lysate and RNA before reverse transcription are carried out swiftly and on ice.
- All primer sequences are listed below:

A	B	C	D	E	F
Oligo	Vendor	Purification	Working Conc.	Sequence ("*" = PTO bond)	Notes
Barcoded Oligo-dT (E3V7NEXT)	Sigma	Cartridge	10 µM	ACACTCTTTCCCTACACGAC GCTCTTCCGATCT[12 bp BC]NNNNNNNNNNNNNNNN VTTTTTTTTTTTTTTTTTTT TTTTTTTTTTVN	
Template Switching Oligo (TSO) (E5V7NEXT)	Sigma	RNase-Free HPLC	100 µM	Biotin- ACACTCTTTCCCTACACGAC GCrGrGrG	
Preamp Primer (SINGV6 PTO)	Sigma	Standard Desalting	10 µM	Biotin- ACACTCTTTCCCTACACGAC *G*C	
3' enrichment primer (P5NEXTPT5)	Sigma	Standard Desalting	5 µM	AATGATACGGCGACCACCG AGATCTACACTCTTTCCCTA CACGACGCTCTTCCGATCT	
i7 Index Primer PTO (Nextera)	IDT	Trugrade	5 µM	CAAGCAGAAGACGGCATA GAGAT[i7]GTCTCGTGGGCT CGGAGATGTGTATAAGAGAC *A*G	
i5 Index Primer PTO (TruSeq)	IDT	Trugrade	5µM	AATGATACGGCGACCACCG AGATCTACAC[i5]ACACTCTT TCCCTACACGACGCTCTTC CGAT*C*T	
blocked prime-seq Adapter AntiSense	IDT	Standard Desalting	1.5 µM	Phos- CTGTCTCTTATACACATCTdd C	Duplexed DNA
blocked prime-seq Adapter Sense	IDT	Standard Desalting	1.5 µM	GTCTCGTGGGCTCGGAGAT GTGTATAAGAGACAGT	Duplexed DNA

Specific barcoded oligodT (E3V7NEXT) sequences:



E3V7_Set1.txt



E3V7_Set2.txt



Materials

MATERIALS

- ✕ DNase I Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0303S**
- ✕ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S**
- ✕ Deoxynucleotide Solution Mix - 40 umol of each **New England Biolabs Catalog #N0447L**
- ✕ Exonuclease I (E.coli) - 3,000 units **New England Biolabs Catalog #M0293S**
- ✕ Quant-it™ PicoGreen® dsDNA Assay Kit **Life Technologies Catalog #P7589**
- ✕ β -mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3148**
- ✕ QuantiFluor(R) RNA System **Promega Catalog #E3310**
- ✕ Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546**
- ✕ 5 M Sodium chloride (NaCl) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S5150-1L**
- ✕ Agilent High Sensitivity DNA Kit **Agilent Technologies Catalog #5067-4626**
- ✕ Buffer RLT Plus **Qiagen Catalog #1053393**
- ✕ Maxima H Minus Reverse Transcriptase (200 U/uL) **Thermo Fisher Scientific Catalog #EP0752**
- ✕ NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns **New England Biolabs Catalog #E6177S**
- ✕ EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889**
- ✕ Ethanol absolute **Carl Roth Catalog #9065.4**
- ✕ Igepal **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I8896**
- ✕ KAPA HiFi 2x RM **Kapa Biosystems Catalog #KR0370**
- ✕ Poly(ethylene glycol) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #89510**
- ✕ UltraPure DNase/RNase Free Distilled Water **Catalog #10977-049**
- ✕ Trizma hydrochloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2694**
- ✕ Aluminium seals for cold storage **Catalog #391-1275**
- ✕ Filter tips 96 low retention 10 uL **Catalog #771265**
- ✕ PCR Seals **Thermo Scientific Catalog #AB0558**
- ✕ twin.tec 96-well DNA LoBind Plates **Eppendorf Catalog #0030129504**
- ✕ Sera-Mag Speed Beads **GE Healthcare Catalog #65152105050250**
- ✕ Sodium Azide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S2002-100G**

Safety warnings

 Please follow all Manufacturer safety warnings and recommendations.

Before start

Wipe bench surfaces with RNase Away and keep working environment clean.



Preparation

12m

- 1 Clean all surfaces and pipettes with RNase Away 5m
- 2 Thaw frozen buffers and primers on ice 10m
- 3 Prepare 80% EtOH (approximately 45 mL for 96 samples) 2m
- 4 **When running the protocol for the first time prepare Cleanup Beads (see end of the protocol)!** 45m
- 5 prime-seq can be used on lysate or extracted RNA. It is essential, however, that the samples either have the same input or that they are normalized after the RNA is extracted, otherwise sequencing depth per sample will be impacted. Based on your starting material, please follow one of the following cases: 

STEP CASE

Lysate (similar input), Direct Lysis 107 steps

Follow this case if you are testing samples that have **similar input** (i.e. the expected RNA amount is the same between samples). The steps here will guide you in digesting residual proteins in your samples, extracting the RNA, digesting DNA, preparing RNA-seq libraries, and finally sequencing.

Example: investigating the genotype effect on transcription in 5,000 neurons

First Time Setup

- 6 When running the direct lysis protocol for the first time, prepare Bead Binding Buffer (see end of the protocol)!

Sample Collection

- 7 Prepare **Lysis Buffer** according to the number of samples. 2m

Reagent	Well	Plate
---------	------	-------



	RLT Plus Buffer	99 μ L	10.89 mL
	β -mercaptoethanol	1 μ L	110 μ L
	Total	100 μL	11 mL

Note

If sample volume exceeds 25 % of total lysate, use 2x TCL buffer (Qiagen, #1070498) + 1 % β -mercaptoethanol

8 Add  100 μ L of **Lysis Buffer** to each well of a semi-skirted 96-well PCR plate

1m

9 Add cells or tissue to wells

**Note****Cells**

Minimum: 100 cells, Optimum: 10,000 cells

Make sure that the same number of cells are used for each sample. Large differences between cells will impact distribution of sequencing reads and can potentially affect normalization.

**Note****Tissue**

If samples are difficult to lyse they should be homogenized using a tissue homogenizer.

Tissue should be a relatively small and not exceed more than 1000 ng of RNA. Tissue samples should be normalized by weight and be the same type of tissue.

Large differences between tissue samples will impact distribution of sequencing reads and can potentially affect normalization.

If you are unsure if the samples will contain the same amount of RNA, it is best to switch to the "Lysate (*variable*)" case in Step 13.

- 10 Transfer  50 μL of **lysate** to a new plate, return one plate immediately to -80°C freezer to save as a backup

1m

Note

Conversely, one can prepare two plates during sorting with 50 μL of lysis buffer.

Proteinase K Digest

30m

- 11 Prepare **Proteinase K Mix**

	A	B	C
Reagent	Well	Plate	
EDTA (50mM)	0.5 μL	55 μL	
H ₂ O	0.5 μL	55 μL	
Proteinase K (20mg/mL)	1 μL	110 μL	
Total	2 μL	220 μL	

- 12 Add  2 μL of **Proteinase K Mix** and mix by pipetting

1m



- 13 Incubate for 00:15:00 at 50 °C and then heat inactivate the Proteinase K for 00:10:00 at 75 °C

25m

Bead Clean Up

20m

- 14 Mix each bulk sample (50 µL per well) with 100 µL of **Cleanup Beads (22% PEG)**

1m

- 15 Incubate for 00:05:00 at 20 °C (Room Temp)

5m

Note

While binding, prepare **DNase I Mix** (Step 28)

- 16 Place on magnet stand until clear (~3 min) and then discard supernatant

3m

- 17 Wash with 100 µL of **80% EtOH** while the plate is on the magnet. Discard the supernatant

2m

Note

After adding EtOH, incubate for 30 s so that all beads are bound to magnet.

- 18 Repeat wash step once more

2m

19 Air dry beads for  00:03:00

3m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

DNase I Digest

1m

20 Add  10 μL H₂O and **resuspend beads by vortexing vigorously**

2m

Note

If you encounter bead clumping at this step, try to resuspend the beads by vigorously pipetting the samples. We have generated high quality prime-seq libraries despite heavy clumping.

21 Prepare **DNase I Mix**

3m

	A	B	C
	Reagent	Well	Plate
	DNase I	1 μL	110 μL
	DNase I Buffer (10x)	2 μL	220 μL
	H ₂ O	7 μL	770 μL
	Total	10 μL	1.1 mL

22 Add  10 μL of **DNase I Mix** and mix by pipetting

2m



- 23 Incubate DNase I Mix and beads for 00:10:00 at 20 °C (Room Temp) 10m
- 24 Heat inactivate the DNase I by adding 1 μL of **EDTA (50 mM)** and 20 μL of **Bead Binding Buffer (2X)** and incubating for 00:05:00 at 65 °C 6m
- 25 Place plate on magnet stand until clear (~3 min) and discard the supernatant. 3m
- 26 Wash with 100 μL of **80% EtOH** while the plate is on the magnet. Discard the supernatant 2m
- 27 Repeat wash step once more 2m
- 28 Air dry beads for 00:05:00 5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

Note

While drying, prepare **Reverse Transcription Mix**.

Reverse Transcription

5m

- 29 Prepare **Reverse Transcription Mix** 5m

	A	B	C
	Reagent	Well	Plate
	Maxima H Minus RT	0.15 µL	16.5 µL
	Maxima RT Buffer (5x)	2 µL	220 µL
	dNTPs (25 mM)	0.4 µL	44 µL
	TSO (E5V7NEXT) (100 µM)	0.1 µL	11 µL
	UltraPure Water	2.35 µL	258.5 µL
	Total	5 µL	550 µL

30 Add  4 µL **H2O** to the beads

1m

Note

The 4 µL of water can be combined with the Reverse Transcription Mix by increasing the water in Row 6 from 2.35 µL to 6.35 µL.

If working with many samples, or if using a stepper pipette or robot, we find that it is better to add some water separately to prevent the beads from drying too much.

31 Add  5 µL **Reverse Transcription Mix**

1m

32 Add  1 µL of **Barcoded oligodT (E3V7NEXT) (10 µM)** per well

2m

33 Incubate for  01:30:00 at  42 °C

1h 30m

cDNA Pooling & Purification

5m

34 Place the plate on a magnet

3m



35 Pool the supernatant of all wells into a 2 mL tube

10m

36 Add  10 μL of **Cleanup Beads (22% PEG)** *per sample* for a 1:1 ratio (e.g. 240 μL for 24 samples)

5m

Note

The EDTA in the **Cleanup Beads (22% PEG)** will inactivate the RT.

37 Incubate for  00:05:00 at  Room temperature to allow binding of the cDNA onto beads

5m

38 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m

39 Wash with  1 mL of **80% EtOH** while the tube is on the magnet, discard the supernatant

1m

Note

Volume of EtOH should be adjusted depending on the number of samples. More samples will require more EtOH to cover the beads completely.

40 Repeat wash step once more

1m

41 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.



42 Elute the beads in  17 μL of **UltraPure Water**

1m

43 Incubate for  00:05:00 at RT and transfer to a new PCR tube or plate

5m

Exonuclease I Treatment

35m

44 Add  2 μL of **ExoI Buffer (10x)** and  1 μL of **Exonuclease I**. Incubate as follows:

35m

	Step	Temperature	Time
	Incubation	37 C	20 min
	Heat Inactivation	80 C	10 min
	Storage	4 C	∞

45 Mix each sample (20 μL per well) with  16 μL of **Cleanup Beads (22% PEG)** for a 1:0.8 ratio

1m

46 Incubate for  00:05:00 at  Room temperature to allow binding of the cDNA onto beads

5m

47 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m

48 Wash with  50 μL of **80% EtOH** while the tube is on the magnet, discard the supernatant

1m

49 Repeat wash step once more

1m

50 Air dry beads for  00:05:00

5m

**Note**

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

51 Elute the beads in  20 μL of **UltraPure Water**

1m

52 Incubate for  00:05:00 at RT and transfer to a new PCR tube or plate

5m

Full length cDNA Amplification

1m

53 Prepare **Pre Amplification Mix**

1m

A	B
Reagent	1x
KAPA HiFi 2x RM	25 μL
Pre-amp Primer (SINGV6-PTO) (10 μM)	3 μL
UltraPure Water	2 μL
Total	30 μL

54 Add  30 μL **Pre Amplification Mix** to sample

1m

55 Incubate the Pre Amplification PCR as follows:

1h 30m

Step	Temperature	Time	Cycles
Initial Denaturation	98 C	3 min	1 cycle
Denaturation	98 C	15 sec	10 cycles*



Annealing	65 C	30 sec	
Elongation	72 C	4 min	
Final Elongation	72 C	10 min	1 cycle
Storage	4 C	∞	

Note

Adjust the number of cycles based on input (sample number, cell number, or concentration).

As a rule of thumb we assume big cells like embryonic stem cells to contain 10 pg of total RNA and small cells like T-cells ~ 1-2 pg

As a general guide we recommend:

Total RNA Input	Cycles
10 ng	16
50 ng	14
100 ng	12
500 ng	10
1000 ng	9

cDNA Bead Purification

1m

56 Mix sample with  40 μ L **Clean Up Beads (22% PEG)** for a ratio of 1:0.8

1m

57 Incubate for  00:05:00 at  20 °C (Room Temp)

5m

58 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m



59 Wash with  100 μL of **80% EtOH** while the tube is on the magnet, discard the supernatant

1m

60 Repeat wash step once more

1m

61 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

62 Elute cDNA in  10 μL **Buffer EB**

1m

63 Incubate for  00:05:00 at RT and transfer to a new PCR tube or plate

5m

Note

Stopping Point. Samples can be safely stored at  -20 °C and protocol can be continued at a later date.



cDNA Quantification and Quality Check

45m

64 Quantify the cDNA using the **QuantiFluor dsDNA System** or equivalent Qubit following the manufacturer's protocol. Use 1 μL of clean cDNA for quantification.

10m




65 Quality check the cDNA using the Agilent 2100 Bioanalyzer with **High Sensitivity DNA Analysis Kits**.

45m

Note

Passing the cDNA quality check does not guarantee that the data will be of high quality, however, if the cDNA fails the quality check it will usually not yield good libraries and will therefore generate lower quality data.

Expected result



Library Preparation

66



Note

Before starting, read the library preparation section carefully as there are a few steps that are very time sensitive.

67 Prepare **Fragmentation Mix**

1m

A	B
Reagent	1x
Ultra II FS Reaction Buffer	1.4 μL
Ultra II FS Enzyme Mix	0.4 μL
cDNA (4-8 ng/μL)	2.5 μL
Buffer TE	1.7 μL
Total	6 μL

Note

Ensure that the Ultra II FS Reaction Buffer is completely thawed. If a precipitate is seen in the buffer, pipette up and down several times to break it up, and quickly vortex to mix. Place on ice until use.

Note

Vortex the Ultra II FS Enzyme Mix for 5-8 seconds prior to use for optimal performance.

68 Vortex the **Fragmentation Mix** for  00:00:05 and immediately proceed to the next step

10s



69 Incubate the Fragmentation reaction as follows:

40m

	A	B	C
	Step	Temperature	Time
	Pre-Cool	4 C	∞
	Fragmentation	37 C	5 min*
	A Tailing and Phosphorylation	65 C	30 min
	Storage	4 C	∞

Note

*Fragmentation time needs to be adjusted for very low (<5ng)/ high (>50ng) cDNA amounts.

**Note**

Set heated lid to 75° C. Make sure the lid is at the correct temperature before you start the reaction.
Skip the first incubation step once you have added your samples.

Adapter Ligation

20m

70 Prepare **Adapter Ligation Mix**

1m

A	B
Reagent	1x
NEBNext Ultra II Ligation Master Mix	6 µL
NEBNext Ligation Enhancer	0.2 µL
blocked prime-seq Adapter (1.5 µM)	0.5 µL
Total	6.7 µL

71 Add  6.7 µL **Adapter Ligation Mix** to each replicate

1m

72 Incubate for  00:15:00 at  20 °C

15m

Note

Turn off heated lid

73 Add  37.3 µL **Buffer EB** to Samples

1m



74 Mix Sample with  26 μL **SPRI select beads**

1m

Note

We use SPRI Select Beads here instead of our home made 22% Clean Up beads for their guaranteed QCed size selection properties.

Note

The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.

75 Incubate for  00:05:00 at  20 °C (Room Temp)

5m

76 Place the plate on the magnet stand until clear and **transfer**  76 μL **supernatant to clean well.**

3m

Note

Be careful not to discard! This is your sample!

77 Mix supernatant with  10 μL **SPRI select beads**

1m

Note

The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.

78 Incubate for  00:05:00 at  20 °C (Room Temp)

5m

79 Place the plate on the magnet stand until clear and discard supernatant

3m



80 Wash with  150 μL of **80% EtOH** while the plate is on the magnet, discard the supernatant

1m

81 Repeat wash step once more

1m

82 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

83 Elute samples in  11 μL **0.1X TE** (dilute 1X TE Buffer 1:10 in water) for  00:05:00

5m

Library PCR

45m

84 Transfer  10.5 μL of samples to clean wells

85 Add  1 μL of **Index Primer (Nextera i7, 5 μM)** to each well

Note

This is the unique index that will be used for demultiplexing libraries.

86 Add  1 μL of **Index Primer (TruSeq i5, 5 μM)** to each well

Note

Alternatively the universal primer P5NEXTPT5 can be used in case the second index will not be sequenced.



87 Prepare **Library PCR Mix** by adding  12.5 μL of **Q5 Master Mix**

Note

Although scaled down, there will not be sufficient Q5 Master Mix (M0544L) in the kit. This item will have to be ordered separately.

88 Incubate the **Library PCR** reaction as follows:

A	B	C	D
Step	Temperature	Time	Cycles
Initial Denaturation	98 C	30 sec	1 cycle
Denaturation	98 C	10 sec	10 cycles*
Annealing/Elongation	72 C	1 min 15 sec	
Final Elongation	72 C	5 min	1 cycle
Storage	4 C	∞	

Note

* Adjust the number of cycles based on cDNA input.

As a general guide we recommend:

	cDNA Input	Cycles
	20 ng	10
	10 ng	11
	5 ng	12



Double Size Selection

25m

89 Add  25 μ L Buffer EB to Index PCR

90 Mix Index PCR with  26 μ L **SPRI select beads**

Note

We use SPRI Select Beads here instead of our home made 22% Clean Up beads for their guaranteed QCed size selection properties.

Note

The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.

91 Incubate for  00:05:00 at  20 °C (Room Temp)

92 Place the plate on the magnet stand until clear and **transfer**  76 μ L **supernatant to clean well.**

Note

Be careful not to discard! This is your library.

93 Mix supernatant with  10 μ L **SPRI select beads**

Note

The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.



94 Incubate for 00:05:00 at 20 °C (Room Temp)

95 Place the plate on the magnet stand until clear and discard supernatant.

96 Wash with 150 µL of **80% EtOH** while the plate is on the magnet, discard the supernatant

97 Repeat wash step once more

98 Air dry beads for 00:05:00

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

99 Elute in 15 µL **Buffer EB.**

100 Incubate for 00:05:00 and then place on magnet until clear. Transfer eluted library to new well.

Note

Stopping point. The libraries can be safely stored at -20 °C until they will be QCed and sequenced.

QC and quantification

45m


101 Quantify and quality control the library using the Agilent 2100 Bioanalyzer with **High Sensitivity DNA Analysis Kits**.

Note

Bulk libraries often yield high concentrations, which should be diluted to get accurate molarity measurements on the Bioanalyzer. Ideally, do not load more than 2 ng onto the chip.

Expected result

Libraries will typically exceed 1-5 ng/ μ l concentration



Sequencing

1m

- 102 Samples should be submitted according to your Sequencing Facility specifications. prime-seq is compatible with Illumina Sequencing.

At least 8 cycles are required for the Index Read (i7) and 28 cycles for the Read 1 (BC+UMI). Dual index sequencing can be done when using patterned flowcells. Read 2 (DNA) should be adjusted based on the quality of the genome being mapped to, but for human and mouse 50 cycles is sufficient.



Some potential sequencing options:

A	B	C	D	E	F
Sequencer	Read 1	Read 2	Index Read (i7)	Index Read (i5)	Kit
NovaSeq q	28	94	8	8	NovaSeq SP v1.5 100 cycle
NextSeq q 500/550	28	56	8	0	NextSeq 500/550 HiOut v3 75 cycle
NextSeq q 1000/2000	28	94	8	8	NextSeq 1000/2000 P2 100 cycle
NextSeq q 2000	28	52	8	0	NextSeq 2000 P3 50 cycles
HiSeq	28	114	8	0	HiSeq 3000/4000 150 cycles

Prepare Cleanup Beads (22% PEG)

10m

103 Prepare **PEG Solution (22%)** by adding all ingredients to a 50 mL falcon tube

10m

Reagent	Amount
PEG 8000	11 g
NaCl (5M)	10 mL
Tris-HCl (1M, pH 8.0)	500 µL
EDTA (0.5M)	100 µL
IGEPAL (10% solution)	50 µL
Sodium Azide (10% solution)	250 µL
UltraPure Water	up to 49 mL
Total	49 mL

**Note**

Do not add the total amount of water until after PEG is completely solubilized

- 104 Incubate at 40 °C and vortex regularly until PEG is completely dissolved 10m
- 105 Resuspend **Sera-Mag Speed Beads** carefully and pipette 1000 µL of bead suspension into a 1.5 mL tube 1m
- 106 Place on magnet stand and remove storage buffer 1m
- 107 Add 1000 µL of **TE Buffer** (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) and resuspend beads 30s
- 108 Place on magnet stand and remove supernatant 30s
- 109 Repeat wash step one more time 1m
- 110 Add 900 µL **TE Buffer** (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) and resuspend beads 30s
- 111 Add the washed **Sera-Mag Speed Beads** to the **PEG Solution (22%)** and mix well 1m 30s

Note

The final **Cleanup Beads (22% PEG)** can be aliquoted and stored at 4 °C for up to six months



Prepare Bead Binding Buffer

10m

112 Prepare **Bead Binding Buffer (2x)**

10m

Reagent	
PEG 8000	1.1 g
NaCl (5 M)	1 mL
Tris-HCl (1 M, pH 8.0)	50 µL
Igepal (10% solution)	5 µL
Sodium Azide (10% solution)	25 µL
H ₂ O	to 5 mL
Total	5 ml

Note

The **Bead Binding Buffer (2x)** can be stored at  Room temperature for up to six months.

Protocol references

Janjic, A., Wange, L.E., Bagnoli, J.W. *et al.* Prime-seq, efficient and powerful bulk RNA sequencing. *Genome Biol* **23**, 88 (2022). <https://doi.org/10.1186/s13059-022-02660-8>