

Oct 20, 2020

Version 1

prime-seq V.1

 [Scientific Reports](#)

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.s9veh66


Aleksandar Janjic¹, Lucas Esteban Wange¹, Johannes JWB Bagnoli¹, Johanna Geuder¹, Phong Nguyen¹, Daniel Richter¹, Christoph Ziegenhain², Wolfgang Enard¹

¹Ludwig-Maximilians-Universität München; ²Karolinska Institute Stockholm



Aleksandar Janjic

Ludwig-Maximilians-Universität München



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.s9veh66>

External link: <https://doi.org/10.1038/s41598-021-82883-0>

Protocol Citation: Aleksandar Janjic, Lucas Esteban Wange, Johannes JWB Bagnoli, Johanna Geuder, Phong Nguyen, Daniel Richter, Christoph Ziegenhain, Wolfgang Enard 2020. prime-seq. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.s9veh66>

Manuscript citation:

Janjic, A., Wange, L.E., Bagnoli, J.W. *et al.* Prime-seq, efficient and powerful bulk RNA sequencing. *Genome Bio* **23**,88 (2022).
<https://doi.org/10.1186/s13059-022-02660-8>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 05, 2018

Last Modified: October 20, 2020

Protocol Integer ID: 15381

Keywords: bulk RNA sequencing, RNA-seq, transcriptomics, SCRB-seq, sensitive single cell rna, open rna, seq method, cell rna, seq2, seq, wide range of model organism, single cell, model organism, cell

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract


prime-seq is a simple and open RNA-seq method, which can be easily established in most research labs or facilities. Based on two of the most sensitive single cell RNA-seq methods available, namely smart-seq2 and mcSCRB-seq (**Picelli et al. 2013; Bagnoli et al. 2018**), prime-seq is an adaptable, affordable, robust, and high-throughput option. Additionally, prime-seq can be used on a wide range of model organisms, from apes to zebrafish, and many in between.


Guidelines

- All reagents and plastic-ware can be found in the 'Materials' section.
- Use only RNase free supplies and clean all surfaces and tools with RNase Away prior to working
- Make sure all steps involving cell lysate and RNA before reverse transcription are carried out swiftly and on ice.
- All primer sequences are listed below:

	Oligo	Vend or	Purificati on	Worki ng Conc.	Sequence	Note s
	Barcoded Oligo-dT (E3V7NEXT)	Sigm a	Cartridge	10 µM	ACACTCTTTCCCTACACGAC GCTCTTCCGATCT[12 bp BC]NNNNNNNNNNNNNNNN VTTTTTTTTTTTTTTTTTTT TTTTTTTTTVN	
	Template Switching Oligo (TSO) (E5V7NEXT)	Sigm a	RNase-Free HPLC	100 µM	Biotin- ACACTCTTTCCCTACACGAC GCrGrGrG	
	Preamp Primer (SINGV6)	Sigm a	Standard Desalting	10 µM	Biotin- ACACTCTTTCCCTACACGAC GC	
	3' enrichment primer (P5NEXTPT5)	Sigm a	Standard Desalting	5 µM	AATGATACGGCGACCACCG AGATCTACACTCTTTCCCTA CACGACGCTCTTCCGATCT	
	i7 Index Primer (N7XX)	IDT	Trugrade	5 µM	CAAGCAGAAGACGGCATAC GAGAT[i7]GTCTCGTGGGCT CGG	
	prime-seq Adapter AntiSense	IDT	Standard Desalting	1.5 µM	/5Phos/CTGTCTCTTATACAC ATCT	Dupl exed DNA
	prime-seq Adapter Sense	IDT	Standard Desalting	1.5 µM	GTCTCGTGGGCTCGGAGAT GTGTATAAGAGACAGT	Dupl exed DNA

Specific barcoded oligodT (E3V7NEXT) sequences:

 E3V7_Set1.txt

 E3V7_Set2.txt



Materials

MATERIALS

- ✕ DNase I Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0303S**
- ✕ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S**
- ✕ Deoxynucleotide Solution Mix - 40 umol of each **New England Biolabs Catalog #N0447L**
- ✕ Exonuclease I (E.coli) - 3,000 units **New England Biolabs Catalog #M0293S**
- ✕ Quant-it™ PicoGreen® dsDNA Assay Kit **Life Technologies Catalog #P7589**
- ✕ β -mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3148**
- ✕ QuantiFluor(R) RNA System **Promega Catalog #E3310**
- ✕ Proteinase K solution, 20 mg ml – 1 **Ambion Catalog #AM2546**
- ✕ 5 M Sodium chloride (NaCl) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S5150-1L**
- ✕ Agilent High Sensitivity DNA Kit **Agilent Technologies Catalog #5067-4626**
- ✕ Buffer RLT Plus **Qiagen Catalog #1053393**
- ✕ Maxima H Minus Reverse Transcriptase (200 U/uL) **Thermo Fisher Scientific Catalog #EP0752**
- ✕ NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns **New England Biolabs Catalog #E6177S**
- ✕ EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889**
- ✕ Ethanol absolute **Carl Roth Catalog #9065.4**
- ✕ Igepal **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I8896**
- ✕ KAPA HiFi 2x RM **Kapa Biosystems Catalog #KR0370**
- ✕ Poly(ethylene glycol) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #89510**
- ✕ UltraPure DNase/RNase Free Distilled Water **Catalog #10977-049**
- ✕ Trizma hydrochloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2694**
- ✕ Aluminium seals for cold storage **Catalog #391-1275**
- ✕ Filter tips 96 low retention 10 uL **Catalog #771265**
- ✕ PCR Seals **Thermo Scientific Catalog #AB0558**
- ✕ twin.tec 96-well DNA LoBind Plates **Eppendorf Catalog #0030129504**
- ✕ Sera-Mag Speed Beads **GE Healthcare Catalog #65152105050250**
- ✕ Sodium Azide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S2002-100G**
- ✕ NEBNext Ultra II FS DNA Library Prep Kit for Illumina **New England Biolabs Catalog #E7805S**

Troubleshooting

Safety warnings

 Please follow all Manufacturer safety warnings and recommendations.

Before start

Wipe bench surfaces with RNase Away and keep working environment clean.



Preparation

12m

- 1 Clean all surfaces and pipettes with RNase Away.
- 2 Thaw frozen buffers and primers on ice.
- 3 Prepare 80% EtOH (approximately 45 mL for 96 samples)

5m

10m

2m

Prepare Cleanup Beads (22% PEG)

10m

- 4 Prepare **PEG Solution (22%)** by adding all ingredients to a 50 mL falcon tube

10m

	Reagent	Amount
	PEG 8000	11 g
	NaCl (5M)	10 mL
	Tris-HCl (1M, pH 8.0)	500 µL
	EDTA (0.5M)	100 µL
	IGEPAL (10% solution)	50 µL
	Sodium Azide (10% solution)	250 µL
	UltraPure Water	up to 49 mL
	Total	49 mL

Note

Do not add the total amount of water until after PEG is completely solubilized



- 5 Incubate at 40 °C and vortex regularly until PEG is completely dissolved 10m
- 6 Resuspend **Sera-Mag Speed Beads** carefully and pipette 1000 µL of bead suspension into a 1.5 mL tube 1m
- 7 Place on magnet stand and remove storage buffer 1m
- 8 Add 1000 µL of **TE Buffer** (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) and resuspend beads 30s
- 9 Place on magnet stand and remove supernatant 30s
- 10 Repeat wash step one more time 1m
- 11 Add 900 µL **TE Buffer** (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) and resuspend beads 30s
- 12 Add the washed **Sera-Mag Speed Beads** to the **PEG Solution (22%)** and mix well 1m 30s

Note

The final **Cleanup Beads (22% PEG)** can be aliquoted and stored at 4 °C for up to six months
- 13 prime-seq can be used on lysate or extracted RNA. It is essential, however, that the samples either have the same input or that they are normalized after the RNA is extracted, otherwise sequencing depth per sample will be impacted. Based on your starting material, please follow one of the following cases:

STEP CASE

Lysate (similar input) 96 steps

Follow this case if you are testing samples that have **similar input** (i.e. the expected RNA amount is the same between samples). The steps here will guide you in digesting residual proteins in your samples, extracting the RNA, digesting DNA,



preparing RNA-seq libraries, and finally sequencing.

Example: investigating the genotype effect on transcription in 5,000 neurons

Prepare Bead Binding Buffer


10m

14 Prepare **Bead Binding Buffer (2x)**

10m

Reagent	
PEG 8000	1.1 g
NaCl (5 M)	1 mL
Tris-HCl (1 M, pH 8.0)	50 µL
Igepal (10% solution)	5 µL
Sodium Azide (10% solution)	25 µL
H ₂ O	to 5 mL
Total	5 ml

Note

The **Bead Binding Buffer (2x)** can be stored at  Room temperature for up to six months.

Sample Collection

2m

15 Prepare **Lysis Buffer** according to the number of samples.

2m


Reagent	Well	Plate
RLT Plus Buffer	99 µL	10.89 mL
β-mercaptoethanol	1 µL	110 µL



	Total	100 μL	11 mL
--	--------------	------------------------------	--------------

Note

If sample volume exceeds 25 % of total lysate, use 2x TCL buffer (Qiagen, #1070498) + 1 % β -mercaptoethanol

16 Add  100 μ L of **Lysis Buffer** to each well of a semi-skirted 96-well PCR plate

1m

17 Add cells or tissue to wells

**Note****Cells**

Minimum: 100 cells, Optimum: 10,000 cells

Make sure that the same number of cells are used for each sample. Large differences between cells will impact distribution of sequencing reads and can potentially affect normalization.

Note**Tissue**


If samples are difficult to lyse they should be homogenized using a tissue homogenizer.

Tissue should be a relatively small and not exceed more than 1000 ng of RNA. Tissue samples should be normalized by weight and be the same type of tissue.

Large differences between tissue samples will impact distribution of sequencing reads and can potentially affect normalization.

If you are unsure if the samples will contain the same amount of RNA, it is best to switch to the "*Lysate (variable)*" case in Step 13.



- 18 Transfer  50 μL of **lysate** to a new plate, return one plate immediately to $-80\text{ }^{\circ}\text{C}$ freezer to save as a backup

1m

Note

Conversely, one can prepare two plates during sorting with 50 μL of lysis buffer.

Proteinase K Digest

30m

- 19 Add  1 μL **Proteinase K** (20 mg/mL) and  1 μL **EDTA** (25 mM) to each well

1m

- 20 Incubate for  00:15:00 at  $50\text{ }^{\circ}\text{C}$ and then heat inactivate the Proteinase K for  00:10:00 at  $75\text{ }^{\circ}\text{C}$



25m

Bead Clean Up

20m

- 21 Mix each bulk sample (50 μL per well) with  100 μL of **Cleanup Beads (22% PEG)**

1m

- 22 Incubate for  00:05:00 at  $20\text{ }^{\circ}\text{C}$ (Room Temp)

5m

Note

While binding, prepare **DNase I Mix** (Step 28)

- 23 Place on magnet stand until clear ($\sim 3\text{ min}$) and then discard supernatant

3m


- 24 Wash with  100 μL of **80% EtOH** while the plate is on the magnet. Discard the supernatant

2m



25 Repeat wash step once more

2m

26 Air dry beads for  00:03:00

3m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

DNase I Digest

1m

27 Add  5 μL H₂O and **resuspend beads by vortexing vigorously**

2m

Note









If you encounter bead clumping at this step, try to resuspend the beads by vigorously pipetting the samples. We have generated high quality prime-seq libraries despite heavy clumping.

28 Prepare **DNase I Mix**

3m

	Reagent	Well	Plate
	DNase I	1 μL	110 μL
	DNase I Buffer (10x)	2 μL	220 μL
	Bead Binding Buffer (2x)	10 μL	1.1 mL
	H ₂ O	2 μL	220 μL
	Total	15 μL	1.65 mL



- 29 Resuspend beads in  15 μL of **DNase I Mix** 2m
- 30 Incubate DNase I Mix and beads for  00:10:00 at  20 °C (Room Temp) 10m
- 31 Heat inactivate the DNase I by adding  1 μL of **EDTA (100 mM)** and incubating for  00:05:00 at  65 °C 6m
- 32 Place plate on magnet stand until clear (~3 min) and discard the supernatant. 3m
- 33 Wash with  100 μL of **80% EtOH** while the plate is on the magnet. Discard the supernatant 2m
- 34 Repeat wash step once more 2m
- 35 Air dry beads for  00:05:00 5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

Note

While drying, prepare **Reverse Transcription Mix** (Step 36).



Reverse Transcription

5m

36 Prepare **Reverse Transcription Mix**

5m

	Reagent	Well	Plate
	Maxima H Minus RT	0.15 μ L	16.5 μ L
	Maxima RT Buffer (5x)	2 μ L	220 μ L
	dNTPs (25 mM)	0.4 μ L	44 μ L
	TSO (E5V7NEXT) (100 μ M)	0.1 μ L	11 μ L
	UltraPure Water	2.35 μ L	258.5 μ L
	Total	5 μL	550 μL

37 Add 4 μ L **H₂O**

1m

Note

The 4 μ L of water can be combined with the Reverse Transcription Mix by increasing the water in Row 6 from 2.35 μ L to 6.35 μ L.

If working with many samples, or if using a stepper pipette or robot, we find that it is better to add some water separately to prevent the beads from drying too much.

38 Add 5 μ L **Reverse Transcription Mix**

1m

39 Add 1 μ L of **Barcoded oligodT (E3V7NEXT) (10 μ M)** per well

2m

40 Incubate for 01:30:00 at 42 $^{\circ}$ C

1h 30m

cDNA Pooling & Purification

5m



41 Place the plate on a magnet

3m

42 Pool all wells of one plate into a 2 mL tube



10m

43 Add  10 μL of **Cleanup Beads (22% PEG)** *for each sample* for a 1:1 ratio (e.g. 240 μL for 24 samples)

5m

Note


The EDTA in the **Cleanup Beads (22% PEG)** will inactivate the RT and make pooling easier due to the color.

44 Incubate for  00:05:00 at  Room temperature to allow binding of the cDNA onto beads

5m

45 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m

46 Wash with  1 mL of **80% EtOH** while the tube is on the magnet. Discard the supernatant


1m

Note

Volume of EtOH should be adjusted depending on the number of samples. More samples will require more EtOH to cover the beads completely.

47 Repeat wash step once more

1m

48 Air dry beads for  00:05:00

5m

**Note**

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

49 Elute the beads in 17 μL of **UltraPure Water**

1m

50 Incubate for 00:05:00 at RT and transfer to a new PCR tube or plate

5m

Exonuclease I Treatment

35m

51 Add 2 μL of **Exol Buffer (10x)** and 1 μL of **Exonuclease I**. Incubate as follows:

35m

Step	Temperature	Time
Incubation	37 C	20 min
Heat Inactivation	80 C	10 min
Storage	4 C	∞

52 Mix each sample (20 μL per well) with 16 μL of **Cleanup Beads (22% PEG)** for a 1:0.8 ration

1m

53 Incubate for 00:05:00 at Room temperature to allow binding of the cDNA onto beads

5m

54 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m


55 Wash with 50 μL of **80% EtOH** while the tube is on the magnet. Discard the supernatant

1m

56 Repeat wash step once more

1m



57 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

58 Elute the beads in  20 μL of **UltraPure Water**

1m

59 Incubate for  00:05:00 at RT and transfer to a new PCR tube or plate

5m

Full length cDNA Amplification

1m

60 Prepare **Pre Amplification Mix**

1m

Reagent	1x
KAPA HiFi 2x RM	25 μL
Pre-amp Primer (SINGV6) (10 μM)	3 μL
UltraPure Water	2 μL
Total	30 μL

61 Add  30 μL **Pre Amplification Mix** to sample

1m

62 Incubate the Pre Amplification PCR as follows:

1h 30m

Step	Temperature	Time	Cycles
------	-------------	------	--------



Initial Denaturation	98 C	3 min	1 cycle
Denaturation	98 C	15 sec	10 cycles*
Annealing	65 C	30 sec	
Elongation	72 C	4 min	
Final Elongation	72 C	10 min	1 cycle
Storage	4 C	∞	

Note

Adjust the number of cycles based on input (sample number, cell number, or concentration).

As a rule of thumb we assume big cells like embryonic stem cells to contain 10 pg of total RNA and small cells like T-cells ~ 1-2 pg

As a general guide we recommend:



Total RNA Input	Cycles
10 ng	16
50 ng	14
100 ng	12
500 ng	10
1000 ng	9

cDNA Bead Purification

1m

63 Mix sample with  40 µL **Clean Up Beads (22% PEG)** for a ratio of 1:0.8

1m

64 Incubate for  00:05:00 at  20 °C (Room Temp)

5m



65 Place the tube on the magnet stand until clear (~3 min) and discard supernatant

3m

66 Wash with  100 μ L of **80% EtOH** while the tube is on the magnet. Discard the supernatant

1m

67 Repeat wash step once more

1m

68 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.


69 Elute cDNA in  10 μ L **UltraPure Water**

1m

70 Incubate for  00:05:00 at RT and transfer to a new PCR tube or plate

5m

Note

Stopping Point. Samples can be safely stored at  -20 °C and protocol can be continued at a later date.

II

cDNA Quantification and Quality Check

45m

71 Quantify the cDNA using the **Quant-iT PicoGreen dsDNA assay kit** or equivalent Qubit following the manufacturer's protocol. Use 1 μ L of clean cDNA for quantification.

10m

*

72 Quality check the cDNA using the Agilent 2100 Bioanalyzer with **High Sensitivity DNA Analysis Kits**.


45m



Note

Passing the cDNA quality check does not guarantee that the data will be of high quality, however, if the cDNA fails the quality check it will usually not yield good libraries and will therefore generate lower quality data.

Expected result



Library Preparation

73



Note

Before starting, read the library preparation section carefully as there are a few steps that are very time sensitive.

74 Prepare **Fragmentation Mix**

1m

Reagent	1x
Ultra II FS Reaction Buffer	1.4 μ L
Ultra II FS Enzyme Mix	0.4 μ L
cDNA (4-8 ng/ μ L)	2.5 μ L



TE	1.7 μL
Total	6 μL

Note

Ensure that the Ultra II FS Reaction Buffer is completely thawed. If a precipitate is seen in the buffer, pipette up and down several times to break it up, and quickly vortex to mix. Place on ice until use.

Note

Vortex the Ultra II FS Enzyme Mix for 5-8 seconds prior to use for optimal performance.

75 Vortex the **Fragmentation Mix** for  00:00:05 and immediately proceed to step 70

10s

76 Incubate the Fragmentation reaction as follows:

40m

Step	Temperature	Time
Pre-Cool	4 C	∞
Fragmentation	37 C	5 min
A Tailing and Phosphorylation	65 C	30 min
Storage	4 C	∞

Note

Set heated lid to 75° C. Make sure the lid is at the correct temperature before you start the reaction.
Skip the first incubation step once you have added your samples.



Adapter Ligation

1m

77 Prepare **Adapter Ligation Mix**

1m

Reagent	1x
NEBNext Ultra II Ligation Master Mix	6 µL
NEBNext Ligation Enhancer	0.2 µL
prime-seq Adapter (1.5 µM)	0.5 µL
Total	6.7 µL

78 Add 6.7 µL **Adapter Ligation Mix** to each replicate

1m

79 Incubate for 00:15:00 at 20 °C

15m

Note

Turn off heated lid

80 Add 37.3 µL Buffer EB to Samples

1m

81 Mix Index PCR with 26 µL **SPRI select beads** (ratio of 1:0.52)

1m

Note

We use SPRI Select Beads here instead of our home made 22% Clean Up beads for their guaranteed QCed size selection properties.

82 Incubate for 00:05:00 at 20 °C (Room Temp)

5m



83 Place the plate on the magnet stand until clear and **transfer 76 μ L supernatant to clean well.**



3m

Note

Be careful not to discard! This is your sample!

84 Mix supernatant with  10 μ L **SPRI select beads** (ratio of 1:0.72)


1m

85 Incubate for  00:05:00 at  20 °C (Room Temp)

5m

86 Place the plate on the magnet stand until clear and discard supernatant.

3m

87 Wash with  150 μ L of **80% EtOH** while the plate is on the magnet. Discard the supernatant

1m

88 Repeat wash step once more

1m

89 Air dry beads for  00:05:00

5m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

90 Elute samples in  10.5 μ L **0.1X TE** (dilute 1X TE Buffer 1:10 in water) for  00:05:00

5m

Library PCR

1m

91 Transfer samples to clean wells

1m



92 Add  1 μL of **Index Primer (i7, 5 μM)** to each well

1m

Note

This is the unique index that will be used for demultiplexing libraries.

93 Prepare **Library PCR Mix**

1m

	Reagent	1x	5x
	NEBNext Ultra II Q5 Master Mix	12.5 μL	62.5 μL
	P5NEXTPT5 primer (5 μM stock)	1 μL	5 μL
	Total	13.5 μL	67.5 μL

Note

Although scaled down, there will not be sufficient Q5 Master Mix (M0544L) in the kit. This item will have to be ordered separately.

Note

Instead of the universal P5NEXTPT5 primer, one can use unique i5 Indices for each sample. This is recommended if using a sequencer with patterned flow cells to account for index hopping.

94 Add  13.5 μL of **Library PCR Mix** to each well

1m

95 Incubate the **Library PCR** reaction as follows:

20m



	Step	Temperature	Time	Cycles
	Initial Denaturation	98 C	30 sec	1 cycle
	Denaturation	98 C	10 sec	10 cycles*
	Annealing/Elongation	65 C	1 min 15 sec	
	Final Elongation	65 C	5 min	1 cycle
	Storage	4 C	∞	

Note


Adjust the number of cycles based on total cDNA input.

As a general guide we recommend:


	cDNA Input	Cycles
	20 ng	10
	10 ng	11
	5 ng	12

Library Double Size Selection

30s

96 Add  25 µL Buffer EB to Index PCR

1m

97 Mix Index PCR with  26 µL **SPRI select beads** (ratio of 1:0.52)

1m

**Note**

We use SPRI Select Beads here instead of our home made 22% Clean Up beads for their guaranteed QCed size selection properties.

98 Incubate for 00:05:00 at 20 °C (Room Temp) 5m

99 Place the plate on the magnet stand until clear and **transfer** 76 μ L **supernatant to clean well.** 30s

Note

Be careful not to discard! This is your library.

100 Mix supernatant with 10 μ L **SPRI select beads** (ratio of 1:0.72) 1m

101 Incubate for 00:05:00 at 20 °C (Room Temp) 5m

102 Place the plate on the magnet stand until clear and discard supernatant. 3m

103 Wash with 150 μ L of **80% EtOH** while the plate is on the magnet. Discard the supernatant 1m


104 Repeat wash step once more 1m

105 Air dry beads for 00:05:00 5m


Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.



106 Elute in  15 μ L **UltraPure Water.**


1m

107 Incubate for  00:05:00 and then place on magnet until clear. Transfer eluted library to new well.

5m



Note

Stopping point. The libraries can be safely stored at  -20 °C until they will be QCed and sequenced.

Library QC

45m

108 Quantify and quality control the library using the Agilent 2100 Bioanalyzer with **High Sensitivity DNA Analysis Kits.**


45m

Note

Bulk libraries often yield high concentrations, which should be diluted to get accurate molarity measurements on the Bioanalyzer. Ideally, do not load more than 2 ng onto the chip.

Expected result

Successful libraries will typically be 2-5 ng/μl concentration



Sequencing

109 Samples should be submitted according to your Sequencing Facility specifications. At least 8 cycles are required for the Index Read (i7) and 28 cycles for the Read 1 (BC+UMI). Read 2 (DNA) should be adjusted based on the quality of the genome being mapped to, but for human and mouse 50 cycles is sufficient.

Some potential sequencing options:

	Read	Read 1	Read 2	Index Read (i7)	Index Read (i5)	Kit
	HiSeq	28	114	8	0	HiSeq 3000/4000 150 cycles
	NextSeq	28	56	8	0	NextSeq 500/

						550 HiO ut v3 75 cycl e
	NovaSeq	28	94	8	8	NovaSeq SP v1.5 100 cycl e