

Oct 18, 2018 Version 4

Primer ID MiSeq Library Prep V.4

DOI


dx.doi.org/10.17504/protocols.io.useewbe

Shuntai Zhou¹

¹University of North Carolina at Chapel Hill



Shuntai Zhou



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.useewbe>

Protocol Citation: Shuntai Zhou 2018. Primer ID MiSeq Library Prep. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.useewbe>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 18, 2018

Last Modified: October 18, 2018

Protocol Integer ID: 16934

Keywords: NGS, Primer ID, qiaamp viral rna extraction kit, primer id miseq library prep, viral rna extraction kit, cdna primer, primer id miseq, viral rna, sequencing library, cdna synthesis, primer id, kapa hifi pcr kit, rna, random nucleotide, pcr, using qiaamp, cdna, block of random nucleotide, miseq, rounds of pcr

Abstract

This is the protocol to prepare Primer ID MiSeq sequencing library. Viral RNA was first extracted using QIAamp viral RNA extraction kit. The block of random nucleotides (Ns) in the cDNA primers served as the Primer ID. The Superscript III kit was used for the cDNA synthesis. We used two rounds of PCR to amplify the cDNA and incorporate Illumina indexed adapters with KAPA2G Robust and KAPA HiFi PCR kits, respectively.

Attachments




[amplicon protocol_v1...](#)

222KB

Guidelines

Viral RNA extraction using QIAamp vRNA mini kit.



Primers (NOTE: HIV-1 ENV REGION AS AN EXAMPLE, SUBSTITUTE 1 in the forward primer and cDNA primer)



V1F(forward) 5'-3'	GCCTCCCTCGCGCCATCAGAGATGTGTATAAGAGACAG AAAGCCTAAAGCCATGTGTA
BV3R <u>Uni</u> (cDNA Primer) 5'-3'	GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTNNNN TGCTCTACTAATGTTACAATGTGC
Universal Adapter	AATGATACGGCGACCACCGAGATCTACACGCCTCCCTC G
Indexed Adapter	CAAGCAGAAGACGGCATACGAGATNNNNNNGTGACT GCTC
ADPT_2a	GTGACTGGAGTTCAGACGTGTGCTC

Note: Primer ID primer and forward primer use random bases. Indexed primers have 24 fixed barcodes.

Table of Indexed Primers



Indexed Primer	Index	Index Sequence	Sequence
PCR Primer, Index 1	1	ATCACGA	CAAGCAGAAGACGGCATAACGAGAT CGTGAT GTG
PCR Primer, Index 2	2	CGATGTA	CAAGCAGAAGACGGCATAACGAGAT ACATCG GTG
PCR Primer, Index 3	3	TTAGGCA	CAAGCAGAAGACGGCATAACGAGAT GCCTAA GTG
PCR Primer, Index 4	4	TGACCAA	CAAGCAGAAGACGGCATAACGAGAT TGGTCA GTG
PCR Primer, Index 5	5	ACAGTGA	CAAGCAGAAGACGGCATAACGAGAT CACTGT GTG
PCR Primer, Index 6	6	GCCAATA	CAAGCAGAAGACGGCATAACGAGAT ATTGGC GTG
PCR Primer, Index 7	7	CAGATCA	CAAGCAGAAGACGGCATAACGAGAT GATCTG GTG
PCR Primer, Index 8	8	ACTTGAA	CAAGCAGAAGACGGCATAACGAGAT TCAAGT GTG
PCR Primer, Index 9	9	GATCAGA	CAAGCAGAAGACGGCATAACGAGAT CTGATC GTG
PCR Primer, Index 10	10	TAGCTTA	CAAGCAGAAGACGGCATAACGAGAT AAGCTA GTG
PCR Primer, Index 11	11	GGCTACA	CAAGCAGAAGACGGCATAACGAGAT GTAGCC GTG
PCR Primer, Index 12	12	CTTGTA	CAAGCAGAAGACGGCATAACGAGAT TACAAG GTG
PCR Primer, Index 13	13	TCCATAA	CAAGCAGAAGACGGCATAACGAGAT TATGGA GTG
PCR Primer, Index 14	14	GTACTAA	CAAGCAGAAGACGGCATAACGAGAT TAGTAC GTG
PCR Primer, Index 15	15	ACAGTAA	CAAGCAGAAGACGGCATAACGAGAT TACTGT GTG
PCR Primer, Index 16	16	CTCATGA	CAAGCAGAAGACGGCATAACGAGAT CATGAG GTG
PCR Primer, Index 17	17	ACGATAA	CAAGCAGAAGACGGCATAACGAGAT TATCGT GTG
PCR Primer, Index 18	18	TGCAGAA	CAAGCAGAAGACGGCATAACGAGAT TCTGCA GTG
PCR Primer, Index 19	19	TTCATAA	CAAGCAGAAGACGGCATAACGAGAT TATGAA GTG
PCR Primer, Index 20	20	TGCTGTA	CAAGCAGAAGACGGCATAACGAGAT ACAGCA GTG
PCR Primer, Index 21	21	TATCACA	CAAGCAGAAGACGGCATAACGAGAT GTGATA GTG
PCR Primer, Index 22	22	TGGATAA	CAAGCAGAAGACGGCATAACGAGAT TATCCA GTG
PCR Primer, Index 23	23	CGCATTA	CAAGCAGAAGACGGCATAACGAGAT AATGCG GTG
PCR Primer, Index 24	24	GCCTTAA	CAAGCAGAAGACGGCATAACGAGAT TAAGGC GTG

Materials

MATERIALS

☒ QIAamp vRNA mini kit **Qiagen Catalog #51304**

☒ Agencourt RNAClean XP **Beckman Coulter Catalog #A63987**

☒ 70% ethanol

☒ DNase-free water

☒ AmpureXP PCR cleanup kits **Beckman Coulter Catalog #A63880**

☒ QIAquick gel extraction kit **Qiagen Catalog #28704**

☒ Qubit dsDNA BR Assay kit **Invitrogen - Thermo Fisher Catalog #Q32850**

Troubleshooting



Prepare Primer Mix (Optional, only for multiplexed Primer ID library prep)

1 For multiplexing sequencing, first, prepare Primer Mix.

Example (For HIV drug resistance pipeline).

	Regions	DR cDNA primer	DR F primer
	PR	R2614_PID	F2163AD
	RT	R3284_PID11	F2620_AD
	IN	R4752_PID11	F4383_AD
	V3	R7209_PID11	V1F_AD

Make 10 μ M primer mix: mix 10 μ L of each primer in one set and 60 μ L of dH₂O.

Primer Tables

primer	sequence (5'-3')
R2614_PID	GTG ACT GGA GTT CAG ACG TGT GCT CTT CCG ATC TNN NNN NNN NCA GTT TAA CTT TTG



		GGC CAT CCA TTC C
	R32 84_ PID1 1	GTG ACT GGA GTT CAG ACG TGT GCT CTT CCG ATC TNN NNN NNN NNN CAG TCA CTA TAG GCT GTA CTG TCC ATT TAT C
	R47 52_P ID11	GTG ACT GGA GTT CAG ACG TGT GCT CTT CCG ATC TNN NNN NNN NNN ATC GAA TAC TGC CAT TTG TAC TGC
	R72 09_P ID11	GTG ACT GGA GTT CAG ACG



		TGT GCT CTT CCG ATC TNN NNN NNN NNN CAG TCC ATT TTG CTY TAY TRA BVT TAC AAT RTG C
	F216 3AD	GCC TCC CTC GCG CCA TCA GAG ATG TGT ATA AGA GAC AGN NNN TCA GAG CAG ACC AGA GCC AAC AGC CCC A
	F26 20_A D	GCC TCC CTC GCG CCA TCA GAG ATG TGT ATA AGA GAC AGN NNN GGC CAT TGA



		CAG AAG AAA AAA TAA AAG C
	F43 83_ AD	GCC TCC CTC GCG CCA TCA GAG ATG TGT ATA AGA GAC AGN NNN AAA AGG AGA AGC CAT GCA TG
	V1F_ AD	GCC TCC CTC GCG CCA TCA GAG ATG TGT ATA AGA GAC AGN NNN TTA TGG GAT CAA AGC CTA AAG CCA TGT GTA

"N" in this Primer table is a random nucleotide.

	P1	AAT GAT ACG
--	----	-------------------



		GCG ACC ACC GAG ATC TAC ACG CCT CCC TCG CGC CAT CAG AGA TGT G
	Indexed Adapter	CAA GCA GAA GAC GGC ATA CGA GAT NNN NNN GTG ACT GGA GTT CAG ACG TGT GCT C
	ADP T ₂ a	GTG ACT GGA GTT CAG ACG TGT GCT C

The 6 "N"s are not random nucleotides. They are a set of 24 pre-designed indexed sequences.

2 Mix well.

cDNA synthesis

3 Pipette the following components into a 0.5 ml RNase-free tube:

	μl/tube		[stock]	[final]	[master mix]
	3.0	dNTP Mix	10 mM each	0.5	
	1.5	cDNA primer	10 μM	0.25 μM	
	34.5	RNA template			
	39.0	Total volume			

🧪 3 μL dNTP Mix

🧪 1.5 μL cDNA primer

🧪 34.5 μL RNA template

- 4 Place tube in 65°C heat block for 3-5'.

🕒 00:05:00 65°C heat block

- 5 Place the tube on ice for 1'.

🕒 00:01:00 on ice

- 6 Add the following components:

	μl/tube		[stock]	[final]	[master mix]
	12.0	5x buffer	5x	1x	
	3.0	DTT	100 mM	5	
	3.0	RNaseOUT	40 u/μl	2	
	3.0	SSIII RT	200 u/μl	10	
	21.0	Per tube			



🧪 12 μ L 5x buffer

🧪 3 μ L DTT

🧪 3 μ L RNaseOUT

🧪 3 μ L SSII RT

7 Mix and incubate at 50°C for 1 hr.

🕒 01:00:00 Incubation at 50°C

8 Increase to 55°C and incubate for 1 hr.

🕒 01:00:00 Incubation at 55°C

9 Inactivate SSIII RT by heating at 70°C for 15'.

🕒 00:15:00 Heating at 70°C

10 To each tube, add 1 μ L RNase H, incubate at 37°C for 20'.

🧪 1 μ L RNase H

🕒 00:20:00 Incubation at 37°C

Purification

11 Purify cDNA using Agencourt RNAClean XP.
Resuspend the beads and take an aliquot out. Keep at room temperature for at least **30 minutes** before use. (Should be in 1ml aliquots)

🕒 00:30:00 Room temperature

12 Transfer the cDNA reactions into 1.7 mL RNase-free tubes.

13 Resuspend the beads (Vortex). Add **42 μ L of beads to 60 μ L cDNA (Ratio: 0.6 – 0.8)**
Agencourt RNAClean XP beads to each cDNA reaction.

🧪 42 μ L beads

🧪 60 μ L cDNA

14 Mix the Agencourt RNAClean XP and sample thoroughly by pipette mixing 15 times. No **vortexing**. Let the tube incubate at room temperature for **20 minutes** before proceeding to the next step.

🕒 00:20:00 Room temperature

15 Place the tube onto the magnetic tube rack for **5 minutes** to separate the beads from solution.



00:05:00 Magnetic tube rack

- 16 Slowly aspirate the cleared solution from the tube and discard. This step should be performed while the tube is situated on the rack. Do not disturb the magnetic beads, which have formed a spot on the side of the tube.
- 17 Dispense 500 μ L of 70% ethanol into the tube and incubate for **30 seconds** at room temperature. Aspirate out the ethanol and discard. (1/3)

It is important to perform this step with the tube situated **on the rack. Do not disturb the separated magnetic beads.** Be sure to remove all of the ethanol from the bottom of the well as it may contain residual contaminants.

00:00:30 Incubation at room temperature

- 18 Dispense 500 μ L of 70% ethanol into the tube and incubate for **30 seconds** at room temperature. Aspirate out the ethanol and discard. (2/3)

It is important to perform this step with the tube situated **on the rack. Do not disturb the separated magnetic beads.** Be sure to remove all of the ethanol from the bottom of the well as it may contain residual contaminants.

00:00:30 Incubation at room temperature

- 19 Let the reaction tube air-dry **10 minutes** on the rack with the cap open. The tube(s) should air-dry until the last visible traces of ethanol evaporate. Over drying the sample may result in a lower recovery.

00:10:00 Air-drying

- 20 **(Optional for clinical samples)** Elute in 60 μ L water. Repeat step 12 to 19.

undefined

- 21 Remove the tube from the rack and resuspend beads in 24 μ L DNase-free water by pipetting up and down. Place tube back on the rack and leave for **3 minutes**.

 24 μ L DNase-free water

00:03:00 Magnetic rack

- 22 Pipette the eluant from the tube while it is situated **on the magnetic tube rack.**

PCR 1

- 23 Complete thaw and vortex KAPA reagents (except for enzyme) before use.

24

	μl/tube		[stock]	[final]	[mastermix]
	10.0	5x Buffer A	5x	1x	
	10.0	Enhancer	5x	1x	
	1.0	dNTPs	10 mM	0.2 mM	
	2.5	Forward primer	10 μM	0.5 μM	
	2.5	ADP T ₂ a	10 μM	0.5 μM	
	0.5	KAP A Robust polymerase	5 U/μl	2.5 U	
	23.5	Template cDNA			

Prepare mastermix in cold box and use repeater pipette to add to each tube in the cold box. Add template cDNA to each tube and pipette up and down to mix.

🧪 10 μL 5x Buffer A

🧪 10 μL Enhancer

🧪 1 μL dNTPs

🧪 2.5 μL Forward primer

🧪 2.5 μL ADPT₂a

🧪 0.5 μL KAPA Robust polymerase

🧪 23.5 μL Template cDNA





25

Cycle

	95°C	1 min

	95°C	15 s
	58°C	1 min
	72°C	30 s
	15–25 cycles	
	72°C	3 min
	4°C	On hold

Purification PCR products

- 26 Purify PCR products using AmpureXP PCR cleanup kits. Vortex the 1 ml aliquot and remove the needed volume. Keep at room temperature for at least **30 minutes** before use.
-  00:30:00 Room temperature
- 27 Transfer the PCR1 reactions into 1.7 mL RNase-free tubes.
- 28 Resuspend the beads. Add **40 µl (Ratio: 0.6 – 0.8: 1, 36µl – 48µl)** Ampure XP beads to each cDNA reaction.
-  40 µL Ampure XP beads
- 29 Mix the Ampure XP and sample thoroughly by **vortexing**. Let the tube incubate at room temperature for **5 minutes** before proceeding to the next step (incube off the rack).
-  00:05:00 Incubation at room temperature
- 30 Place the tube onto the magnetic tube rack for **5 minutes** to separate the beads from solution.
-  00:05:00 Incubation magnetic rack
- 31 Slowly aspirate the cleared solution from the tube and discard. This step should be performed while the tube is situated on the rack. Do not disturb the magnetic beads, which have formed a spot on the side of the tube.
- 32 Dispense 500 µL of 70% ethanol into the tube and incubate for **30 seconds** at room temperature. Aspirate out the ethanol and discard. (1/2)

**Note**

It is important to perform this step with the tube situated **on the rack**. **Do not disturb the separated magnetic beads**. Be sure to remove all of the ethanol from the bottom of the well as it may contain residual contaminants.

🧪 500 µL 70% ethanol

⌚ 00:00:30 Incubation at room temperature

- 33 Dispense 500 µL of 70% ethanol into the tube and incubate for **30 seconds** at room temperature. Aspirate out the ethanol and discard. (2/2)

Note

It is important to perform this step with the tube situated **on the rack**. **Do not disturb the separated magnetic beads**. Be sure to remove all of the ethanol from the bottom of the well as it may contain residual contaminants.

🧪 500 µL 70% ethanol

⌚ 00:00:30 Incubation at room temperature

- 34 Let the reaction tube air-dry **10 minutes** on the rack with the cap open. The tube(s) should air-dry until the last visible traces of ethanol evaporate. Over drying the sample may result in a lower recovery.

⌚ 00:10:00 Air-drying

- 35 Remove the tube from the rack and resuspend beads in 50 µL DNase-free water by pipetting up and down. Place tube back on the rack and leave for **3 minutes**.

🧪 50 µL DNase-free water

⌚ 00:03:00 Magnetic rack

- 36 Pipette the 45 µL eluant from the tube while it is situated **on the magnetic tube rack**.

PCR 2

- 37 Complete thaw and vortex KAPA reagents (except for enzyme) before use.

38

	µl/tu be		[sto ck]	[fina l]	[mas term ix]
--	-------------	--	-------------	-------------	---------------------

	5.0	5x KAP A HiFi Fidel ity Buff er	5x	1x	
	1.0	dNT P Mix	10 mM	0.4 mM	
	1.0	Uni Ada pter (AD PT_P 1)	10 μM	0.4 μM	
	0.5	KAP A HiFi poly mer ase	1 U/ μl	0.5 U	
	1.0	Inde xed Ada pter	10 μM	0.4 μM	
	2.0	Purif ied tem plate DNA			
	14.5	Wat er			

Prepare mastermix in cold box and use repeater pipette to add to each tube in the cold box. Add Indexed Adapter to each tube.

🧪 5 μL 5x KAPA HiFi Fidelity Buffer

🧪 1 μL dNTP Mix

🧪 1 μL Uni Adapter (ADPT_P1)

🧪 0.5 μL KAPA HiFi polymerase

🧪 1 μL Indexed Adapter

🧪 2 μL Purified template DNA



🧪 14.5 μL Water

39 **Cycle** (PCR machine #5 SZ → ILM2):



	95°C	2 min
	98°C	20 s
	63°C	15 s
	72°C	30 s
	25 - 35 cycles	
	72°C	3 min
	4°C	On hold

Gel Purification and quantification

- 40 Before gel purification, run 2 µl products on 1% agarose gel to check the bands.
- 41 **Gel purification.** (Qiagen QIAquick gel extraction kit)
Run 2nd round PCR products on 1.2% agarose gel. E = 4 V/cm, T = 60 min.
- 42 Excise DNA fragment.
- 43 Weight the gel; add 3 volume of Buffer QG to 1 volume of gel.
- 44 Incubate at 50 °C for **10 minutes** to completing dissolve. Vortex every **2-3 minutes** to help dissolve.
 00:10:00 Incubation at 50 °C
 00:02:30 Vortex
- 45 Check the color of gel solution (should be yellow, otherwise add 10 µl 3M sodium acetate).
- 46 Place **MinElute** column, apply the sample to the column and centrifuge for **1 minute**.



00:01:00 Centrifugation

47 Add 500 µl buffer QG and centrifuge for **1 minute**.

500 µL Buffer QG

00:01:00 Centrifugation

48 Add 0.75 ml buffer PE, **incubate for 5 minutes** at room temperature, centrifuge for **1 minute**.

00:05:00 Incubation at room temperature

00:01:00 Centrifugation

49 Discard the fluid, centrifuge for additional **3 minutes**.

00:03:00 Centrifugation

50 Put the column in a new 1.7 ml tube, add **10 µl** buffer EB. Stand for **4 minutes**, centrifuge for **2 minutes**.

00:04:00 Stand

51 Quantification using Invitrogen Qubit dsDNA BR Assay kit. See Qubit dsDNA BR assay protocol. **Don't use Nanodrop to quantify!**

52 After quantification, pool libraries in equal amount. Use AMPure XP beads to purify pooled libraries (2:3, two washes), elute in 20 to 30 µL Elution Buffer.

53 Quantify the pooled library using Qubit dsDNA BR assay kit. Check the quality of the library by Agilent Bioanalyzer or Bio-rad Experion. Repeat purification if primer dimer is present (~200bp).