

Nov 20, 2020

🌐 Primers for bark beetles and associated fungi

DOI

<https://dx.doi.org/10.17504/protocols.io.bnvcme2w>



Andrew J. Johnson¹, Demian F Gomez¹, Jiri Hulcr¹

¹University of Florida

Protocols Bark Beetle M...



Bark Beetle Mycobiome Research Coordination Network Bark

Beetle

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bnvcme2w>

Document Citation: Andrew J. Johnson, Demian F Gomez, Jiri Hulcr 2020. Primers for bark beetles and associated fungi.
protocols.io <https://dx.doi.org/10.17504/protocols.io.bnvcme2w>

License: This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: October 23, 2020

Last Modified: November 20, 2020

Document Integer ID: 43652

Keywords: primers for bark beetle, bark beetle mycobiome, part of the bark beetle mycobiome, bark beetle, fungi this protocol, fungus symbiosis, associated fungi, beetle, widespread insect, fungi, fungus, symbiosis, more information on the bbm, primer, bark

Abstract

This protocol describes the primers used for bark beetles and associated fungi.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Seltmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. Symbiosis 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

Diluting primers

Reconstituting lyophilized primers to create 100X stock:

Resuspend to 100 uM/L concentration: for 28.5 nM of lyophilized primer, add 285 ul water (10 times more, in uL).

Dilute to make working solution :

To 10uL of the working stock above, add 90 uL of water, to make 10 uM/L (or 10 pM/uL) solution.

Working concentration: Use 1uL of the working solution above in 25uL reaction ($10/25=0.4\text{pM}$ per rxn).

Primers we use

	GEN E	PRI MER COD E	SEQ	REF ERE NCE
	16S	1492 R	GGT TAC CTT GTT ACG ACT T	
	16S	27F	AGA GTT TGA TCC TGG CTC AG	
	16S bact rev	338 R	TGC TGC CTC CCG TAG GAG T	Ama nn et al. 1995 . Micr obiol . Rev. 59, 143– 169.
	16S bact +arc haea for	515F	GTG CCA GC MG CCG CGG TAA	Berg man n et al 2011 (Fier er lab)
	16S bact +arc	806 R	GGA CTA CVS	Berg man n et



	haea rev		GGG TAT CTA AT	al 2011 (Fier er lab)
	28S	A40 53 (alt_ nam e: A1)	TCC KGT KTT CAA GAC GGG GTC	Whit ing et al. 1997
	28S	A42 21 (alt_ nam e: A160)	CGC CTC TTC TCG CAA TGA GA	Seq ueira et al, 200 0
	28S	A42 85 (alt_ nam e: A24 7)	CCT GAC TTC GTC CTG ACC AGG C	Seq ueira et al, 200 0
	28S	A43 94 (alt_ nam e: A33 5)	TCG GAR GGA ACC AGC TAC TA	Whit ing et al. 1997
	28S	D2F1	ACT GTT GGC GAC GAT GTT CT	Jord al etal. (200 8)
	28S	D3R 2	TCT TCG CCC CTA TAC CC	Jord al etal. (200 8)
	28S	f109 4	GGA TCG TCG CAA GAC GGA CAG AAG	Seq ueira et al, 200 0



	28S	f116 (alt_ name: e)	CTG GTT GAT CCT GCC ACG T	Ham by & Zim mer 1988
	28S	f140 3	CGG AAG GAT TGA CAG ATT GAG AG	Seq ueira et al, 200 0
	28S	f420	GGC GAC GCA TCT TTC AAA TGT CTG	Seq ueira et al, 200 0
	28S	r109 4 (alt_ name: P)	CGC TTT CGT AAA CGG TT	Ham by & Zim mer 1988
	28S	r113 8	CGC CTT CGA ACC TCT AAC	Seq ueira et al, 200 0
	28S	r162 6	GGC ATC ACA GAC CTG TTA TTG CTC AAT CTC	Seq ueira et al, 200 0
	28S	r185 6 (alt_ name: Q)	CAG CGA GGA TGG CTA ACT TA	Ham by & Zim mer 1988
	28S	r803	CCA CCG GCA	Seq ueira et al,



			GGA CGT CTC	200 0
	28S	S36 90F	GAG AGT TMA ASA GTA CGT GAA AC	Dow ton & Aust in 1998
	28S	S40 72 (alt_ nam e: S1)	GAC CCG TCT TGA AMC AMG GA	Whit ing et al. 1997
	Ap10 98re v2		ATA TTR TTK GGC ARY TGD CCK CCC	[not e: TmD 61.1 °C] [pri mer_ refer ence :Dan forth et al. 200 6]
	Ap7 87fo r2		TGC TTY GAR CCD AGY CTH GAT TAY TG	[not e: TmD 60.0 °C] [pri mer_ refer ence :Dan forth et al. 200 6]
	Ap8 35re v1		GCA THA CYT CHC CCA CRC TYT TC	[not e: TmD 59.8 °C] [pri mer_ refer ence :Dan forth et al.



				2006]
	ApC ADfor1		GG WTA TCC CGT DAT GGC BMG WGC	[not e: TmD 62.1 °C] [primer_ reference :Dan forth et al. 2006]
	ApC ADfor2mod		GAT GGG AYC TNR GNA ART TYC	[not e: TmD 53.1 °C] [primer_ reference :Dan forth et al. 2006]
	ApC ADfor3		CTC HGT KGA RTT YGA TTG GTG YGC	[not e: TmD 60.9 °C] [primer_ reference :Dan forth et al. 2006]
	ApC ADfor4		TGG AAR GAR GTB GAR TAC GAR GTG GTY CG	[not e: TmD 63.1 °C] [primer_ reference :Dan forth et al. 2006]



	ApC ADfor4		TGG AAR GAR GTB GAR TAC GAR GTG GYC G	This is not good, do not use.
	ApC ADrev1mod		GCC ATY RCY TCB CCY ACR CTY TTC AT	[note: TmD 62.2 °C] [primer_ reference :Dan forth et al. 200 6]
	ApC ADrev1mod		GCC ATY RCT CBC CTA CRC TYT TCA T	This is not good, do not use.
	ApC ADrev4a		GGC CAY TGN GCN GCC ACY GTG TCT ATY TGY TTN ACC	[note: TmD 68.6 °C] [primer_ reference :Dan forth et al. 200 6]
	ArgK	ArgK - forB 2	GAY TCC GG WAT YGG WAT CTA YGC TCC	Dole et al. 2010



	ArgK	ArgK - revB 1	TCN GTR AGR CCC ATW CGT CTC	Danf orth et al. 200 5
	ArgK	ArgK - revB 2	GTA TGY TCM CCR CGR GTA CCA CG	Dole et al. 2010
	cello bioh ydrol ase	fung CBH 1f		
	cello bioh ydrol ase	fung CBH 1r		
	CO1	1495 b	AAC AAA TCA TAA AGA TAT TGG RAC	Smit h and Cog nato, 2014
	CO1	A22 37R	CCG AAT GCT TCT TTT TTA CCT CTT TCT TG	
	CO1	A241 1R	GCT AAT CAT CTA AAA ACT TTA ATT CC WGT WG	
	CO1	HCO 2198	TAA ACT TCA	Heb ert et al.



			GGG TGA CCA AAA AAT CA	200 3
	CO1	LCO 1490	GGT CAA CAA ATC ATA AAG ATA TTG G	Heb ert et al. 200 3
	CO1	LepF 1	ATT CAA CCA ATC ATA AAG ATA TTG G	Heb ert et al. 200 4
	CO1	Lep R1	TAA ACT TCT GGA TGT CCA AAA AAT CA	Heb ert et al. 200 4
	CO1	r750	GAA ATT ATN CCA ATT CCT GG	Smit h and Cog nato, 2014
	CO1	S171 8F	GGA GGA TTT GGA AAT TGA TTA GTT CC	
	EF1a	A104 3	GTA TAT CCA TTG GAA ATT TGA CCN	Nor mar k et al. 1999



			GGR TGR TT	
	EF1a	S149	ATC GAG AAG TTC GAG AAG GAG GCY CAR GAA ATG GG	Nor mar k et al. 1999
	efTu	52F		Lud wig et al
	ITS (end of SSU)	ITS1	TCC GTA GGT GAA CCT GCG G	
	ITS (in 5.8S , for)	ITS3	GCA TCG ATG AAG AAC GCA GC	
	ITS (in 5.8S , rev)	ITS2	GCT GCG TTC TTC ATC GAT GC	
	ITS (rev, start of LSU)	ITS4	TCC TCC GCT TAT TGA TAT GC	
	ITS fung al (like ITS1, basi dio)	ITS1 -F	CTT GGT CAT TTA GAG GAA GTA A	



	LSU (28S)	JH- LSU -91	AGT AAC GGC GAG TGA AG	prob ably univ ersal for fung i, not for inse ct, posit ion 91 in S. cere visia e 28S
	LSU (28S)	LR0 R	ACC CGC TGA ACT TAA GC	ident ical in inse ct and fung i
	LSU (28S) rev	JH- LSU -377	TTT TCA AAG TGC TTT TC	almo st ident ical to LR2 but with one bp mor e spec ific to filam ento us asco s, posit ion 377 in S. cere visia e; does n't work very well
	LSU (28S)	JH- LSU	GGT CCG	reve rse



	) rev	-638	TGT TTC AAG AC	com plem ent to Vilga lis's LR3 R (whi ch is not the sam e as reve rse com plem ent of LR3)
	LSU (28S) rev	JH- LSU -643	CTC CTT GGT CCG TGT TTC	teste d with Raff aele a spp., inco nsist ent prod uct
	LSU (28S) rev	JH- LSU -830	GTT TCC TCT GGC TTC ACC C	teste d with Raff aele a spp., inco nsist ent prod uct
	LSU (28S) rev	LR2	TTT TCA AAG TTC TTT TC	Whit e et al., from yeas t, 1bp diffe renc e from filam ento us



				asco s
	LSU (28S) rev	LR3	CCG TGT TTC AAG ACG GG	ident ical for beetl es and fung i
	SSU (18S) gen eral	NS1	GTA GTC ATA TGC TTG TCT C	euka ryot es?
	SSU (18S) gen eral	NS4	CTT CCG TCA ATT CCT TTA AG	euka ryot es?
	SSU fung al	AU2	TTT CGA TGG TAG GAT AGD GG	Van denk oorn huys e et al. 200 2, Scie nce 295: 2051
	SSU fung al	AU4	RTC TCA CTA AGC CAT TC	Van denk oorn huys e et al. 200 2, Scie nce 295: 2051
	SSU fung al	JH- SSU -304	CCC TAT CAA CTT TCG ATG	upst rea m of AU2, 2bp diffe renc e from



				Curc ulion idae, can use with NS2 for 240 bp cons erve d prod uct or with JH- SSU -822 for ~50 0bp of varia ble prod uct
	SSU fung al (rev)	JH- SSU -164 6	CTA AGC CAT TCA ATC GGT A	6bp upst rea m of AU4, 2bp diffe renc e from Curc ulion idae
	SSU fung al (rev)	JH- SSU -822	CTA GAA ACC AAC AAA ATA G	2bp upst rea m of SR4, 4bp diffe renc e from Curc ulion idae
	SSU rev	NS2	GGC TGC TGG CAC CAG	550 in S. cere visia e; euka



			ACT TGC	ryot es?
--	--	--	------------	-------------

R= A G (purine) D= G A T

Buying primers

ITS1 example CType: DNA Oligos – modified & non-modified

Oligo Name: ...

Synthesis Scale: 0.025 μ mol

Purification: Desalt

Sequence (5' to 3'): ...

Format: Lab-ready

Number of Tubes per Oligo: 1

Concentration (μ M): NA