

Sep 09, 2020

## Priming and loading a MinION flowcell

DOI

[dx.doi.org/10.17504/protocols.io.7q5hmy6](https://dx.doi.org/10.17504/protocols.io.7q5hmy6)

Josh Quick<sup>1</sup>

<sup>1</sup>University of Birmingham

Diaz-Munoz Lab

Coronavirus Method De...

1 more workspace



Josh Quick

University of Birmingham

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7q5hmy6>

External link: <https://doi.org/10.1016/j.remle.2020.05.007>

**Protocol Citation:** Josh Quick 2020. Priming and loading a MinION flowcell. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.7q5hmy6>

**Manuscript citation:**

Fernández-Rodríguez A, Casas I, Culebras E, Morilla E, Cohen MC, Alberola J, COVID-19 and post-mortem microbiological studies. Spanish Journal of Legal Medicine doi: [10.1016/j.remle.2020.05.007](https://doi.org/10.1016/j.remle.2020.05.007)

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 27, 2019



**Last Modified:** September 09, 2020

**Protocol Integer ID:** 28157




**Keywords:** minion flowcell, minion

## Troubleshooting



- 1 Thaw the following reagents at room temperature before placing on ice:  
  
Sequencing buffer (SQB)  
Loading beads (LB)  
Flush buffer (FLB)  
Flush tether (FLT)
- 2 Add  30 µL FLT to the FLB tube and mix well by vortexing.
- 3 If required place a new MinION flowcell onto the MinION by flipping open the lip and pushing one end of the flowcell under the clip and pushing down gently.
- 4 Rotate the inlet port cover clockwise by 90° so that the priming port is visible.
- 5 Take a P1000 pipette and tip and set the volume to  800 µL . Place the tip in the inlet port and holding perpendicularly to the plane of the flowcell remove any air from the inlet port by turning the volume dial anti-clockwise.  

**Note**


Be careful not to remove so much volume that air is introduced onto the rectangular array via the outlet.
- 6 Load  800 µL of FLB (plus FLT) into the flow cell via the inlet port, dispense slowly and smoothly trying to avoid the introduction of any air bubbles.
- 7 Wait for  00:05:00 .
- 8 Gently lift the SpotON cover to open the SpotON port.
- 9 Load another  200 µL of FLB (plus FLT) into the flow cell via the inlet port, this will initiate a siphon at the SpotON port to allow you to load the library dilution.
- 10 In a new tube prepare the library dilution for sequencing:

**Component****Volume**

SQB

 37.5 µL

LB


 25.5 µL

Final library

 12 µL**Total** 75 µL**Note**

Mix LB immediately before use as they settle quickly.

Dilute library in EB if required.

- 11 Mix the prepared library gently by pipetting up and down just prior to loading.
- 12 Add the  75 µL library dilution to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop siphons into the port before adding the next.
- 13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the inlet port and close the MinION lid.