

Apr 07, 2025

Version 3

Processing frozen cells for population-scale SQK-LSK114 Oxford Nanopore long-read DNA sequencing SOP V.3

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr347bvmk/v3>



Jackson Mingle¹, Kimberly Paquette¹, Guillaume Cogan¹, Pilar Alvarez Jerez¹, Laksh Malik¹, Cornelis Blauwendraat¹, Kimberley J Billingsley¹, on behalf of the CARD Long-read Team¹

¹Center for Alzheimer's and Related Dementias, National Institute on Aging, Bethesda, Maryland, USA

NIH Center for Alzheime...



Guillaume Cogan

NIH

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr347bvmk/v3>

Protocol Citation: Jackson Mingle, Kimberly Paquette, Guillaume Cogan, Pilar Alvarez Jerez, Laksh Malik, Cornelis Blauwendraat, Kimberley J Billingsley, on behalf of the CARD Long-read Team 2025. Processing frozen cells for population-scale SQK-LSK114 Oxford Nanopore long-read DNA sequencing SOP . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr347bvmk/v3>Version created by **Jackson Mingle**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: February 28, 2025

Last Modified: April 07, 2025

Protocol Integer ID: 123600

Keywords: Long-read sequencing, Oxford Nanopore sequencing, High molecular weight DNA extraction, DNA extraction, Cultured cell extraction, DNA size selection, DNA shearing, read sequencing of cultured cell, read sequencing, Isk114 oxford nanopore, sequencing data, nanopore team, frozen cells for population, processing frozen cell, genome, nih center for alzheimer, sequencing, read dna, sequencing sop, patients with alzheimer, cell, alzheimer, lewy body dementia, dna, related dementia, dementia, cultured cell, genome sequencing

Disclaimer

Still in development.

Abstract

At the NIH Center for Alzheimer's and Related Dementias (CARD) (<https://card.nih.gov/research-programs/long-read-sequencing>), we will generate long-read sequencing data from thousands of patients with Alzheimer's disease, frontotemporal dementia, and Lewy body dementia, as well as healthy subjects. With this research, we will build a public resource consisting of long-read genome sequencing data from a large number of confirmed patients with Alzheimer's disease and related dementias and healthy individuals. To generate this large-scale nanopore sequencing data, we have developed a protocol for processing and long-read sequencing of cultured cells, targeting an N50 of ~25-30 kb and ~30X coverage.

Acknowledgements:

We would like to thank the Nanopore team (Jade Bartolo, Olivor Holman, Androo Markham, and Jessica Anderson), PacBio team (Jeffrey Burke, Michelle Kim, Duncan Kilburn, and Kelvin Liu), and the whole CARD long-read team. Workflow graphic by Paige Jarreau.

†Correspondence to: Kimberley Billingsley billingsleykj@nih.gov




Figure 1. Overview of the HMW DNA extraction and ONT sequencing protocol

Materials

Materials:

A	B
Material	Vendor (Part Number)
Nanobind HT CBB kit	PacBio (102-301-900)
KingFisher 96 deep-well plates, barcoded	Thermo Fisher Scientific (95040450B)
KingFisher 96 tip comb, barcoded	Thermo Fisher Scientific (97002534B)
1X PBS	Any major lab supplier (MLS)
Ethanol (96-100%)	Any MLS
Isopropanol (100%)	Any MLS
1.5 mL DNA LoBind tubes	Eppendorf (022431021)
Qubit 1X dsDNA BR assay kit	Thermo Fisher Scientific (Q33266)
Qubit Flex Assay Tube Strips	Thermo Fisher Scientific (Q33252)
Femto Pulse gDNA 165kb analysis kit	Agilent Technologies, Inc. (FP-1002-0275)
Megaruptor 3 DNAFluid+ kit	Diagenode (07020001)
TE Buffer	Thermo Fisher Scientific (12090)
2X Loading Solution	Sage Science (190192)
High Pass Plus Cassette	Sage Science (BPLUS10)
Ligation Sequencing Kit XL V14	Oxford Nanopore Technologies (SQK-LSK114-XL)
NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing	New England Biolabs (E7180L)
Nuclease-free water	Any MLS
0.2 mL thin-walled PCR tubes	Thermo Fisher Scientific (AB-0620B)
AMPure XP reagent	Beckman Coulter (A63882)
Qubit 1X dsDNA HS assay kit	Thermo Fisher Scientific (Q33231)
PromethION flow cell	Oxford Nanopore Technologies (FLO-PRO114M)

A	B
Flow Cell Wash Kit	Oxford Nanopore Technologies (EXP-WSH004-XL)

Table 1. Required materials

Equipment:

A	B
Equipment	Vendor (Part Number)
KingFisher Apex system	Thermo Fisher Scientific (5400940)
KingFisher Apex 96 DW head	Thermo Fisher Scientific (24079930)
KingFisher Apex 96 DW heating block	Thermo Fisher Scientific (24075930)
Qubit Flex fluorometer	Thermo Fisher Scientific (Q33327)
Vortex mixer	Any major lab supplier (MLS)
Minicentrifuge	Any MLS
Femto Pulse system	Agilent Technologies, Inc (M5330AA)
Megaruptor 3 system	Diagenode (B06010003)
BluePippin Size Selection System	Sage Science (BLU0001)
ThermoMixer	Eppendorf (5382000023)
Thermal cycler	Bio-Rad (1851197)
Platform rocker	Any MLS
DynaMag-2 Magnet	Thermo Fisher Scientific (12321D)
PromethION 24 or 48 sequencing unit	Oxford Nanopore Technologies (PRO-SEQ024 or PRO-SEQ-048)

Table 2. Required equipment

Materials for NGS STAR:

A	B
Material	Vendor (Part Number)
2 mL tubes	Sarstedt (72.693.005)

	A	B
	PCR ComfortLid	Hamilton (814300)
	MIDI plate	Thermo Fisher Scientific (AB-0859)
	HSP plate	Bio-Rad (HSP9601)
	50 µL tips	Hamilton (235979)
	300 µL tips	Hamilton (235903)
	1000 µL tips	Hamilton (235940)
	60 mL reservoir	Hamilton (56694-01)
	20 mL reservoir	Roche (3004058001)

Table 3. Additional materials required for automated library preparation on Hamilton NGS STAR

Equipment:

	A	B
	Equipment	Vendor (Part Number)
	NGS STAR	Hamilton (STAR)
	Magnetic Stand-96	Thermo Fisher Scientific (AM10027)

Table 4. Additional equipment required for automated library preparation on Hamilton NGS STAR

Part 1: Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on the KingFisher Apex system

- 1 Obtain frozen cell pellet from -80 °C freezer. Thaw on ice for  00:05:00 .

Note: Cell input amount should be 1×10^6 – 2.5×10^6 diploid human cells or equivalent. Cell counts should be accurately determined using a hemocytometer or cell counter. Cell pellets should be frozen dry with as much liquid removed as possible. No cryoprotectant is needed.

Note: For non-diploid or non-human cells, the cell input should be scaled appropriately to contain 5-25 µg of DNA. Warning: >25 µg DNA inputs can cause Nanobind disks to be "dropped" in the Lysis/Binding solution and/or cause well-to-well contamination.

Note: For fresh cell samples, harvest cells and centrifuge at 500 x *g* for 3-5 minutes at  4 °C to pellet cells. Remove the supernatant.

- 2 Prepare the KingFisher Apex 96 deep-well plates:
 - **Plate 1:** Lysis/Binding Plate: Sample and reagents from step 4
 - **Plate 2:** Nanobind Storage Plate: One 3 mm Nanobind disk per well
 - **Plate 3:** Wash Plate 1:  700 µL Buffer CW1 per well
 - **Plate 4:** Wash Plate 2:  700 µL Buffer CW2 per well
 - **Plate 5:** Wash Plate 3:  700 µL Buffer CW2 per well
 - **Plate 6:** Elution Plate :  100 µL Buffer EB per well
 - **Plate 7:** Tip Plate: KingFisher Apex 96 deep-well tip comb

Note: Nanobind disks do not need to be perfectly centered in the wells, but ensure they are at the bottom of the wells and not stuck to the side walls.

Note: Buffer CW1 and CW2 are supplied as concentrates. This kit uses CW1 with a 60% final ethanol concentration and CW2 with a 60% final ethanol concentration. Before using, add the appropriate amount of ethanol (96–100%) to Buffer CW1 and Buffer CW2 as indicated on the bottles.

- 3 Add  50 µL of 1X PBS to each cell pellet. Pipette-mix 10x with a standard P200 pipette to resuspend cells.

Note: Mix until cell pellet is fully resuspended without visible lumps. Sticky cell types may require additional pipette-mixing or vortexing.

Note: Aggressive mixing at this step will not affect DNA size. However, incomplete resuspension will result in inefficient lysis and digestion which will lead to low yield, low purity, and high heterogeneity.

4 Prepare the Lysis/Binding Plate:

- Add  50 μL resuspended cells per well
- Add  20 μL Proteinase K per well
- Add  5 μL Buffer CLE3 per well
- Add  20 μL RNase A per well

Note: The sample and reagents MUST be added to the wells in the order described above.

5 Let rest for 00:05:00 at room temperature.

6 Add 150 μL Buffer BL3 per well.

Note: Add Buffer BL3 gently against the side of the well into the Lysis/Binding solution. Adding Buffer BL3 directly to the Lysis/Binding solution may affect extraction performance.

7 Ensure the KingFisher Apex instrument is set up with the 96 deep-well magnetic head and the 96 deep-well heating block.

8 Select the Cell_Nanobind_HT_APEX script on the KingFisher Apex instrument. Insert plates into the KingFisher Apex instrument as indicated on the display and press 'Next' after every plate to confirm position. The protocol will start when the final plate is loaded and the 'Next' button is pressed.

9 When prompted by the instrument (~21 minutes after start), remove the Lysis/Binding Plate from the instrument and add 250 μL of isopropanol to each well. Re-insert the plate and press 'Next' to resume the protocol.

Note: Add isopropanol gently against the side of the well into the Lysis/Binding solution. Adding isopropanol directly to the Lysis/Binding solution may affect extraction purity.

10 When prompted by the instrument, remove the Elution Plate from the instrument. Transfer eluate from each well to a new 1.5 mL Eppendorf DNA LoBind tube.

Note: Residual sample volume may be present in the tip comb plate. Transfer residual eluate to the 1.5 mL Eppendorf DNA LoBind tube.

- 11 Pipette-mix the sample 10x with a standard P200 pipette to homogenize and disrupt any unsolubilized "jellies" that may be present.
- 12 Let eluate rest overnight at room temperature to allow DNA to solubilize.
- 13 Follow overnight rest, pipette-mix 10x with a standard P200 pipette.
- 14 Quantify by taking a single measurement on the Qubit Flex Fluorometer with the Qubit 1x dsDNA BR Assay. Use  2 μ L of DNA per measurement.

Note: Triplicate measurements (top, middle, and bottom) are advised for new cohorts or cell types. Additional hand-shearing with a syringe may be necessary for samples with exceptionally variable measurements.

- 15 **Optional:** Quantify by taking a single measurement on the NanoDrop 8000 spectrophotometer. Use  2 μ L of DNA per measurement.

Note: Measurement on NanoDrop 8000 is advised for new cohorts or brain regions. Additional hand-shearing with a syringe may be necessary for samples with exceptionally variable measurements.

- 16 **Optional:** Size on the Agilent Femto Pulse System with the Genomic DNA 165 kb kit. The expected size range for samples post-extraction is >100kb.




Figure 2. Femto Pulse trace post-extraction



Note: Sizing on Agilent Femto Pulse System is advised for new cohorts or cell types. Additional hand-shearing with a syringe may be necessary for samples with exceptionally variable measurements.

- 17 Store at  4 °C until ready for shearing.

Part 2: Shearing with Diagenode Megaruptor 3 DNAFluid+ Kit

- 18 Prepare DNA in  100 µL Buffer EB, at a concentration of up to 150 ng/µL in a Megaruptor 3 shearing tube.

Note: Shearing conditions are optimized for sequencing N50s ~25-35 kb. Shearing with Megaruptor 3 is concentration- and volume-dependent. Samples with a concentration > 150 ng/µL may be susceptible to under-shearing. Do not shear at a concentration > 150 ng/µL.

Note: For samples with concentrations > 150 ng/µL, prepare DNA in > 100 µL, and ensure sample concentration is under 150 ng/µL. Ensure the volume is input on MegaRuptor 3 accurately in step 21.

Note: For samples with low yields, shear the entire sample in 100 µL, then pool with additional samples after shearing to reach the required input for size selection, if desired.

- 19 Remove the Megaruptor 3 DNA Fluid+ shearing syringe from the package. Tighten the assembly and ensure the plunger is completely compressed. Attach the syringe to the shearing tube and ensure the syringe is snug on the cap of the tube.
- 20 Load samples and syringes onto the Megaruptor 3. If running an odd number of samples, samples can be balanced with an empty corresponding tube.
- 21 Run samples on the Megaruptor 3 at speed 45. Enter the appropriate sample volume (100 µL) and sample concentration.

Note: Sample volume entered into MegaRuptor 3 protocol must accurately reflect the sample volume. If running samples with different volumes (due to high concentrations), samples must be run on separate MegaRuptor 3 protocols.

Note: Ensure samples of different volumes and concentrations are run on separate MegaRuptor 3 protocols. For different volumes, keep samples within +/- 10 µL of each other. For different concentrations, keep samples within +/- 15 ng/µL of each other.

Note: Shearing speed may need to be optimized for new cohorts or cell types to achieve targeted size.

- 22 Repeat the previous step for a total of 2 runs on the Megaruptor 3 at speed 45.
- 23 Remove the samples from the Megaruptor 3. Remove the syringe from the tube. Make sure the plunger is fully depressed to avoid losing sample volume. Disassemble the needle and remove residual sample from the needle and syringe to avoid losing sample volume.
- 24 Quantify by taking a single measurement on the Qubit Flex Fluorometer with the Qubit 1x dsDNA BR Assay. Use 1 μL of DNA per measurement.
- 25 Size on the Agilent Femto Pulse System with the Genomic DNA 165 kb kit. The expected size range for samples post-shearing is ~25-40 kb.




Figure 3. Femto Pulse trace post-shearing

- 26 Store at 4 °C until ready for size selection. For long-term storage, store at -80 °C .

Part 3: Size Selection with Sage Science BluePippin

- 27 Prepare up to 10.5 μg sheared DNA in 77 μL TE buffer or Buffer EB.

Note: If samples are low volume, low concentration, or limited due to lack of starting material, samples can be prepared in 32 μL or 62 μL TE buffer.

28 Remove 2X Loading Solution from 4 °C fridge and allow it to reach room temperature.

Note: If preparing samples at a lower volume of  32 µL or  62 µL , standard (1X) Loading Solution should be used and allowed to come to room temperature.

29 Create a new protocol using the BluePippin software and select the cassette type:

- Click the "Protocol Editor" tab
- Click the "Cassette Type" folder icon
- Click "+" icon next to "0.75% Agarose Dye-Free"
- Select "10kb High Pass Plus Marker U1"
- Click "SELECT"


Note: The standard protocol used is "10kb High Pass Plus Marker U1". If samples require more aggressive size selection following analysis of post-shearing Femto Pulse traces, "15kb High Pass Plus Marker U1" or "20kb High Pass Plus Marker U1" protocols can be run instead.

30 Enter the sample IDs into the Sample ID fields. Enter Marker for the Marker well in Lane 4.

Note: Lane 4 is the preferred lane for the Marker, but other lanes can be used.

31 Make sure the check box for "End Run when Elution is Completed" is selected.

32 In the "Range" column, ensure the rectangle boxes next to sample lanes are light gray, or selected. Ensure the rectangle box in the Marker lane is dark gray, or not selected.



Select Range Mode and de-select Lane 4 (for loading Marker U1)

Figure 4. Range mode selection

33 Enter the lane number to which the DNA marker will be loaded into the "Reference Lane" field, and click the "APPLY REFERENCE TO ALL LANES" button.

34 Enter the "BP Start" value of 10000 for all sample wells.


Note: The "BP End" and "BP Target" values are automatically populated.

Note: If using the "15kb High Pass Plus Marker U1" or "20kb High Pass Plus Marker U1" protocols, "BP Start" values should be entered as 15000 and 20000, respectively.

35 Click "Save As" to name and save the run.

36 Remove the gel cassette from the foil packaging and inspect the levels of buffer in all buffer reservoirs.

Note: Reservoirs should be nearly full. If the buffer level in any reservoir appears lower, or less than 50% full, fill with spare electrophoresis buffer.



If a buffer reservoir is low, fill with spare buffer.

Figure 5. Buffer reservoir volume inspection

37 Inspect for bubbles due to delamination of agarose from the bottom of the cassette in the region used for optical detection of DNA. To inspect, turn the cassette upside down and gently tilt the cassette under a light source. If a bubble is detected, do not use the affected lane to run the DNA reference marker. Affected lanes can be used for size selection of sample DNA.

Note: Sometimes, flat bubbles can occur between the top surface of the gel column and the plastic top of the channel. This will not affect run quality or optical detection and can be used for markers or samples without any adverse effects. If it is difficult to tell if a bubble will be problematic, it may be better to not use the affected lane for the marker. If

a gel piece is detached, do not use that particular lane at all. Sage Science may be able to provide a replacement cassette.



If the agarose has delaminated in the optical region, do not use that lane to run the reference marker .

Figure 6. High Pass Plus cassette inspection

- 38 Dislodge air bubbles from behind the elution wells. Tilt the cassette sample well side down, with the lower buffer chambers facing up, and tap gently to release the bubbles.




Figure 7. High Pass Plus cassette diagram

- 39 Place the cassette into the optical nest. Keep the cassette slightly tilted down so that the air bubbles do not return to the area behind the elution modules. Make sure the cassette is fully seated into the bottom of the nest to ensure proper optical alignment.
- 40 Remove the white tabbed adhesive strips from the cassette. Place one hand on the cassette, and hold it firmly in the nest. Grab the white tabs and pull the strips firmly and slowly toward the front of the instrument until they are removed.
- 41 If necessary, fill buffer reservoirs with spare electrophoresis buffer if any reservoir appears low or less than 50% full. For any buffer reservoirs that appear full, remove  750 μL buffer to reduce leakage at the end of the run.
- 42 Completely remove all buffer from each elution module using a P100 pipette tip. Move the pipette tip slowly to avoid piercing the gel at the corner of the elution modules.
- 43 Add  80 μL of fresh electrophoresis buffer to each elution module. Change pipette tips between elution modules. When refilling, dispense the fresh buffer slowly, raising the pipette tip slowly from the bottom of the elution module to avoid trapping air bubbles in the module.
- 44 Seal the elution modules with an adhesive tape strip. Place the tape over the elution wells and rub firmly to seal the elution ports using a straight edge. The port should be tightly sealed without any wrinkles around the edges of the port. If necessary, use two overlapping tape strips.
- 45 Add  11 μL 2X Loading Solution to the prepared DNA samples. Pipette-mix or flick mix, and spin down.

Note: The 2X Loading Solution is very viscous, so care should be taken to aspirate and dispense fully.


Note: If preparing samples at  32 μL or  62 μL DNA,  10 μL or  20 μL Standard (1X) Loading Solution should be added, respectively.

46 Gently remove the cassette from the optical nest.

47 Calibrate the optical sensors:

- Click the "Main" tab
- Place the calibration fixture onto the optical nest — the dark side of the fixture must be down, and completely cover all LED detectors
- Close the lid of the BluePippin instrument
- Press the "CALIBRATE" button; a calibration sub-window will launch
- Check that the calibration target box ("Target I pH, mA") at the top of the window reads "0.60"; if not, manually change the target value to "0.60"
- Press the "CALIBRATE" button within the calibration window
- If the status window reads "Calibration OK", press "EXIT" to return to the main screen

Note: Calibration should be performed no more than 30 minutes before the start of a run. If more than 30 minutes elapses between calibration and the start of the run, a new calibration should be performed.



Correct placement of the calibration fixture.

Figure 8. Calibration fixture placement

- 48 Place the cassette into the optical nest. Keep the cassette slightly tilted down so that the air bubbles do not return to the area behind the elution modules. Make sure the cassette is fully seated into the bottom of the nest to ensure proper optical alignment.
- 49 Close the lid and perform the Continuity Test:
- Press the "Test" button — the test sub-window will open, and the test routine will automatically measure the current in each separation and elution channel
 - A successful Continuity Test will return a "PASS" indication
 - The cassette temperature must be above 17 °C (62 °F); if a lane fails due to low elution current (the affected lane will be highlighted with an orange color), refilling the elution module will usually resolve the problem
 - Press "Return" to close the Continuity Test sub-window
- 50 Remove  85 µL of buffer from sample wells and load  85 µL of sample into the sample well. Take care not to pierce the gel with the pipette tip.

Note: When removing buffer, it is useful to immerse the pipette tip just below the surface of the buffer and follow the liquid level down with the tip as the buffer is removed. When buffer removal is completed, there will be ~70 µL of buffer left in the well. When adding sample, place the tip of the pipette just below the surface of the buffer, and slowly eject the sample. Don't be concerned if the sample well slightly overfills. The density of the sample will allow it to sink before it can flow out of the well.

Note: Cassettes may be reused if 2-3 lanes are used (1-2 samples). Reseal the cassette with the adhesive tape and store at room temperature. However, one remaining lane must again be loaded with the DNA marker.

Note: If preparing samples at  32 µL or  62 µL DNA, remove  40 µL or  80 µL of buffer from sample wells and load  40 µL or  80 µL of sample into the sample well, respectively.

- 51 Remove DNA marker from 4 °C fridge. Do not keep out any longer than needed. Return to 4 °C fridge promptly after use.
- 52 Remove  40 µL of buffer from the marker well and load  40 µL of marker into the marker well. Take care not to pierce the gel with the pipette tip.

Note: When removing buffer, it is useful to immerse the pipette tip just below the surface of the buffer and follow the liquid level down with the tip as the buffer is removed. When buffer removal is completed, there will be ~70 µL of buffer left in the well. When adding sample, place the tip of the pipette just below the surface of the buffer, and slowly eject

the sample. Don't be concerned if the sample well slightly overfills. The density of the sample will allow it to sink before it can flow out of the well.

53 Close the lid. Ensure the proper protocol is loaded in the "Protocol Name" field in the "Main" tab.

54 Press "START". The run will automatically stop when every collection is complete.

Note: Samples can be left in the BluePippin following elution overnight. Remove samples from the elution modules first thing the following morning.

55 For best recovery, wait at least 30-45 minutes after the end of the run before removing samples from the elution modules.

56 Remove samples from the elution modules using a standard P200 pipette. Aspirate slowly to minimize breakage of large DNA fragments. Samples will be ~80 µL.

Note: Additional 10-30% yield can often be achieved by rinsing the elution well with TE+0.1% Tween 20 (provided by Sage Sciences). Add 80 µL Tween solution to the elution well, wait 1 minute, then remove solution and pool with original extracted sample, or process separately.

57 Remove the cassette and dispose of properly. Do not keep used cassettes in the BluePippin with the cover closed. Humidity from the cassette may cause damage to the electrodes.

58 Quantify by taking a single measurement on the Qubit Flex Fluorometer with the Qubit 1x dsDNA BR Assay. Use  1 µL of DNA per measurement.

59 Size on the Agilent Femto Pulse System with the Genomic DNA 165 kb kit. The expected size range for samples post-size selection is ~24-40 kb, with a DQN > 9.8.




Figure 9. Femto Pulse trace post-shearing and size selection

- 60 Store at 4 °C until ready for library preparation. At least 2.5 µg DNA are required for library preparation. If necessary, repeat shearing and/or size selection to achieve required mass for library preparation.

Part 4a: Ligation Sequencing DNA V14 (SQK-LSK114) library preparation

- 61 **Note:** For automated library preparation, see Part 4b: Ligation sequencing DNA V14 (SQK-LSK114) library preparation on Hamilton NGS STAR.

Thaw the following reagents at room temperature, then store on ice:

- NEBNext FFPE DNA Repair Buffer
- Ultra II End-Prep Reaction Buffer
- Ligation Buffer (LNB)
- Elution Buffer (EB)

Place the following reagents at room temperature:

- AMPure XP beads
- Short Fragment Buffer (SFB)
- Qubit 1X HS dsDNA working solution and standards **NOTE:** Qubit 1X HS dsDNA working solution must be stored away from light

- 62 Prepare 2.5 µg DNA in 48 µL nuclease-free water in a 0.2 mL thin-walled PCR tube.

Note: If necessary, input volume may be greater than 48 μL to accommodate 2.5 μg DNA input mass. Input volume should not exceed 168 μL .

63 Prepare the following DNA repair enzyme master mix, adjusting the reagent volumes for the number of samples being processed:

-  3.5 μL NEBNext FFPE DNA Repair Buffer (vortex and spin down)
-  3.5 μL Ultra II End-Prep Reaction Buffer (vortex and spin down)
-  2 μL NEBNext FFPE DNA Repair Mix (do not vortex, spin down)
-  3 μL Ultra II End-Prep Enzyme Mix (do not vortex, spin down)

Pipette gently 10x to mix. Avoid pipetting bubbles.

Note: Keep NEB enzyme mixes in the freezer until use and return promptly.

Note: Prepare the DNA repair enzyme master mix with a 10% overage.

64 Add  12 μL DNA repair enzyme master mix. Pipette-mix or flick mix 10x, and spin down.

65 Incubate samples at  20 $^{\circ}\text{C}$ for  00:30:00, followed by  65 $^{\circ}\text{C}$ for  00:05:00 in a thermocycler.

Note: Start and pause thermocycler to allow lid to come to 85 $^{\circ}\text{C}$ before putting samples in.

66 Transfer samples to new 1.5mL Eppendorf DNA LoBind tubes.

67 Resuspend AMPure XP beads by vortexing.

68 Add  60 μL or equal sample volume of resuspended AMPure XP beads. Flick mix 10x and spin down.

Note: AMPure XP bead volume should be adjusted to match the input sample volume and enzyme mix post-incubation.

69 Incubate on platform rocker for  00:05:00 at room temperature.

70 Prepare  500 μL of fresh 80% ethanol per sample in nuclease-free water.



- 71 Spin down samples and place on magnetic rack. Wait until supernatant is clear and colorless, about  00:02:00 .
- 72 Keep the tube on the magnetic rack and pipette off the supernatant.
- 73 With the samples remaining on the magnetic rack, add  200 μL of 80% ethanol. Pipette on the opposite wall to avoid disturbing the pellet. After  00:00:05 , remove the ethanol. Do not resuspend the beads in ethanol.
- Note:** If initial volume of beads was significantly higher than 60 μL , more ethanol may be required to keep the beads fully covered.
- 74 Repeat the previous step.
- 75 Spin down and place tubes on the magnetic rack. Pipette off any residual ethanol.
- 76 Allow to dry for ~  00:00:30 , but do not overdry to the point of cracking.
- 77 Remove tubes from the magnetic rack. Add  60 μL nuclease-free water. Flick mix 10x and spin down.
- 78 Incubate for  00:03:00 at  37 $^{\circ}\text{C}$ and 450 rpm in Thermomixer.
- 79 Spin down and place the samples on the magnetic rack until eluate is clear and colorless, about  00:02:00 .
- 80 Transfer  60 μL of eluate into a new 1.5 mL Eppendorf DNA LoBind tube.
- Note:** It is possible to store samples at  4 $^{\circ}\text{C}$ overnight at this step.
- 81 Prepare the following ligation enzyme master mix, adjusting the reagent volumes for the number of samples being processed:
-  25 μL Ligation Buffer (pipette-mix, spin down)
 -  10 μL Quick T4 DNA Ligase (do not vortex, spin down)

-  5 μL Ligation Adapter (do not vortex, spin down)

Pipette gently 10x to mix. Avoid pipetting bubbles.

Note: Keep Quick T4 DNA Ligase and Ligation Adapter in the freezer until use and return promptly.

Note: Prepare the ligation enzyme master mix with a 10% overage.

82 Add  40 μL ligation enzyme master mix. Pipette-mix or flick mix 10x, and spin down.

83 Incubate the samples for  00:30:00 at room temperature.

84 Resuspend AMPure XP beads by vortexing.

85 Add  45 μL of resuspended AMPure XP beads. Flick mix 10x and spin down.

86 Incubate on platform rocker for  00:05:00 at room temperature.

87 Spin down samples and place on magnetic rack. Wait until supernatant is clear and colorless, about  00:02:00 .

88 Keep the tube on the magnetic rack and pipette off the supernatant.

89 With the samples remaining on the magnetic rack, add  250 μL of Short Fragment Buffer. Remove samples from the magnetic rack, flick mix until beads are fully resuspended, and spin down.

Wait until supernatant is clear and colorless, about  00:02:00 , then remove the supernatant.

90 Repeat the previous step.



- 91 Spin down and place tubes on the magnetic rack. Pipette off any residual Short Fragment Buffer.
- 92 Allow to dry for ~ 00:00:30 , but do not overdry to the point of cracking.
- 93 Remove tubes from the magnetic rack. Add 26 μL Elution Buffer. Flick mix 10x and spin down.
- 94 Incubate for 00:20:00 at 37 °C and 450 rpm in Thermomixer.
- 95 Spin down and place the samples on the magnetic rack until eluate is clear and colorless, about 00:02:00 .
- 96 Transfer 26 μL of eluate into a new 1.5 mL Eppendorf DNA LoBind tube.
- 97 Quantify by taking a single measurement on the Qubit Flex Fluorometer with the Qubit 1X dsDNA HS Assay. Use 1 μL of DNA per measurement.
- 98 Store at 4 °C until ready for loading. At least 750 ng DNA is required for loading for a 72 hour run. If necessary, repeat library preparation with additional sheared and size selected DNA to achieve required mass for loading. For long-term storage, store libraries at -80 °C .

Note: Recommended DNA input per load is 20 fmol, calculated to be 300 ng of 24kb dsDNA using <https://nebiocalculator.neb.com/#!/dsdnaamt>. For 3 full loads, 900 ng DNA is required. At a minimum, 750 ng is required for 2 full loads of 300 ng each, followed by an incomplete load comprised of the recovered DNA library from the first load and the remaining DNA library under 300 ng.

Note: If greater than 900 ng DNA library is available, PromethION flow cell may be loaded 4 times in a 96 hour run for improved data output.

Part 4b: Ligation Sequencing DNA V14 (SQK-LSK114) library preparation on Hamilton NGS STAR

- 99 **Note:** For manual library preparation, see Part 4a: Ligation sequencing DNA V14 (SQK-LSK114) library preparation.

Thaw the following reagents at room temperature, then store on ice:

- NEBNext FFPE DNA Repair Buffer
- Ultra II End-Prep Reaction Buffer
- Ligation Buffer (LNB)
- Elution Buffer (EB)

Place the following reagents at room temperature:

- AMPure XP beads
- Short Fragment Buffer (SFB)
- Qubit 1X HS dsDNA working solution and standards **NOTE:** Qubit 1X HS dsDNA working solution must be stored away from light

- 100 Prepare  2.5 µg DNA in  48 µL nuclease-free water in a BioRad Hard-Shell PCR Plate.

Note: DNA input mass should be  2 µg at a minimum. Input volume may not be greater than  48 µL .

Note: Samples must be processed in multiples of 8. If preparing less samples than a multiple of 8, fill remaining wells with  48 µL nuclease-free water.

- 101 Launch the LSK-114 protocol in the Hamilton Method Manager.

Note: Ensure the deck layout of the NGS STAR matches the protocol being run.

Note: Daily and weekly maintenance must be performed prior to beginning the protocol. Follow manufacturer guidelines for preventive maintenance and cleaning procedures.

- 102 Select "Process01: DNA repair and end-prep" as the start process.



Figure 10. Select start process

- 103 Select "Process02: DNA repair and end-prep clean-up" as the stop process.



Figure 11. Select stop process

104 Select the input file worklist.




Figure 12. Select input file worklist

Note: Samples must be processed in multiples of 8. If preparing less samples than a multiple of 8, fill remaining wells with  48 μ L nuclease-free water.

- 105 Prepare the DNA repair enzyme mix in 2 mL Sarstedt tubes according to the instructions in the method. Avoid pipetting bubbles.




Figure 13. Prepare DNA repair enzyme mix

- 106 Load PCR ComfortLids onto the NGS STAR deck according to the instructions in the method.




Figure 14. Load PCR ComfortLids

- 107 Load MIDI plates and HSP plates onto the NGS STAR deck according to the instructions in the method.




Figure 15. Load MIDI and HSP plates

- 108 Load 50 μ L tips onto the NGS STAR tip carriers according to the instructions in the method. Use the tip count screen to input loaded 50 μ L tips.



Figure 16. Load 50 μ L tips




Figure 17. Select loaded 50 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

- 109 Load 300 μ L tips onto the NGS STAR tip carriers according to the instructions in the method. Use the tip count screen to input loaded 300 μ L tips.




Figure 18. Load 300 μ L tips




Figure 19. Select loaded 300 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

- 110 Load 80% ethanol in a 60 mL reagent trough in the reagent trough carrier according to the instructions in the method.




Figure 20. Load 80% ethanol

- 111 Load nuclease-free water in a 60 mL reagent trough and AMPure XP beads in a 20 mL reagent trough in the reagent trough carrier according to the instructions in the method.




Figure 21. Load nuclease-free water and AMPure XP beads

Note: AMPure XP beads must be resuspended by vortexing prior to use.

- 112 Load 1000 μ L tips and the prepared sample plate with  48 μ L DNA per well onto the NGS STAR deck according to the instructions in the method. Use the tip count screen to

input loaded 1000 μ L tips.




Figure 22. Load 1000 μ L tips and DNA sample plate




Figure 23. Select loaded 1000 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

Note: Ensure the magnetic stand is placed on the carrier in the 5th position.

- 113 Load the DNA repair enzyme mix tubes onto the CPAC cooler according to the instructions in the method.




Figure 24. Load DNA repair enzyme mix

Note: Enzyme mixes should be spun down and placed on the NGS STAR deck with the caps removed. Take care to avoid pipetting bubbles in the tubes to avoid inaccurate dispensing.

- 114 Close the NGS STAR door and press 'Ok.' The protocol will now begin.
- 115 When prompted at the end of the run, remove the sample plate from the NGS STAR. Discard the reagent troughs and enzyme mix tubes.




Figure 25. Run completion dialog

Note: It is possible to store samples at  4 °C overnight at this step.

116 Launch the LSK-114 protocol in the Hamilton Method Manager.

Note: Ensure the deck layout of the NGS STAR matches the protocol being run.

Note: Daily and weekly maintenance must be performed prior to beginning the protocol. Follow manufacturer guidelines for preventive maintenance and cleaning procedures.

117 Select "Process03: Adapter ligation" as the start process.



Figure 26. Select start process

- 118 Select "Process04: Adapter ligation clean-up" as the stop process.



Figure 27. Select stop process

- 119 Select the input file worklist.




Figure 28. Select input file worklist

Note: Samples must be processed in multiples of 8. If preparing less samples than a multiple of 8, fill remaining wells with  60 μL nuclease-free water.

- 120 Prepare the ligation enzyme mix in 2 mL Sarstedt tubes according to the instructions in the method. Avoid pipetting bubbles.




Figure 29. Prepare ligation enzyme mix

- 121 Load MIDI plates and HSP plates onto the NGS STAR deck according to the instructions in the method.




Figure 30. Load MIDI and HSP plates

- 122 Load 50 μ L tips onto the NGS STAR tip carriers according to the instructions in the method. Use the tip count screen to input loaded 50 μ L tips.




Figure 31. Load 50 μ L tips




Figure 32. Select loaded 50 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

- 123 Load 300 μ L tips onto the NGS STAR tip carriers according to the instructions in the method. Use the tip count screen to input loaded 300 μ L tips.




Figure 33. Load 300 μ L tips




Figure 34. Select loaded 300 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

- 124 Load AMPure XP beads in a 20 mL reagent trough, Short Fragment Buffer in a 60 mL reagent trough, and Elution Buffer in a 60 mL reagent trough in the reagent trough carrier according to the instructions in the method.




Figure 35. Load AMPure XP beads, Short Fragment Buffer, and Elution Buffer

Note: AMPure XP beads must be resuspended by vortexing prior to use.

- 125 Load 1000 μ L tips and the prepared sample plate with  60 μ L DNA per well onto the NGS STAR deck according to the instructions in the method. Use the tip count screen to input loaded 1000 μ L tips.




Figure 36. Load 1000 μ L tips and DNA sample plate




Figure 37. Select loaded 1000 μ L tips

Note: Method should be started with full tip racks to avoid unnecessary errors.

Note: Ensure the magnetic stand is placed on the carrier in the 5th position.

- 126 Load the ligation enzyme mix tubes onto the CPAC cooler according to the instructions in the method.




Figure 38. Load ligation enzyme mix

Note: Enzyme mixes should be spun down and placed on the NGS STAR deck with the caps removed. Take care to avoid pipetting bubbles in the tubes to avoid inaccurate dispensing.

- 127 Close the NGS STAR door and press 'Ok.' The protocol will now begin.
- 128 When prompted at the end of the run, remove the sample plate from the NGS STAR. Discard the reagent troughs and enzyme mix tubes.




Figure 39. Run completion dialog

- 129 Transfer  26 μL of eluate into a new 1.5 mL Eppendorf DNA LoBind tube.
- 130 Quantify by taking a single measurement on the Qubit Flex Fluorometer with the Qubit 1X dsDNA HS Assay. Use  1 μL of DNA per measurement.
- 131 Store at  4 °C until ready for loading. At least  750 ng DNA is required for loading. If necessary, repeat library preparation with additional sheared and size selected DNA to achieve required mass for loading. For long-term storage, store libraries at  -80 °C .

Note: Recommended DNA input per load is 20 fmol, calculated to be  300 ng of 24kb dsDNA using <https://nebiocalculator.neb.com/#!/dsdnaamt>. For 3 full loads,  900 ng DNA is required. At a minimum,  750 ng DNA is required for 2 full loads of  300 ng each, followed by an incomplete load comprised of the recovered DNA library from the first load and the remaining DNA library under 300 ng.

Note: If greater than  900 ng DNA library is available, PromethION flow cell may be loaded 4 times in a 96 hour run for improved data output.

Part 5: Loading the PromethION Flow Cell

- 132 Thaw the following reagents at room temperature:
- Sequencing Buffer (SB)
 - Library Solution (LIS)
 - Flow Cell Tether (FCT)
 - Flow Cell Flush (FCF)
 - Elution Buffer (EB)

Mix by vortexing, spin down, and place on ice.

- 133 Take flow cells out of the fridge. Allow to sit at room temperature for  00:20:00 .

Note: Condensation can form on the flow cell in humid environments. Inspect the gold connector pins on the top and underside of the flow cell for condensation and wipe off with a lint-free wipe if any is observed. Ensure the heat pad (black pad) is present on the underside of the flow cell.

- 134 Load flow cells into the PromethION docking ports. Perform flow cell check prior to flow cell priming.

Note: Only flow cells with a starting pore count greater than 6500 pores should be used for sequencing runs. Flow cells with starting pore counts less than 6500 pores can be used to generate additional sequencing data after an initial run.

135 Prepare the flow cell priming mix by combining  1170 μL Flow Cell Flush and  30 μL Flow Cell Tether per sample. Mix by vortexing.

136 Prepare  300 ng DNA library in  32 μL Elution Buffer in a 1.5 mL Eppendorf DNA LoBind tube.

Note: Recommended DNA input is 20 fmol, calculated to be 300 ng of 24kb dsDNA using <https://nebiocalculator.neb.com/#!/dsdnaamt>.

137 Add  100 μL Sequencing Buffer (SB).

138 Add  68 μL Library Solution (LIS).

139 Slide the flow cell inlet port cover clockwise to open. Draw back a small volume to remove any air bubbles:

- Set a P1000 pipette tip to 200 μL
- Insert the tip into the inlet port
- Turn the wheel until dial shows 220-230 μL , until a small volume of buffer enters the pipette tip

Note: Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 μL , and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.

140 Load  500 μL of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles.

141 Wait  00:05:00 .

142 Load  500 μL of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles.

143 Mix the prepared library gently by pipetting up and down 5x just prior to loading.

- 144 Load  200 µL of library into the inlet port using a P1000 pipette.
- 145 Close the valve to seal the inlet port. Install the light shield. Close the PromethION door.
- 146 Wait  00:10:00 before initiating sequencing run in MinKNOW:
- Navigate to the start page and click "Start sequencing"
 - Fill in the experiment name and sample ID, select the flow cell position, and load run configuration preset (if applicable)
 - Select the sequencing kit used in the library preparation (SQK-LSK114-XL) on the Kit page
 - Configure sequencing and output parameters - Fast basecalling model, 72 or 96 hour run time, basecalled output off, raw reads output .POD5, and minimum Q score of 8
 - Click "Start" to initiate the sequencing run
- 147 Within the first hour of sequencing, pay attention to pore occupancy and pore scan results. If necessary to resolve sequencing issues related to sample or flow cell quality:
1. Perform pore scan on flow cell and note if it resolves pore count issues.
 2. Stop the sequencing run, remove the flow cell from the PromethION, and insert it again into a different position. Restart the run and note if it resolves pore count issues.
 3. Stop the sequencing run. Recover DNA library and load onto a new flow cell, following protocol from step 139.
- 148 Following 72 hours of sequencing and 3 loads of prepared DNA libraries, the sample should yield a data output ~90-100 Gb with an N50 ~20-30 kb.

^ READ LENGTHS · OUTLIERS REMOVED

The read length graph shows the total number of bases vs the read length. The longest 1% of strands are classified as outliers, and excluded to allow focus on the main body of data.

N50*

29.05 kb

% Basecalled

100 %

Basecalled | Estimated




Figure 40. Read length distribution for cell sample with a 72 hour sequencing run

^ PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- Pore available
Pore in channel available for sequencing
- Reserved pore
Pore in reserve, will return to available when required
- Unavailable
Pore inhibited from sequencing
- Saturated
Possible contamination in the sample
- Zero
No current is passing through this pore, possibly due to bubbles on the membrane
- Inactive
Pore no longer suitable for further sequencing




Figure 41. PromethION pore scan results for cell sample with a 72 hour sequencing run and 3 loads of DNA library

Part 6: Washing and reloading the PromethION Flow Cell

149 Thaw the following reagents at room temperature:

- Sequencing Buffer (SB)
- Library Solution (LIS)

- Flow Cell Tether (FCT)
- Flow Cell Flush (FCF)
- Elution Buffer (EB)
- Wash Diluent (DIL)

Mix by vortexing, spin down, and place on ice.

- 150 Prepare the flow cell wash mix by combining  398 μL Wash Diluent and  2 μL Wash Mix per sample. Mix by pipetting. Do not vortex.

Note: Keep Wash Mix (WMX) in the freezer until use and return promptly. Do not vortex.

Note: Prepare the flow cell wash mix with a 10% overage.

- 151 Pause the sequencing run in MinKNOW.

- 152 If necessary, recover the already loaded DNA library:

- Slide the inlet port cover clockwise to open
- Set a P1000 pipette tip to 200 μL
- Insert the tip into the inlet port
- Turn the wheel until dial shows 400 μL

Note: Do not attempt to recover the loaded library after removing waste from the flow cell to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.

Note: DNA library recovery is necessary for samples with less than  900 ng of available DNA library for loading, and can be used for the third load in place a complete  300 ng DNA library load.

- 153 Remove waste buffer from flow cell:
- Close the inlet port
 - Insert a P1000 pipette into waste port 2 or 3 and remove the waste buffer

Note: It is vital that the inlet port is closed before removing waste to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.

- 154 Slide the inlet port cover clockwise to open. Draw back a small volume to remove any air bubbles:
- Set a P1000 pipette tip to 200 μL
 - Insert the tip into the inlet port

- Turn the wheel until dial shows 220-230 μL , until a small volume of buffer enters the pipette tip

Note: Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 μL , and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.

155 Load  400 μL of the flow cell wash mix into the flow cell via the inlet port, avoiding the introduction of air bubbles. Close the inlet port and ensure the light shield is installed.

156 Wait  01:00:00 .

157 Prepare the flow cell priming mix by combining  1170 μL Flow Cell Flush and  30 μL Flow Cell Tether per sample. Mix by vortexing.

158 Prepare  300 ng DNA library in  32 μL Elution Buffer in a 1.5 mL Eppendorf DNA LoBind tube.

Note: Recommended DNA input is 20 fmol, calculated to be 300 ng of 24kb dsDNA using <https://nebiocalculator.neb.com/#!/dsdnaamt>.

Note: If available sample mass for the third load is under 300 ng, DNA library should be recovered for the first load and used for the third load. Additional DNA library can be spiked in to the recovered DNA library.

159 Add  100 μL Sequencing Buffer (SB).

160 Add  68 μL Library Solution (LIS).

161 Slide the inlet port cover clockwise to open. Draw back a small volume to remove any air bubbles:

- Set a P1000 pipette tip to 200 μL
- Insert the tip into the inlet port
- Turn the wheel until dial shows 220-230 μL , until a small volume of buffer enters the pipette tip

Note: Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 μL , and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.



- 162 Load  500 μL of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles.
- 163 Wait  00:05:00 .
- 164 Load  500 μL of the priming mix into the flow cell via the inlet port, avoiding the introduction of air bubbles.
- 165 Mix the prepared library gently by pipetting up and down 5x just prior to loading.
- 166 Load  200 μL of library into the inlet port using a P1000 pipette.
- 167 Close the valve to seal the inlet port. Install the light shield. Close the PromethION door.
- 168 Wait  00:10:00 before resuming sequencing run in MinKNOW.
- 169 Within the first hour of sequencing, pay attention to pore occupancy and pore scan results. If necessary to resolve sequencing issues related to sample or flow cell quality:
1. Perform pore scan on flow cell and note if it resolves pore count issues.
 2. Stop the sequencing run, remove the flow cell from the PromethION, and insert it again into a different position. Note if it resolves pore count issues.
 3. Stop the sequencing run. Recover DNA library and load onto a new flow cell, following protocol from step 161.
- 170 Following 72 hours of sequencing and 3 loads of prepared DNA libraries, the sample should yield a data output ~90-100 Gb with an N50 ~20-30 kb.

^ READ LENGTHS · OUTLIERS REMOVED

The read length graph shows the total number of bases vs the read length. The longest 1% of strands are classified as outliers, and excluded to allow focus on the main body of data.

N50*

29.05 kb

% Basecalled

100 %

Basecalled | Estimated




Figure 42. Read length distribution for cell sample with a 72 hour sequencing run

^ PORE SCAN

A Pore scan is performed at configurable time intervals to determine the current status of pores within channels on a Flow Cell. For this run a Pore scan is performed every 1.5 hrs.

Legend

- Pore available
Pore in channel available for sequencing
- Reserved pore
Pore in reserve, will return to available when required
- Unavailable
Pore inhibited from sequencing
- Saturated
Possible contamination in the sample
- Zero
No current is passing through this pore, possibly due to bubbles on the membrane
- Inactive
Pore no longer suitable for further sequencing




Figure 43. PromethION pore scan results for cell sample with a 72 hour sequencing run and 3 loads of DNA library



Protocol references

Extracting HMW DNA using the Nanobind HT CBB kit for mammalian cultured cells on the KingFisher Apex system (PacBio): <https://www.pacb.com/wp-content/uploads/Procedure-checklist-Extracting-HMW-DNA-using-the-Nanobind-HT-CBB-kit-for-mammalian-cultured-cells-on-KingFisher-Apex-system.pdf>

Megaruptor 3 DNAFluid+ Kit Guide (Diagenode): <https://www.diagenode.com/files/products/kits/protocol-DNAFluid-V1.pdf>

High-Pass DNA Size Selection User Guide (Sage Science): <https://sagescience.com/wp-content/uploads/2019/01/BPLUS10-user-guide-460045-Rev-E.pdf>

Ligation sequencing DNA V14 (SQK-LSK114) library preparation (Oxford Nanopore): <https://nanoporetech.com/document/genomic-dna-by-ligation-sqk-lsk114?device=PromethION>

Washing and reloading the PromethION Flow Cell (Oxford Nanopore): <https://nanoporetech.com/document/flow-cell-wash-kit-exp-wsh004>