

Sep 06, 2023

## Processing of LRRK2-RCKW:GZD-824:E11

DOI

[dx.doi.org/10.17504/protocols.io.81wgbxz81lpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbxz81lpk/v1)

Amalia Villagran Suarez<sup>1</sup>

<sup>1</sup>Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093



Amalia Villagran Suarez

Department of Cellular and Molecular Medicine, University of...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbxz81lpk/v1>

**Protocol Citation:** Amalia Villagran Suarez 2023. Processing of LRRK2-RCKW:GZD-824:E11 . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.81wgbxz81lpk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** September 05, 2023



**Last Modified:** September 06, 2023

**Protocol Integer ID:** 87409

**Keywords:** protocol for processing, protocol, e11, lrrk2, rckw, lrrk2rckw, processing

**Funders Acknowledgements:**

ASAP

Grant ID: ASAP-000519

## Abstract

Protocol for processing of LRRK2<sup>RCKW</sup> bound to GZD-824 and DARPin-E11. This protocol covers everything from preprocessing to refinement.

## Materials

Requires the use of CryoSparc v.4.3.1 and Topaz 0.2.5

## Troubleshooting

## Preprocessing data

- 1 Use your preferred software. The original publication used MotionCor 2 and CTFFIND4. As a note, all of the data processed was collected on a UltrAuFoil Holey Gold 2/2 200 mesh grids and we used CryoSparc v.4.3.1 to process the data.

## Particle picking

### 2 **Blob picker**

We used cryoSparc-live blob picker option with a minimum particle diameter of 90Å to 250Å, lowpass filter to 20Å.

#### Note

Depending on the LRRK2<sup>RCKW</sup> variant, trimer particles might appear, we process them parallel using the following:

- Monomer particles were picked with a 320pixel box size with a pixel of 0.935Å. This were bin by 4 to a box size of 80pixel with a final pixel size of 3.74Å
- Trimer particles were picked with a 400pixel box size with a pixel of 0.935Å. This were bin by 5 to a box size of 80pixel with a final pixel size of 4.68Å.

### 3 **2D Classification and Topaz picking**

Using the extracted particles, we run 2D classification jobs in cryoSparc, using 40 iterations to clean our particle dataset. Selected particles then were use to train a Topaz model, and model was used to pick particles, followed by rounds of 2D classification. Selected particles can then be re-extracted to original pixel size (in this case 0.935Å, box size of 320pixel for monomers and 400pixel for dimers)

#### Note

Selecting 2D classes depends on the user. Take into account shape of the particles (J-shaped) and the visualization of secondary structure elements.

## Refinement for monomer

### 4 **Ab-initio reconstruction and Heterogeneous Refinement**



Select monomer particles then are subject to an Ab-initio reconstruction in cryoSparc asking 3 classes. Followed by a Heterogeneous Refinement.

#### Note

In this section you could try more than 3 classes, in order to sort bad particles out. This can be coupled with the final 2D classification jobs from **step 3**.

## 5 **Nu-Refinement and Local refinement**

Selected refinement from **step 4**, is subject to a Nu-refinement at C1 symmetry. Followed by a local refinement with a continuous mask surrounding the Kinase, WD40 and DARPin-E11.

## 6 **3D Variability**

To account for the heterogeneity of the ROC-COR domains, create a mask using the Nu-refinement of **step 5**, and using the same amount of particles of the same job, run a 3DVA.

#### Note

Recommendation: use a filter resolution of 7Å, and we ask for 3 components.

### 6.1 **3D Variability Display**

- Use this jobs with the output from Step 6, to visualize the movement of the ROC-COR domain. Input the following: for **output-mode "intermediates"**, **"11: frames"**, and **filter resolution of "7Å"**.
- To visualize the movement of each component download file to chimera and use volume morph. Select component with the most movement in the ROC-COR domain.
- Clone 3D Variability Display job and now select: **skip reconstruction, only used this component** (write which component you picked), and select **Intermediates: output particle subsets**. The output of this job should be the 11 particles stacks.

## 7 **Local Refinement of the ROC-COR domain**

- Particles of the component obtained from step 6.1, are split in half, selecting the particles that account for the movement of the ROC-COR closest to the Kinase domain.



- This are subject to a local refinement. The Nu-refinement volume from **step 5** is used as the input for the local refinement, a continuous mask surrounding the ROC-COR domains is created and the selected particles of the component.

## Refinement for trimer

### 8 **Ab-initio and Nu-Refinement**

- Select trimer particles were use to run an Ab-initio reconstruction in cryoSparc. Followed by a Nu-Refinement C1 and at C3 symmetry.
- We compare both the C1 and C3 maps to check for density for the inhibitors in the protomers and overall resolution.

#### Note

In our case, C3 symmetry improved resolution of Kinase active site.

### 9 **Particle expansion**

Particles from Nu-Refinement in **step 8** were used for a **Symmetry Expansion** job in cryoSparc, based on the volume symmetry.

### 10 **Local Refinement**

A local refinement was carried using the following inputs: mask of one protomer of the trimer, expanded particles from **step 9**, Nu-refinement from **step 8** and C1 symmetry.