

Feb 25, 2025

Processing plasma from Whole Blood

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1d19kvr2/v1>

Sara Lucas Del Pozo^{1,2}, Kai-Yin Chau^{1,2}, Jane Macnaughtan^{1,2}, Giuseppe Uras^{1,2}, Roxana Mezabrovschi^{1,2}

¹University College London; ²Aligning Science Across Parkinson's



Sofia Koletsi

University College London, University of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1d19kvr2/v1>

Protocol Citation: Sara Lucas Del Pozo, Kai-Yin Chau, Jane Macnaughtan, Giuseppe Uras, Roxana Mezabrovschi 2025. Processing plasma from Whole Blood. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov1d19kvr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 05, 2025

Last Modified: February 25, 2025

Protocol Integer ID: 119993

Keywords: ASAPCRN, plasma from whole blood, processing plasma, processing of plasma, whole blood sample, plasma, whole blood, processing, blood, protocol detail, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000420

Abstract

This protocol details the processing of plasma from the whole blood sample.

Guidelines

The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.

Safety warnings

 The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.

Before start

The protocol needs prior approval by the users' Institutional Review Board (IRB), Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s) as applicable.



Procedure

40m

- 1 Gently invert the plasma tubes at 180° and back 8-10 times immediately after collection. Allow 00:20:00 - 00:30:00 to cool down prior to sample processing. 30m
 - 2 Centrifuge the green cap tube(s) in at 1800 x g, Room temperature, 00:10:00 **NO BRAKE**. 10m
 - 3 After centrifugation there should be a clear separation between the red blood cells at the bottom and the plasma at the top.
 - 4 Use the Pasteur pipette provided or a micropipette to carefully aliquot 500 µL of the plasma from the green cap tube (Lithium Heparin) into the 2-3 × 1ml cryovials (500 µL into each tube) labelled as LH-PLA.
- Note**

Do not to disturb the red blood cell layer. If the sediment is disturbed the centrifugation step must be repeated.
- 5 Place plasma samples in -80 °C .