

Apr 24, 2023

Promega Wizard DNA extraction - Drosophila whole body protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwjrzdlmk/v1>

Finley Grover-Thomas¹

¹UCL



Finley Grover-Thomas

University College London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjrzdlmk/v1>

Protocol Citation: Finley Grover-Thomas 2023. Promega Wizard DNA extraction - Drosophila whole body protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwjrzdlmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 21, 2023



Last Modified: April 24, 2023

Protocol Integer ID: 80886

Keywords: Promega, DNA extraction, Wizard, Drosophila, Insect, promega wizard genomic dna purification, promega wizard dna extraction, drosophila whole body protocol, body fruit fly sample, drosophila, dna extraction, genomic dna, dna, purification, whole body protocol, extraction

Abstract

Edited version of the 'Promega Wizard Genomic DNA Purification' protocol optimised to work on pooled whole-body fruit fly samples.



Lysis & Incubation

- 1 Add **10 flies into a 1.5ml Eppendorf tube** and cover with  300 µL Nuclei Lysis Solution . 2m
- 2 Use a **micro-tube pestle to crush** the flies until they are fully ground. 5m


- 2.1 For the initial crushing, use your hands and the pestle only.

Once the flies are partially disintegrated an electronic homogeniser can be used to increase efficiency.

Take care not to froth the solution up too much.
- 3 Once the flies are thoroughly homogenised, add another  300 µL of Nuclei Lysis Solution . 2m

By this point, the solution should be opaque and red-brown in colour.
- 4 **Vortex** the homogenised solution then **incubate** at  65 °C for 20 minutes . 20m
- 5 Remove the samples from the heat and lower the set temperature to  55 °C .
- 6 Allow for the **samples to cool** for a few moments, then add  17.5 µL Proteinase K . 2m

(This is not included in the Promega Wizard kit).
- 7 **Vortex the samples** thoroughly then **incubate** at  55 °C for 3 hours . 3h

Vortex the samples regularly throughout this time.



Protein Precipitation

4m

- 8 Remove the samples from the heat block and **allow to cool** to room temperature.
- 9 Add  200 μ L Protein Precipitation Solution to each sample & **vortex** for **20 seconds** (until thoroughly mixed).
- 10 Leave the samples  On ice for 5 minutes .
- 11 **Centrifuge** the samples at high speed:  13200 rpm , for 4 minutes .
- 12 **Carefully transfer the supernatant** (liquid) to a new (labelled) micro-centrifuge/Eppendorf tube.

Take care not to disturb the pellet.
- 13 *If any tissue or protein precipitate (white mass) remains in the supernatant, **repeat steps 11 & 12.***
- 14 Once the supernatant is clear of tissue mass and protein precipitate, add  600 μ L isopropanol (at room temperature).

5m

4m

2m

2m

DNA elution and cleansing

31m

- 15 **Gently invert** the samples to mix the isopropanol and supernatant.

White threads of DNA may or may-not form; if no threads are visible after ~5 minutes, continue on with the protocol.
- 16 **Centrifuge** the samples at  13200 rpm , for 1 minute .
- 17 Taking care not to disturb the pellet (of DNA), **remove the supernatant.**

If the supernatant does not contain the pellet, it can be discarded.
- 18 Add  600 μ L 70% ethanol , then **gently invert** the tube several times (to wash the DNA).

5m

1m

4m

3m



19 **Centrifuge** at  13200 rpm , for 1 minute .

1m

20 Making sure that the DNA pellet is not disturbed, **remove the ethanol** (this can be discarded).

2m

21 **Invert the (open) tube** onto clean absorbent paper; **leave for 10-15 minutes**.

15m

Towards the end of this time, set the heating bloc to  65 °C .

22 Add  100 µL DNA Rehydration Solution **(Nuclease Free Water)** and **incubate** at  65 °C for 1 hour .

Mix the solution by gently tapping and shaking the tubes.

23 Store the DNA at  2-8 °C .

Protocol references

Adapted from:

https://www.promega.co.uk/-/media/files/resources/protocols/technical-manuals/0/wizard-genomic-dna-purification-kit-protocol.pdf?rev=776154d180fe4b41a0a3dc8a9e328c8a&sc_lang=en