

Jul 21, 2025

Protein Digestion and Mass Spectrometry Analysis

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9o67qv3e/v1>

Mara monetti¹, Zhuoning Li¹

¹Proteomics Core Facility, Sloan Kettering Institute, Memorial Sloan Kettering Cancer Center



Michael G Hanna

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9o67qv3e/v1>

Protocol Citation: Mara monetti, Zhuoning Li 2025. Protein Digestion and Mass Spectrometry Analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9o67qv3e/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2025



Last Modified: July 21, 2025

Protocol Integer ID: 219694

Keywords: ASAPCRN, mass spectrometry analysis this protocol, mass spectrometry analysis

Abstract

This protocol describes ...

Protein Digestion

- 1 The beads were resuspended in 80 μL of 2M Urea, 50 mM ammonium bicarbonate (ABC) and treated with DL-dithiothreitol (DTT) (final concentration 1 mM) for 30 minutes at 37°C with shaking at 1100 rpm on a Thermomixer (Thermo Fisher).
- 2 Free cysteine residues were alkylated with 2-iodoacetamide (IAA) (final concentration 3.67 mM) for 45 minutes at 25°C with shaking at 1100 rpm in the dark.
- 3 The reaction was quenched using DTT (final concentration 3.67 mM), and LysC (750 ng) was added, followed by incubation for 1h at 37°C at 1150 rpm.
- 4 Trypsin (750 ng) was added, and the mixture was incubated for 16 hours at 37°C with shaking at 1150 rpm.
- 5 After incubation, an additional 500 ng of trypsin was added to the sample, followed by a 2-hour incubation at 37°C at 1150 rpm.
- 6 The digest was then acidified to pH <3 by adding 50% trifluoroacetic acid (TFA), and the peptides were desalted on C18 stage tips (Empore C18 extraction disks).
- 7 Briefly, the stage tips were conditioned with sequential additions of:
 - 7.1 i) 100 μL methanol,
 - 7.2 ii) 100 μL 70% acetonitrile (ACN)/0.1% TFA,
 - 7.3 iii) 100 μL 0.1% TFA twice.
- 8 After conditioning, the acidified peptide digest was loaded onto the stage tip, and the stationary phase was washed with 100 μL 0.1% formic acid (FA) twice.
- 9 Peptides were eluted with 50 μL 70% ACN/0.1% FA twice.

- 10 Eluted peptides were dried under vacuum in a Speed-Vac centrifuge, reconstituted in 12 μL of 0.1% FA, sonicated and transferred to an autosampler vial.
- 11 Peptide yield was quantified using a NanoDrop (Thermo Fisher).

Mass Spectrometry Analyses

- 12 Peptides were separated on a 25 cm column with a 75 μm diameter and 1.7 μm particle size, composed of C18 stationary phase (IonOpticks Aurora 3 1801220) using a gradient from 2% to 35% Buffer B over 90 minutes
- 12.1 This was followed by an increase to 95% Buffer B for 7 minutes (Buffer A: 0.1% FA in HPLC-grade water; Buffer B: 99.9% ACN, 0.1% FA) with a flow rate of 300 nL/min on a NanoElute2 system (Bruker).
- 13 MS data were acquired on a TimsTOF HT (Bruker) with a Captive Spray source (Bruker) using a data-independent acquisition PASEF method (dia-PASEF). The mass range was set from 100 to 1700 m/z, and the ion mobility range from 0.60 V.s/cm² (collision energy 20 eV) to 1.6 V.s/cm² (collision energy 59 eV), a ramp time of 100 ms, and an accumulation time of 100 ms. The dia-PASEF settings included a mass range of 400.0 to 1201.0 Da, mobility range 0.60-1.60, and an estimated cycle time of 1.80 seconds. The dia-PASEF windows were set with a mass width of 26.00 Da, mass overlap 1.00 Da, and 32 mass steps per cycle.

DIA Data Analysis

- 14 Raw data files were processed using Spectronaut version 17.4 (Biognosys) and searched with the PULSAR search engine against the Homo sapiens UniProt protein database (226,953 entries, downloaded on 2022/09/23). Cysteine carbamidomethylation was set as fixed modifications, while methionine oxidation, acetylation of the protein N-terminus, and deamidation (NQ) were defined as variable modifications. A maximum of two trypsin missed cleavages was allowed.
- 15 Searches used a reversed sequence decoy strategy to control the peptide false discovery rate (FDR), with a 1% FDR threshold for identification. An unpaired t-test was used to calculate p-values for differential analysis, and volcano plots were generated based on log₂ fold change and q-value (multiple testing corrected p-value). A q-value of ≤ 0.05 was considered the statistically significant cut-off.