

Oct 28, 2020

## 🌐 Protein expression in *Komagataella phaffii* (formerly *Pichia pastoris*)

📁 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bd32i8qe](https://dx.doi.org/10.17504/protocols.io.bd32i8qe)


Maja Rennig<sup>1</sup>, Kristoffer Bach Falkenberg<sup>1</sup>, Cristina Hernandez Rollan<sup>1</sup>, Andreas Birk Bertelsen<sup>1</sup>, Morten Norholm<sup>1</sup>

<sup>1</sup>Technical University of Denmark



**Maja Rennig**

Technical University of Denmark



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bd32i8qe](https://dx.doi.org/10.17504/protocols.io.bd32i8qe)

**Protocol Citation:** Maja Rennig, Kristoffer Bach Falkenberg, Cristina Hernandez Rollan, Andreas Birk Bertelsen, Morten Norholm 2020. Protein expression in *Komagataella phaffii* (formerly *Pichia pastoris*). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bd32i8qe>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 23, 2020

**Last Modified:** October 28, 2020

**Protocol Integer ID:** 34650

**Keywords:** LyGo, LPMO, *K. phaffii*, *P. pastoris*, protein production,

## Abstract

The methylotrophic yeast *Komagataella phaffii* (formerly *Pichia pastoris*) is the most commonly used yeast species in the production of recombinant proteins. This is most likely due to its ability to grow to high cell density, to express recombinant genes in a tightly controlled manner and to efficiently secrete proteins. Despite its biotechnological importance and wide use in industry, relatively few genetic tools are readily available for academic research. Here, we present a protocol for production of proteins in *K. phaffii*.

The protocol found here is modified from the protocol provided in the *Pichia* Expression Kit (Catalog Number K1710-01, Invitrogen).

## Guidelines


Transformation and integration into *K. phaffii* GS115 can be done utilizing the *Pichia* EasyComp Kit (Catalog no. K1730-01, Invitrogen).

Recipes for all media and solutions can also be found in the manual for the *Pichia* Expression Kit (Catalog Number K1710-01, Invitrogen).

## Materials

### MATERIALS

 Yeast Extract **Catalog #Y1625**

 GS115, *Pichia pastoris* Yeast Strain **Thermo Fisher Catalog #C18100**

 YNB without Amino Acids **Sigma Aldrich Catalog #Y0626-1KG**

 Biotin **Sigma Aldrich Catalog #B4639-1G**

 Peptone from soybean enzymatic digest **Sigma Aldrich Catalog #70178-500G**


 Potassium phosphate monobasic **Sigma Aldrich Catalog #P0662-1KG**

 Potassium phosphate dibasic **Sigma Aldrich Catalog #P3786-1KG**

 Methanol **Sigma Aldrich Catalog #34860-2.5L-M**

 Glycerol **Sigma Aldrich Catalog #15523-1L-M**

## Safety warnings

 This protocol describes the use of GMO classified organisms. Make sure that the local GMO and safety legislations are respected.

This protocol involves the use of 100 % methanol. Make sure to store and handle 100 % methanol accordingly.

## Before start

Make sure to have your expression strain streaked on an YPD agar plate

## Prepare media and stock solutions

1 **Prepare** [M] 1 Molarity (M) **potassium phosphate buffer,**  $\text{pH } 6$

1.1 Combine  $\text{132 mL}$  of [M] 1 Molarity (M)  $\text{K}_2\text{HPO}_4$  and  $\text{868 mL}$  of [M] 1 Molarity (M)  $\text{KH}_2\text{PO}_4$

1.2 Confirm that the  $\text{pH } 6.0$  and adjust using phosphoric acid or KOH

1.3 Sterilize by autoclaving for  $\text{00:20:00}$  on liquid cycle program and store at  $\text{Room temperature}$

### Note

The shelf life of this solution is greater than one year

2 **Prepare** [M] 13.4 Mass / % volume **YNB with ammonium sulfate without amino acids**

2.1 Dissolve  $\text{134 g}$  of yeast nitrogen base (YNB) with ammonium sulfate and without amino acids in  $\text{1.000 L}$  of water.

The solution can be heated to dissolve YNB completely in water

### Note



Alternatively, use  $\text{34 g}$  of YNB without ammonium sulfate and amino acids and  $\text{100 g}$  of ammonium sulfate.

2.2 Sterilize by filtration and store at  $\text{4 }^\circ\text{C}$

Note

The shelf life of this solution is approximately one year

3 **Prepare** [M] 0.2 Mass / % volume **biotin solution**

3.1 Dissolve  20 mg biotin in  100 mL of water

3.2 Sterilize by filtration and store at  4 °C

Note

The shelf life of this solution is approximately one year

4 **Prepare BMGY and BMMY media**

**Ingredients (final concentrations):**

[M] 1 Mass / % volume yeast extract

[M] 2 Mass / % volume peptone

[M] 100 millimolar (mM) mM potassium phosphate, pH 6.0

[M] 1.34 Mass / % volume YNB


[M] 0.0004 Mass / % volume biotin

[M] 1 % volume glycerol (BMGY) or [M] 0.5 % volume methanol (BMMY)




Note



*K. phaffii* cells exhibit optimal growth with higher YNB concentrations than *S. cerevisiae*

4.1 Dissolve  10 g of yeast extract,  20 g of peptone in  700 mL of water



4.2 Autoclave for  00:20:00 on liquid a cycle program

4.3 Let cool down to room temperature, then add the following:

 100 mL  1 Molarity (M) potassium phosphate buffer,  6 (autoclaved stock solution)

 100 mL  13.4 Mass / % volume YNB with ammonium sulfate without amino acids (filter sterilized stock solution)

 2 mL  0.02 Mass / % volume biotin (filter sterilized stock solution)


4.4 For BMGY add  100 mL  10 % volume glycerol

For BMMY add  100 mL  5 % volume methanol

#### Note

If not used directly, store the media at  4 °C for up to 2 months

## Inoculation - Day 1




5 Inoculate  25 mL of BMGY in a 250 mL shake flask with a single colony of *K. phaffii* GS115 with the desired construct integrated on the genome.

6 Incubate the culture at  28 °C with 250RPM shaking until saturation

#### Note

It is also possible to grow the cells at  30 °C

## Pre-culture - Day 2

7 Transfer  2.5 mL of the saturated culture into  100 mL BMGY in a 500 mL shake flask with low baffles and grow overnight at  28 °C

Note

It is also possible to grow the cells at  $30\text{ }^{\circ}\text{C}$

## Expression - Day 3

8 Measure OD<sub>600</sub> of the overnight culture

9 Pellet  $50\text{ mL}$  of the culture in a 50 mL falcon tube at

$4000\text{ x g}$ , Room temperature, 00:10:00

10 Discard the supernatant and resuspend the pellet in BMMY medium to a final OD<sub>600</sub> of 50. The appropriate amount of media can be calculated as:

$$V_{BMMY} = \frac{OD_{600,measured} \cdot V_{pre-culture}}{50}$$

11 Inoculate  $100\text{ mL}$  BMMY medium with  $2\text{ mL}$  of the normalized pre-culture, giving a start OD<sub>600</sub> of 1.

12 Grow the cells at  $28\text{ }^{\circ}\text{C}$  and 250 RPM of shaking




Note

It is crucial that the temperature does not exceed  $28\text{ }^{\circ}\text{C}$  when producing protein. It is advised to double-check the incubator's temperature reading with a thermometer.




## Feeding - Day 4

13 Add 1% methanol to the expression culture ( $1\text{ mL}$  of 100 % volume methanol to  $100\text{ mL}$  of BMMY culture) around 24 hours after inoculation





## Feeding - Day 5

- 14 Add 1% methanol to the expression culture (  1 mL of  100 % volume methanol to  100 mL of BMMY culture) around 24 hours after last addition of methanol

## Feeding - Day 6

- 15 Add 1% methanol to the expression culture (  1 mL of  100 % volume methanol to  100 mL of BMMY culture) around 24 hours after last addition of methanol

## Harvesting - Day 7

- 16 Harvest the cells by centrifugation at  4000 x g, 4°C, 00:10:00 and discard the pellets
- 17 Keep the sample  On ice when working with it and at  4 °C or  -20 °C for storage