

Oct 20, 2019

Protein expression in OnePot PURE cell-free system

DOI

dx.doi.org/10.17504/protocols.io.8avhse6



Konstantinos Ragios¹

¹EPFL - EPF Lausanne

iGEM EPFL



Konstantinos Ragios

EPFL - EPF Lausanne

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.8avhse6

External link: <https://pubs.acs.org/doi/10.1021/acssynbio.8b00427>

Protocol Citation: Konstantinos Ragios 2019. Protein expression in OnePot PURE cell-free system. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8avhse6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2019

Last Modified: October 20, 2019

Protocol Integer ID: 28725

Keywords: iGEM, cell-free, OnePot, PURE, protein, expression,

Abstract

In this protocol we explain the procedure of protein expression in the OnePot PURE cell-free system.

Materials

Material/Consumables:

- OnePot Protein solution
- OnePot Ribosome solution
- OnePot Energy solution
- DNA template
- Nuclease free water

	Com pone nt	Volu me (μL)	Final conc entra tion of the reacti on (nM)
	Protei n Soluti on	0.65	-
	Ribos ome Soluti on	0.9	-
	Energ y Soluti on	2	-
	DNA Temp late	x	5
	Water	1.45- x	-

Before start

In order to create the components needed refer to the following protocols:

Protocol



NAME

Protein Purification for OnePot PURE cell-free system

CREATED BY

Konstantinos Ragios

PREVIEW

Protocol



NAME

Ribosome Purification for OnePot PURE cell-free system

CREATED BY

Konstantinos Ragios

PREVIEW

Protocol



NAME

Energy solution preparation for OnePot PURE cell-free system

CREATED BY

Konstantinos Ragios

PREVIEW

Protein expression

- 1 For a 5 μ l reaction add to a tube 2 μ l of Energy solution, 0.65 μ l of Protein solution, 0.9 μ l of Ribosome solution, the DNA template (5nM final concentration in the reaction) and if needed add water to reach the final volume.

Note

The minimum reaction volume is 5 μ l while the suggested one is 10 μ l.
For any final reaction volume (e.g. X μ l) you just need to multiply the volumes of the components needed by a factor of X/5.

Note

If you haven't produce ribosomes then any store-bought product can be used.

- 2 Centrifuge for a few seconds and then your solution will start reacting.

- 3

Note

If you want to measure it with a plate reader, ideally you should add each component in different corners of the well, so the reaction will not initiate before adding the plate on the reader. The centrifuge for a few seconds.